Test your health

PRIMA LAB SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
+41 91 605 1030
support@primalabsa.ch
primalabsa.ch

EC REP Qarad EC-REP BV
Pas 257
2440 Geel - Belgium

800025_IFU_11_5.0 06/2022

REF 800025

Registered trademarks.

INSTRUCTIONS FOR USE GEBRAUCHSANWEISUNG NOTICE D'UTILISATION ISTRUZIONI PER L'USO

IRON FER TEST

Self-test for the detection of ferritin levels
in whole blood samples

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ENGLISH IRON AND ANAEMIA

Iron is an essential metal for our bodies and is crucial for transporting oxygen in the blood, for cell multiplication and to build the structure of tissues and organs. However, in excessively high levels it is toxic for the body. For this reason, every one of us has a system for taking up iron from the external environment (e.g. via a diet rich in iron-containing foods) and storing it in cells in a way that is not excessive (and therefore not toxic). Ferritin is the protein responsible for this storage function. The level of ferritin is an excellent indicator of the amount of iron available to the body. Low levels of this protein in the blood are an indication of depleted iron stores, a condition that precedes the development of anaemia. A decrease can be caused by anaemia, pregnancy, haemorrhages, alterations in iron uptake, tuberculosis.

TEST PRINCIPLE


Iron FER Test is an immunochromatographic assay which detects the protein Ferritin thanks to special monoclonal gold-conjugate antibodies embedded to test strip.

REFERENCES

- Wick M, Pingerra W, Lehmann P. *Iron metabolism: diagnosis and therapy of anemias*, 5th ed, Vienna, New York: Springer Verlag, 2003; p. 151.
- Worwood M. *The laboratory assessment of iron status – an update. Clin Chim Acta* 1997; 259: 3-23.
- Kaltwasser JP, Werner E. *Diagnosis and clinical evaluation of iron overload. Baillieres Clin Haematol* 1989; 2; 363-89.
- Baynes RD, Cook JD. *Current issues in iron deficiency. Curr Opin Hematol* 1996; 3:145-9.
- Lee MH, Means RT Jr. *Extremely elevated serum ferritin levels in a university hospital: associated diseases and clinical significance. Am J Med* 1996; 98: 566-71.

CONTENT

* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the external box.

* Material required but not supplied: cotton wool, a device to measure time (i.e. timer, watch).

* Do not open the sealed aluminium bag until just before performing the test. Take care to open it as marked.

* The desiccant packet must not be used. Dispose of it with household waste without opening it.

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PRECAUTIONS

- Read these instructions for use carefully before performing the test. The Test is reliable only if all the instructions are followed correctly.
- Keep the Test out of the reach of children.
- Do not use the Test after the expiry date or if the package has been damaged.
- Follow the procedure exactly, using only the specified quantities of blood and diluent.
- Store the Test components at a temperature between +4°C and +30°C. Do not freeze.
- Use the test and lancing device once only.
- The test is for external use only. DO NOT SWALLOW.
- In vitro diagnostic device for individual use.
- Not recommended for people who take anti-coagulant medications (blood thinners) or people suffering from haemophilia.
- After using, please dispose of all components according to your local waste disposal laws. Ask your pharmacist for advice.

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