## Test your health

**PRIMA LAB SA**  
Via Antonio Monti 7  
CH-6828 Balerna - SWITZERLAND  
support@primallabs.ch  
primallabs.ch

EC REP Qarad EC-REP BV  
Pas 257  
2440 Geel - Belgium

100077\_IFU\_11\_5.3\_01/2023 REF 100077

Registered trademarks

### INSTRUCTIONS FOR USE GEBRAUCHSANWEISUNG NOTICE D'UTILISATION ISTRUZIONI PER L'USO

## CELIAC TEST

*Self-test to check for the presence of anti deamidated gliadin (anti-DGP) antibodies, IgA and IgG, in the blood, for celiac disease screening.*

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### ENGLISH GLIADIN AND CELIAC DISEASE


Celiac disease is an autoimmune disease that causes chronic inflammation of the small intestine. Disease symptoms include diarrhoea, weight loss, abdominal bloating, malnutrition and skin reactions. This inflammation is caused by gliadin, the protein component of gluten, contained in wheat, barley, rye, spelt and many other cereals. Recent studies have shown that subjects affected by celiac disease produce specific antibodies (class IgA and IgG) against deamidated gliadin peptides (anti-DGP). These antibodies have high specificity and sensitivity in celiac disease, which permits very accurate diagnosis. Early diagnosis can reduce the progression of intestinal damage and prevent more serious complications. The Celiac Test is a rapid immunochromatographic test to reveal the presence of class IgA and IgG anti-DGP antibodies in whole blood for celiac disease screening.

### REFERENCES

1. P.D. Mooney, S.H.Wong, A.J.Johnston, M.Kurien, A.Averginos and D.S.Sanders "Increased detection of celiac disease with measurement of Deamidated Gliadin Peptide Antibody before endoscopy" *Cli.Gastroenterology and Hepatology* 2015;13:1278-1284
2. N. R. Lewis & B. B. Scott "Meta-analysis: deamidated gliadin peptide antibody and tissue transglutaminase antibody compared as screening tests for coeliac disease" *Alimentary Pharmacology & Therapeutics*
3. A. Lammi, P. Arikoski, S. Simell, T. Kinnunen, V. Simell, S. Paavanan-Huhtala, A. Hinkkanen, R. Veijola, M. Knip, J. Toppari, O. Vaarala, O. Simell, and J. Ilonen "Antibodies to Deamidated Gliadin Peptide in Diagnosis of Celiac Disease in Children" *Gastroenterology - JPN* 2015;60: 626-631

### CONTENT

\* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the external box.

\* Material required but not supplied: cotton wool, a device to measure time (i.e. timer, watch).  
\* Do not open the sealed aluminium bag until just before performing the test. Take care to open it as marked.  
\* The desiccant packet must not be used. Dispose of it with household waste without opening it.

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
Vitrex Medical A/S  
Vasekaer 6-8  
DK-2730 Herlev  
Denmark

SterilLance Medical (Suzhou) Inc.  
No. 168 PutuoShan Road  
New District  
215153 Suzhou, Jiangsu - P. R. China

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### PRECAUTIONS

1. Read these instructions for use carefully before performing the test. The Test is reliable only if all the instructions are followed correctly.
2. Keep the Test out of the reach of children.
3. Do not use the Test after the expiry date or if the package has been damaged.
4. Follow the procedure exactly, using only the specified quantities of blood and diluent.
5. Store the Test components at a temperature between +4°C and +30°C. Do not freeze.
6. Use the test and lancing device once only.
7. The test is for external use only. DO NOT SWALLOW.
8. In vitro diagnostic device for individual use.
9. Not recommended for people who take anti-coagulant medications (blood thinners) or people suffering from haemophilia.
10. After using, please dispose of all components according to your local waste disposal laws. Ask your pharmacist for advice.

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