SPECIMEN COLLECTION AND PREPARATION

Before testing, read the instructions carefully and completely. Before testing, it is important that you refrain from any sexual activity for 3-7 days.

- 1) Before testing, it is important that you refrain from any sexual activity for 3-7 days. This ensures that the volume and quality of sperm is at its peak and the test will then be an accurate determination of sperm concentration.
- 2) Using masturbation, the semen should be collected directly into the collection cup. **ATTENTION: Do not remove the liquefying powder from the cup since it will help the semen to liquefy.**
- 3) Care should be taken that collected semen is not contaminated by touch of hands or tissues or any other materials.
- 4) Shake the semen evenly in the collection cup and leave it to stand for 15 minutes at room temperature until the semen liquefies. Do not use semen stored for more than 12 hours after liquefaction.

TEST PROCEDURE:

- 1) Remove the test cassette from the foil pouch and lay it horizontally on a flat surface. Using the pipette provided in the foil pouch, dispense one drop of semen into sample well A. (FIG. A-B)
- 2) Once the semen is soaked into Well A, add three drops of the blue staining solution to Test Well A. Let it soak for 1-2 minutes. (FIG. C)
- 3) Now add two drops of the transparent washing solution to Test Well A (FIG. D), and let it soak for 1-2 minutes, and then read the results immediately.
- 4) Read the colour of Test Well A, comparing the colour of Test Well A to Reference Well B. The darker the colour of Well A, the higher the sperm concentration.


RESULTS INTERPRETATION

Compare the colour of test well A to reference Well B

 2 min.


NORMAL

The colour of test well A is darker than or the same as the standard colour of reference Well B. It means that the sperm concentration is greater than or equivalent to 15 million/ml. The likelihood of conception is high with this sperm concentration, subject to other conditions, such as ovulation time of the partner being favorable.

ABNORMAL

The colour of test well A is lighter than the standard colour of reference Well B. It means that the sperm concentration is less than 15 million/ml. This is known as oligospermia (a range that is normally between 5 million/ml and 15 million/ml). The likelihood of conception is less with this sperm concentration and further medical consultation is recommended.

NOTE: If well A is colourless, it means the sperm concentration is less than 5 million/ml or zero. This condition is known as azoospermia. If you are unsure of the result or you feel the result is abnormal you should repeat the test using the second test that is included in the pack but make sure you do not ejaculate through any sexual activity for 6 days before carrying out the second test. If the second test is also abnormal, you should discuss the results with your doctor or medical professional.

F.A.Q. - QUESTIONS AND ANSWERS


WHY THE TIME IS 15 MINUTES BEFORE THE SEMEN BEING TAKEN OUT FROM THE SEMEN COLLECTING CUP FOR TEST, AND THE STORAGE TIME AFTER SAMPLING DOES NOT EXCEED 12 HOURS? The fresh semen is viscous, and normal semen needs to be incubated for 30-60 minutes at 37 °C to liquefy completely. Only in liquid state, the semen can be used for test, because the viscous semen cannot completely pass the membrane of the test well. The sheet at bottom of the collection cup can make the semen liquefy quickly within 15 minutes. If the storage time of sperm specimen is too long, it may cause lysis of sperms, which may affect accuracy of the results.

HOW LONG DO THE SEMEN AND STAINING SOLUTION NEED TO BE IN THE TEST WELL? In general, they may pass membrane of the test well within several seconds, if they cannot pass the membrane completely after 5 minutes, and which shows that the semen has not liquefied completely, or density of the sperm is too high, and we must repeat the test. The reasons of semen non-liquefaction may be that the enzyme in the collection cup becomes invalid or the user does

not use the semen cup correctly or there is some other cup used for semen collection. **IS THE SOLUTIONS USED FOR TEST SAFE?** The solution is safe, and which is a synthetic dye, intake of which in concentration less than 5mg/kg does not harm to human body. The concentration of the staining fluid is less than 10ug/ml. **CAN THE ABNORMAL RESULTS SHOW THAT THE SUBJECT HAS NO ABILITY TO HAVE CHILDREN?** Sperm concentration is one of the several semen analysis tests. There are other factors that should be taken into account such as motility. Therefore, a medical expert's opinion in such cases would be highly recommended. **WHAT IS THE REASON THAT MAY RESULT IN WRONG TEST RESULTS?** Any mistake at any point of time from sample collection to test timing to non-compliance to abstinence may result in erroneous test results. **HOW ACCURATE IS THE SPERM TEST?** The test is very accurate. Assessment reports show a concordance rate above 99% (95%: 97.0-99.9%) with reference methods.

not use the semen cup correctly or there is some other cup used for semen collection. **IS THE SOLUTIONS USED FOR TEST SAFE?** The solution is safe, and which is a synthetic dye, intake of which in concentration less than 5mg/kg does not harm to human body. The concentration of the staining fluid is less than 10ug/ml. **CAN THE ABNORMAL RESULTS SHOW THAT THE SUBJECT HAS NO ABILITY TO HAVE CHILDREN?** Sperm concentration is one of the several semen analysis tests. There are other factors that should be taken into account such as motility. Therefore, a medical expert's opinion in such cases would be highly recommended. **WHAT IS THE REASON THAT MAY RESULT IN WRONG TEST RESULTS?** Any mistake at any point of time from sample collection to test timing to non-compliance to abstinence may result in erroneous test results. **HOW ACCURATE IS THE SPERM TEST?** The test is very accurate. Assessment reports show a concordance rate above 99% (95%: 97.0-99.9%) with reference methods.

not use the semen cup correctly or there is some other cup used for semen collection. **IS THE SOLUTIONS USED FOR TEST SAFE?** The solution is safe, and which is a synthetic dye, intake of which in concentration less than 5mg/kg does not harm to human body. The concentration of the staining fluid is less than 10ug/ml. **CAN THE ABNORMAL RESULTS SHOW THAT THE SUBJECT HAS NO ABILITY TO HAVE CHILDREN?** Sperm concentration is one of the several semen analysis tests. There are other factors that should be taken into account such as motility. Therefore, a medical expert's opinion in such cases would be highly recommended. **WHAT IS THE REASON THAT MAY RESULT IN WRONG TEST RESULTS?** Any mistake at any point of time from sample collection to test timing to non-compliance to abstinence may result in erroneous test results. **HOW ACCURATE IS THE SPERM TEST?** The test is very accurate. Assessment reports show a concordance rate above 99% (95%: 97.0-99.9%) with reference methods.

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