## Your partner for rapid testing

**PRIMA Lab SA**  
Via Antonio Monti 7  
CH-6828 Balerna - SWITZERLAND

support@primalabsa.ch  
primalabsa.ch

REF 800060-xP

800060-xP\_IFU\_10\_5.5 06/2023

EC REP Qarad EC-REP BV  
Pas 257, 2440 Geel  
Belgium

PRIMA®

Registered trademarks

## STREP A RAPID TEST

Rapid Test for the detection of group A β-hemolytic Streptococcus in pharyngeal swab.  
IVD for professional use only.

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### ENGLISH STREPTOCOCCUS A

The throat infection with group A β-hemolytic Streptococcus (also known as *Streptococcus pyogenes*) is a Gram-positive bacterium that constitutes the most frequent bacterial cause of acute pharyngitis and it is also responsible for some skin infections such as impetigo and erysipelas. Most sore throats, on the other hand, are caused by viral infections, which are cured without the need of antibiotic treatments, unnecessary in these cases. Streptococcus A infections can resolve in a few days without treatment, although doctors generally prefer to prescribe antibiotics to prevent infection-related complications that can also be serious, such as rheumatic fever or acute glomerulonephritis. Group A Streptococcus continues to be a focus of interest not only because of their causal role in acute streptococcal pharyngitis and other pyogenic infections but also because of its association with poststreptococcal sequelae, specifically, acute rheumatic fever and acute glomerulonephritis. In order to properly treat the disease using antibiotic therapy, it is important to use an accurate diagnostic method to identify the pathological agent. For the screening of Group A Streptococcus infection several methods are currently used including Bacitracin sensitivity tests on a sheep blood agar plate, latex agglutination and enzyme immunoassay. STREP A RAPID TEST is a rapid immunological test for the qualitative and visual detection of group A Streptococcus bacterium in throat infections. The test helps to quickly identify whether a sore throat is caused by Strep A rather than by other pathogens (usually viruses) that may not require antibiotic treatment.

#### PRINCIPLE OF THE TEST

STREP A RAPID TEST detects specific Streptococcus A antigens thanks to special gold conjugate antibodies embedded to a test strip. The specimen, which is pharyngeal swab, is taken with a disposable sterile rayon swab. The specimen is placed in a plastic collection tube where the user dispenses a specific volume of 2 extraction reagents (in order to extract specific Streptococcus A antigens). After an incubation time, the obtained sample-solution is dispensed onto the test cassette sample well.

The provided reagents are:

- R1: solution mainly containing Sodium Nitrite and Phenol Red.

- R2: solution mainly containing Acetic Acid.

Once mixed, the two solutions create Nitrous Acid, acting as Streptococcus A specific antigen extractor.

In case of a Strep A concentration equal or higher than 10<sup>4</sup> CFU/mL in the patient's sample, the specific antigens present in the sample solution are firstly bound by the anti-Strep A polyclonal antibodies labelled with gold molecules. The immunocomplexes are selectively captured by the test line (T line, made of purified anti-Strep A polyclonal antibodies), thus showing a red-purple band appearance in correspondence of the T line, pointing out a positive result. In case of absence of Strep A or when its levels in patient's sample are lower than the cut-off value (10<sup>4</sup> CFU/mL), not enough specific immunocomplex is formed. Consequently, only one band, the C line, is visible indicating a negative result. To serve as a procedural control, a colored line will always appear in the Control line region, indicating that the test procedure has been properly performed and that the test components and reagents have operated as intended.

#### REFERENCES

- Anjos LM, Marcondes MB, Lima MF, Mondelli AL, Okoshi MP, Streptococcal acute pharyngitis, Rev Soc Bras Med Trop 2014; 47(4):409-13;
- Castro SA, Dorfmueller HC A brief review on Group A Streptococcus pathogenesis and vaccine development, R Soc Open Sci 2021;
- Centers for Disease Control and Prevention. Adult appropriate antibiotic use summary; <http://www.cdc.gov/getsmart/community/materials-references/print-materials/hcp/adult-approp-summary.pdf>. Accessed March 8, 2016;
- Banerjee S, Ford C. "Rapid Tests for the Diagnosis of Group A Streptococcal Infection: A Review of Diagnostic Test Accuracy, Clinical Utility, Safety, and Cost-Effectiveness". Ottawa (ON): Canadian Agency for Drugs and Technologies in Health; 2018 May. CADTH Rapid Response Reports;
- Llor C, Bjerrum L, Munck A, Cots JM, Hernández S, Moragas A; HAPPY AUDIT Investigators. "Access to point-of-care tests reduces the prescription of antibiotics among antibiotic-requesting subjects with respiratory tract infections." Respir Care. 2014 Dec;
- Georgers N, Nakhoul MD and John Hickner, MD, MSc, Management of Adults with Acute Streptococcal Pharyngitis: Minimal Value for Backup Strep Testing and Overuse of Antibiotics, J Gen Inter Med. 2013 Jun, 28(6):830-834;
- Emily H. Stewart, Brian Davis, B. Lee Clemans-Taylor, Benjamin Littenberg, Carlos A. Estrada, and Robert M. Centor, Rapid Antigen Group A Streptococcus Test to Diagnose Pharyngitis: A Systematic Review and Meta-Analysis, PLoS One. 2014; 9(11);
- Parks T, Barrett L, Jones N. Invasive streptococcal disease: a review for clinicians. British Medical Bulletin, Volume 115, Issue 1, September 2015, Pages 77-89.

### PRECAUTIONS

- Read these instructions for use carefully before performing the test. The test is reliable if the instructions are followed correctly (reaction times, how the sample is collected and the extraction instructions shown).
- Do not use the test after the expiry date or if the package has been damaged.
- Store the test components at a temperature between +4 °C and +30 °C. DO NOT FREEZE.
- In vitro diagnostic device for professional use.
- The test is for external use only. Do not swallow solutions 1 and 2 and the other items in the box! If swallowed, seek medical advice immediately. Show the remaining parts of the box, the instructions for use and the package.
- If liquids come into contact with the eyes, skin or other mucosae, rinse immediately and thoroughly with water and contact a physician. Show the vials' labels.
- After using, dispose of all components according to local laws. Liquid wastes should not be disposed of down the drain.
- Do not open the sealed elements (test case, swab and tongue depressor) until just before performing the test. Take care to open as marked.
- Do not switch the caps of the reagent vials or mix components from different test kits.
- All kit components are single-use, with the exception of the workstation (which can be re-used), the instructions for use and solutions 1 and 2, which must be stored, closed with their caps, at a temperature between +4 °C and +30 °C.
- The operator is required to comply with the local regulations in force and to operate according to good laboratory practices (GLP). Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.

### LIMITATIONS

- The test procedure and the interpretation of the test result shall be carefully followed. Correct sample collection is essential for optimal test performance. Failure to follow the procedure may result in inaccurate results.
- STREP A RAPID TEST is for the professional in vitro diagnostic use only. The test should be used for the determination of Streptococcus A antigen on pharyngeal swabs only. It cannot determine either the quantitative value or the percentage increase in Streptococcus A antigen concentration.
- This test only indicates the presence of Streptococcus A antigen in the specimen due to viable or non-viable Group A Streptococcus bacteria or it does not differentiate asymptomatic carriers of group A Streptococcus from those with a symptomatic infection.
- The accuracy of the test depends on the quality of the swab specimen. False negative results may occur due to improper specimen collection or storage. A negative result may also be obtained from patients at the onset of the disease due to low antigen concentration.
- An excessive amount of blood or mucus on the swab may interfere with the test results and lead to a false positive result. Avoid contact of the swab with the tongue, cheeks and teeth and any oral areas with open wounds at the time of sampling.
- As with all diagnostic tests, results should be interpreted in the light of the patient's overall clinical history.
- If the reading time (5 minutes) is not strictly respected, wrong results may be obtained.
- Respiratory infections, including pharyngitis, can be caused by streptococci of serogroups other than group A as well as other pathogens. A negative test result for group A Streptococcus does not exclude infection with other pathogenic microorganisms. Positive results do not preclude co-infections with other pathogens.
- If clinical signs and symptoms are not consistent with test results, additional testing using other clinical methods (e.g. bacterial culture) is recommended.
- The use of antibiotics or over-the-counter medication may suppress the growth of group A Streptococcus in culture despite the presence of microorganisms detectable by rapid antigen tests. This might lead to discrepant results.

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
## ENGLISH PERFORMANCE

Two swabs (test and reference method) were obtained from each patient and collected by physicians participating in the study. The comparator swab was analysed according to its instruction for use. The results were compared to those of the reference method and reported as positive and negative in Table 1. The performance of the test are reported in terms of Sensitivity, Specificity and Accuracy.

Table 1. Table of results.

STREP A RAPID TEST	REFERENCE METHOD		
	POSITIVE	NEGATIVE	
POSITIVE	34	2	36
NEGATIVE	3	199	202
	37	201	238

SPECIFICITY =  $(199/201) * 100 = 99,00\%$  (CI 95%: 97,63 – 100%)  
 SENSITIVITY =  $(34/37) * 100 = 91,89\%$  (CI 95%: 78,70 – 97,20%)  
 ACCURACY =  $(233/238) * 100 = 97,90\%$  (CI 95%: 96,08 – 99,72%)

## CROSS-REACTIVITY

Each sample was tested as indicated in the table below, according to STREP A RAPID TEST instruction for use in force, at the indicated concentration.

BACTERIA	ATCC	LOT	CONCENTRATION ASSAYED
<i>Bordetella pertussis</i>	9340	BP-02	10 <sup>7</sup> CFU/mL
<i>Campylobacter jejuni</i>	29428	CJ-02	10 <sup>7</sup> CFU/mL
<i>Candida albicans</i>	MUT6298	CA-01	10 <sup>7</sup> CFU/mL
<i>Clostridium perfringens</i>	13124	CP-02	10 <sup>7</sup> CFU/mL
<i>Corynebacterium diphtheriae</i>	11913	CD-02	10 <sup>7</sup> CFU/mL
<i>Enterobacter cloacae</i>	MUT6313	EC-02	10 <sup>7</sup> CFU/mL
<i>Helicobacter pylori</i>	43504	HP-01	10 <sup>7</sup> CFU/mL
<i>Hemophilus influenzae</i>	49766	HI-01	10 <sup>7</sup> CFU/mL
<i>Klebsiella pneumoniae</i>	MUT6314	KP-02	10 <sup>7</sup> CFU/mL
<i>Proteus vulgaris</i>	MUT6312	PV-02	10 <sup>7</sup> CFU/mL
<i>Pseudomonas aeruginosa</i>	MUT6310	PA-02	10 <sup>7</sup> CFU/mL
<i>Shigella flexneri</i>	12022	SF-02	10 <sup>7</sup> CFU/mL
<i>Staphylococcus aureus</i>	MUT6311	SA-02	10 <sup>7</sup> CFU/mL
<i>Streptococcus agalactiae</i>	12386	SAG-02	10 <sup>7</sup> CFU/mL
<i>Streptococcus dysgalactiae</i>	12388	SD-02	10 <sup>7</sup> CFU/mL
<i>Streptococcus gallolyticus</i>	9808	SG-02	10 <sup>7</sup> CFU/mL
<i>Yersinia enterocolitica</i>	51871	YE-02	10 <sup>7</sup> CFU/mL

None of the bacteria tested in this study, at the concentration given in the table, have shown cross reactions or interferences with STREP A RAPID TEST.

## MEASURING RANGE

The limit of detection study was conducted on Certified Positive ATCC *Streptococcus pyogenes*. The LOD of STREP A RAPID TEST is 10<sup>4</sup> CFU/mL. No adverse effect on T-line formation (prozone effect) was observed for group A *Streptococcus* concentrations up to 5 x 10<sup>8</sup> CFU/mL.

## CONTENT

THE PACKAGE CONTAINS:

- 20\* hermetically sealed aluminum pouches each containing: 1 STREP A TEST cassette and 1 desiccant bag;
- 20\* sterile swabs;
- 20\* sterile wooden tongue depressors;\*\*
- 1 workstation;
- 20\* empty plastic extraction tubes with dropper
- 1\* dropper vial containing enough STREP A TEST-R1 solution to perform 20\* tests

R1 contains 20.7% Sodium Nitrite  
 H302: Harmful if swallowed  
 P264: Wash carefully your hands after use  
 P270: Do not eat, drink or smoke while using this product.

- 1\* dropper vial containing enough STREP A TEST- R2 solution to perform 20\* tests
- R2 contains 0.55% Acetic Acid
- 1\* Instructions for use leaflet

\* the number of devices and accessories inside the kit may vary depending on the commercial reference purchased, for the exact number please consult the "contents" section on the outside of the package.

\*\*optional accessory only present in some STREP A RAPID TEST formats.

- Required material not supplied: a timer (e.g. timer, watch).  
 - Do not open the sealed elements (test case, swab and tongue depressor) until just before performing the test. Take care to open as marked.

## SPECIMEN COLLECTION AND PREPARATION

Bring tests, specimens, reagents to room temperature (15-30°C) prior to testing. To avoid cross contamination, do not allow the tips of the reagent bottles to come into contact with specimen material.

- 1) Place the workstation on a clean, dry, flat surface. Place the empty plastic test tube in one of the compartments of the workstation. – FIG. A
- 2) Open the package containing the sterile swab. Avoid touching the cotton tip and remove the swab using the plastic handle. Tilt the patient's head backwards and let open the mouth as much as possible. Pick up the tongue depressor with one hand and flatten the tongue. Use the other hand to place the swab close to the throat. Touch the back of the throat - the area around the tonsils and any reddened or painful part (soft palate, uvula) with the cotton tip. Rotating the swab is recommended because it increases the amount of collected sample. – FIG. B

It is recommended that swab specimens be processed as soon as possible after collection. If swabs are not to be processed immediately, they should be placed in a sterile, dry, tightly capped tube or bottle and refrigerated. Do not freeze swabs. Swabs can be stored at room temperature (15-30°C) for up to 8 hours or refrigerated (2-8°C) for up to 72 hours. All specimens should be brought to room temperature (15-30°C) prior to testing.

## TEST PROCEDURE

- 1) After collecting the sample, insert the swab's cotton tip into the supplied plastic test tube that you had previously placed in the workstation. – FIG. C
- 2) Open the vial STREP A RAPID TEST - R1 by unscrewing only the blue cap (leave the white part firmly sealed), then add 4 drops into the test tube and close it with the cap. – FIG. D
- 3) Next, open the vial STREP A RAPID TEST - R2 by unscrewing only the green cap (leave the white part firmly sealed), then add 4 drops into the test tube and close it with the cap. – FIG. E  
**WARNING:** the addition of R2 to R1 leads to the formation of a solution for which the indications described in point 8 and 9 of the "precautions" apply.
- 4) Holding the swab's plastic handle, rotate the swab against the sides of the test tube about 10 times to thoroughly mix the solution. Leave the swab to incubate for 2-5 minutes. – FIG. F
- 5) At the end of the incubation period, remove the test tube (with the swab and extraction fluid still in it) from the workstation. Using your thumb and index finger, press the sides of the test tube to release as much fluid as possible from the swab's cotton tip and collect it in the test tube. Remove the swab – FIG. G. Dispose of the swab in compliance with guidelines for the handling of infectious agents and put the test tube back into one of the compartments in the workstation.
- 6) Add the supplied dropper to the plastic test tube. – FIG. H
- 7) Dispense 1 drop into the well in the test case. – FIG. I **NOTE:** If the dispensed drop contains air bubbles, add another drop to the well.
- 8) Read the results after 5 minutes (some positive results can be seen after just 1 minute). Do not read after 10 minutes.

## RESULTS INTERPRETATION

READ THE RESULTS AFTER 5 MINUTES. DO NOT READ THE RESULTS AFTER 10 MINUTES  
 The intensity of the line colors is not relevant for the purposes of interpretation of the Test's results.

### POSITIVE

Two colored lines appear. A colored line should always appear in the Control (C) region and another line should be in the Test (T) region.

This result means that the test has detected the presence of *Streptococcus A* specific antigens in the sample.

### NEGATIVE

One colored line appears in the Control region (C). No apparent colored line appears in the Test line region (T).

This result is negative and means that no *Streptococcus A* antigens have been detected in the sample from the pharynx or that their amount is too low to be detected.

### INVALID

No bands appear or there is a line only under the T (Test) sign and not under the C (Control) sign.

In this case it is not possible to interpret the result of the test, which must be considered not valid. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor. Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

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