Your partner for rapid testing

PRIMA Lab SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
support@primabsa.ch
primabsa.ch

REF 800087
800087\_IFU\_11\_1\_0 12/2023

Registered trademarks

KI SCREENING TEST

The first rapid test for semi-quantitative detection of NGAL in urine for early screening of kidney damage in dogs. For professional use only.

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SYMBOLS / SYMBOLE / SIMBOLES / SIMBOLOGIA

Read the instructions for use carefully... Sufficient for <math>\lt; n \gt;</math> tests... Do not reuse... Temperature limits... Manufacturer

ENGLISH

NGAL PROTEIN

Neutrophil Gelatinase Associated Lipocalin (NGAL), a protein synthesised by many tissues, is freely filtered by the glomeruli and reabsorbed by renal tubules.

PRINCIPLE AND USE OF THE TEST

KI SCREENING TEST is a lateral flow immunochromatographic device, for professional use, for semi-quantitative detection of NGAL in dog urine samples.

PRECAUTIONS AND LIMITATIONS

- 1. Read these instructions carefully before testing. The test is reliable if the instructions are carefully followed.
2. This kit is for canine veterinary use only, do not use it on other animals.
3. Do not use the test after the expiry date or if the packaging is damaged.
4. Do not mix kit components from different batches.
5. Store kit components at a temperature of +4 °C to +30 °C. Do not freeze.
6. Do not open the foil package unless it is to be used.
7. After use, dispose of all components in compliance with local regulations.
8. All components of the kit are to be considered single-use, with the exception of the instructions for use.
9. The presence of a high number of leukocytes in the urine could cause an increase in urinary NGAL, since the biomarker is produced by neutrophil granulocytes.
10. The presence of significant leukocytosis could cause an increase in urinary NGAL.
11. Do not eat, drink or smoke in the area where samples or kits are handled.
12. The test is for external use only.

PERFORMANCE STUDIES

ANALYTICAL PERFORMANCE

The test was evaluated on 180 urine samples from different clinical settings of urinary tract disease by comparing it with the gold standard (ELISA, BioPorto Diagnostics inc).

SENSITIVITY: 99.21% (CI 95%: 95.67-99.86%)

SPECIFICITY: 98.15% (CI 95%: 90.23-99.67%)

ACCURACY: 98.90% (CI 95%: 96.06-99.70%)

CLINICAL PERFORMANCE

When used in the setting of a clinical trial on dogs, the device showed that a urinary NGAL concentration of 20 ng/mL was the cut-off capable of distinguishing the presence of AKI, with a sensitivity of 97.6%.

INTERFERENCE

No interference was observed in the presence of the following substances: ethanol, ultrasound gel, haemoglobin, bilirubin, glucose, ketones, leukocytes, nitrites, total protein, urobilinogen, and under different urinary pH conditions, even non-physiological (pH = 3 and pH = 9).

CONTENTS

\*for the exact number of components listed, please see the label on the outside of the kit.

- \*rapid tests
\*single-use pre-dosed bottles
\*pipettes
\*colour chart
1 instructions for use

Materials required but not provided: timer, urine sample container.


PREPARATION AND COLLECTION OF THE SAMPLE

- 1) Collect the patient's urine in a clean, dry container using the collection method deemed most appropriate.
2) Centrifuge the patient's urine (5000 rfc for 5 minutes). Centrifugation of the sample is recommended but not mandatory.

If not used immediately, storage at +2 °C/+8 °C is recommended. If not used within 24 hours at this temperature, freezing at -20 °C is recommended.

TEST PROCEDURE

If stored cool, allow the kit components and/or sample to reach room temperature (15-30 °C) before testing.

- 1) Open the thinner bottle by unscrewing the white portion of the cap and prepare the necessary material as follows: open the aluminium bag, take out only the cassette and place it on a flat, dry surface. Discard the desiccant bag.
2) Transfer 2 drops of sample into the diluent bottle.
3) Close the diluent bottle with the cap and mix by gently turning the bottle upside down 3-5 times.
4) Open the dropper of the bottle by unscrewing the blue portion of the cap and dispense 3 drops of solution into the well (S) of the test cassette.
5) Wait 10 minutes and place the cassette on the dotted section of the colour chart provided.

INTERPRETATION OF THE RESULTS

Results must be read after 10 minutes and no later than 12 minutes.

The intensity of the test line (T) is directly proportional to the NGAL concentration present in the sample analysed.

The intensity of the control line (C), as long as it is visible, is of no importance for the reading of the result.

RESULT ASSOCIATED CLINICAL SETTING

Table with 2 columns: RESULT and ASSOCIATED CLINICAL SETTING. Rows include NGAL ~90 ng/mL or greater, NGAL ~20 ng/mL, NGAL ~4 ng/mL, and NGAL <4 ng/mL.

Note: The limit of detection (LoD) of the test is 2.4 ng/mL. Therefore, samples with a concentration below 4 ng/mL but equal to or above the LoD may give a positive result with a test line (T) of less intensity than that shown on the colour chart yet still visible.

INVALID

The test is considered invalid if the control line (C) does not appear. The most common reasons for the absence of the control line (invalid result) are insufficient sample volume or performing the procedure incorrectly.

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BIBLIOGRAPHY / BIBLIOGRAFIE / BIBLIOGRAPHIE / BIBLIOGRAFIA

1. Davis J, Rasis AL, Miller DW, Hosgood GL, Rossi G. Analytical validation and reference intervals for a commercial multiplex assay to measure five novel biomarkers for acute kidney injury in canine urine. Res Vet Sci. 2021 Oct;139:78-86.
2. Lee YJ, Hu YY, Lin YS, Chang CT, Lin FY, Wong ML, Kuo-Hsuan H, Hsu WL. Urine neutrophil gelatinase-associated lipocalin (NGAL) as a biomarker for acute canine kidney injury. BMC Vet Res. 2012 Dec 28;8:248.
3. Nabby M, Hokamp J. Urinary Biomarkers of Kidney Disease in Dogs and Cats. Vet Clin North Am Small Anim Pract. 2023 Jan;53(1):53-71. doi: 10.1016/j.cvsm.2022.07.006. Epub 2022 Oct 19. PMID: 36270837.
4. Steinbach S, Weis J, Schweighäuser A, Francey T, Neiger R. Plasma and urine neutrophil gelatinase-associated lipocalin (NGAL) in dogs with acute kidney injury or chronic kidney disease. J Vet Intern Med. 2014 Mar-Apr;28(2):264-9. doi: 10.1111/jvim.12282.

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