Your partner for rapid testing

PRIMA Lab SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
support@primalabs.ch
primallabs.ch

Qarad EC-REP BV
Pas 257, 2440 Geel - Belgium

800034IP-Z0P 800034IP-Z0P_IFU_90_50 03/2024

Table 1. Performance / Tabelle 1. Leistung / Tableau 1. Performances / Tabella 1. Prestazioni / Tabla 1. Rendimiento / Tabela 1. Desempenho / Tabela 1. Działanie

Table with 5 columns: RESULTS, Negative, Positive, >80, TOTAL. Rows include Negative, Positive, and TOTAL results for C-REACTIVE PROTEIN RAPID TEST.

Relative Specificity = 96.0% (CI 95%:86.5%–98.9%)
Relative Sensitivity = 98.7% (CI 95%: 93.0%–99.8%)
Relative Accuracy = 97.6 % (CI*: 93.3%–99.2%)

Table 2. Interfering substances / Tabelle 2. Beeinträchtigende Substanzen / Tableau 2. Substances interférentes / Tabella 2. Sostanze interferenti / Tabla 2. Sustancias que interfieren / Tabela 2. Substâncias interferentes / Tabela 2. Substancje zakłócające

Table with 2 columns: SUBSTANCE, CONCENTRATION ASSAYED. Lists various substances like Bilirubin, Cholesterol, Triglycerides, etc.

SYMBOLS / SYMBOLE / SIMBOLOS / SIMBOLOGIA / SIMBOLS / SIMBOLGIA / SYMBOLE

IVD, I, CE, REF, LOT, EC REP icons and their corresponding text descriptions in multiple languages.

REFERENCES / BIBLIOGRAFIE / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA

- 1. C-reactive protein concentrations as a marker of inflammation or infection for interpreting biomarkers of micronutrient status.
2. Calarco et al. 2023. Analytical Performance of 17 Commercially Available Point-of-Care Tests for CRP to Support Patient Management at Lower Levels of the Health System.
3. Dittrich et al. (2016). Target Product Profile for a Diagnostic Assay to Differentiate between Bacterial and Non-Bacterial Infections and Reduce Antimicrobial Overuse in Resource-Limited Settings.
4. Martínez-González et al. 2020. Point-of-Care C-Reactive Protein Testing to Reduce Antibiotic Prescribing for Respiratory Tract Infections in Primary Care.

C-REACTIVE PROTEIN RAPID TEST

Rapid test for the semi-quantitative detection of C-Reactive Protein (CRP) in whole blood or serum samples. Schnelltest für den semi-quantitativen Nachweis von C-reaktivem Protein (CRP) in Vollblut- oder Serumproben.

LIMITATIONS

- 1. As with any diagnostic method, the physician must evaluate the test result in light of other clinical evidence.
2. People suffering from haemophilia or taking anticoagulants could have an excessive bleeding after the fingerprick.
3. Despite the reliability of the test, false positive or negative results can not be excluded.
4. For CRP values above 750 mg/L, the C-REACTIVE PROTEIN TEST may give unreliable results.
5. The test results should only be interpreted after 5 minutes. If the reading time is not strictly respected, wrong results may be obtained.

CONTENT

- 20 hermetically sealed aluminum bags containing 1 C-REACTIVE PROTEIN TEST cassette and 1 desiccant bag;
• 20 vials with dropper tip containing enough C-REACTIVE PROTEIN TEST diluent for 20 tests;
• 20 capillary pipette (10 µL each) for collecting the blood sample;
• 1 instructions for use leaflet.

Materials required but not provided: a time instrument (i.e. stopwatch, clock), centrifuge (for serum/plasma), sample collection containers.

PERFORMANCE CHARACTERISTICS

To assess the performance of C-REACTIVE PROTEIN RAPID TEST, a clinical evaluation was performed by comparing the results obtained against the laboratory's reference method. 127 whole blood samples were tested with both methods.

LIMIT OF DETECTION

The study concerning the Limit Of Detection (LOD) was conducted with the C-Reactive Protein (CRP) calibrators Sentinel. The LOD of C-REACTIVE PROTEIN RAPID TEST is 10 mg/L.

INTERFERING SUBSTANCES
The following compounds were tested with C-REACTIVE PROTEIN RAPID TEST and no interference was observed (Table 2).

TEST PROCEDURE

Allow the kit components and/or samples to reach room temperature (15-30 °C) before testing. Prepare the necessary material as follows: open the aluminum bag, remove the cassette and discard the desiccant bag.

WHOLE BLOOD SAMPLING PROCEDURE

- 1) Prick the finger with the lancet and obtain a drop of blood. Fill the capillary pipette to the black line - fig. A1.
2) Transfer 10 µl of whole blood into the vial previously opened. Press the bulb 2 or 3 times to ensure that all the blood has moved into the diluent solution.
3) Replace the screw cap onto the vial and mix well for at least 10 seconds.
4) Start the timer. Read the results at exactly 5 minutes.

PRECAUTIONS

- 1. Read these instructions for use carefully before performing the test.
2. Do not open the sealed aluminum pouch until just before performing the test.
3. Humidity and temperature can adversely affect the results.
4. Handle all specimens as if they contain infectious agents.
5. The test is for external and single use only.
6. Do not use the kit after the expiry date or if the packaging and/or internal components are damaged.
7. Store the buffer correctly closed with its cap.

RESULTS INTERPRETATION

READ THE RESULTS AT EXACTLY 5 MINUTES. The intensity of the line colors is not relevant for the purposes of interpretation of the test's results.

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POSITIVE, C-REACTIVE PROTEIN BETWEEN 10 AND 40 mg/L

Two coloured lines appear. One coloured line should always appear in the Control (C) line region and another line should appear in the Test line region 1 (T1).

POSITIVE, C-REACTIVE PROTEIN BETWEEN 40 AND 80 mg/L

Three coloured lines appear in the reading window, one at the C sign and two more at the T1 and T2 signs, as shown in the figure opposite.

POSITIVE, C-REACTIVE PROTEIN OVER 80 mg/L

Four coloured lines appear in the reading window, one at sign C plus three more at signs T1, T2 and T3 as shown in the figure opposite.

NEGATIVE, C-REACTIVE PROTEIN BETWEEN BELOW 10 mg/L

A coloured line appears in the Control (C) line region. No line appears in the Test (T1, T2, T3) regions.

INVALID RESULT

The Control (C) line does not appear. An insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the Control line (invalid result).

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