TEST PROCEDURE

- Wash hands with soap and warm water, rinse with clean water and allow to dry. Note: The use of warm water facilitates capillary blood collection as it induces vasodilation.
- Layout the material needed as follows: open the aluminium bag, take out only the test cassette and dispose of the desiccant packet. Open the plastic packet containing the pipette. Open the vial by unscrewing the white cap, making sure it does not fall to the ground. **FIG. A**
- Carefully rotate the protective cap of the sterile lancet 360° without pulling it, then extract and discard the cap. **FIG. B**
- Carefully massage the finger chosen for the puncture (the side of the ring finger is recommended). It is important that the massage is done from the palm of the hand to the fingertip, to improve blood flow. Press the open end of the lancet (the side the cap has been extracted from), against the fingertip. **FIG. C**. The tip of the lancet automatically retracts after use. If the lancet does not work properly, discard it and use the second one supplied. If the second one is not required, it can be disposed of without special precautions.
- Holding the hand downwards, massage the finger until a large drop of blood forms. It is important to massage from the palm of the hand to the fingertip to improve blood flow. **FIG. D**
- Take the pipette without pressing the bulb. Two sampling methods are suggested:
 - FIG. E1:** hold the pipette horizontally without pressing the bulb place it in contact with the drop of blood, it will enter the pipette by capillarity. Move the pipette away when the black line is reached. If there is not enough blood, continue to massage the finger until the black line is reached.
 - FIG. E2:** place the pipette on a clean, flat surface with the tip protruding from the surface, then place the drop of blood in contact with the pipette, it will enter by capillarity. If blood is not sufficient, continue massaging the finger until the blood has reached the black line. Bring the pipette in contact with the finger, it will enter by capillarity. Avoid, as far as possible, to continuously move the tip of the pipette away from the finger in order to prevent the formation of air bubbles.
- Deposit the blood collected with the pipette into the previously opened vial. Press the pipette bulb 2 or 3 times to ensure that all the blood has moved into the diluent solution. Replace the screw cap onto the vial and mix well for at least 10 seconds. **FIG. F**
- Unscrew the blue cap from the dropper vial (leave the white cap tightly screwed on). **FIG. G**
- Deposit 3 drops into the well shown on the cassette (S). **FIG. H**
- Wait 5 minutes and then read the results as indicated.

RESULTS INTERPRETATION

READ THE RESULT AT EXACTLY 5 MINUTES.
DO NOT READ THE RESULT AFTER 6 MINUTES

5 min.

The result must be interpreted independently of the colour intensity of the bands that may appear in the reading window.

NEGATIVE, C-REACTIVE PROTEIN BELOW 10 MG/L

A coloured band appears only under the C (Control) sign, as shown in the figure on the right. This result indicates that the concentration of CRP is below 10 mg/L and suggests the absence of infection and inflammation.

POSITIVE, C-REACTIVE PROTEIN BETWEEN 10 AND 40 MG/L

Two coloured bands appear in the reading window, one by the C (Control) sign, the other by the T1 sign, as shown in the figure on the right. It is recommended to consult a physician, as CRP values between 10 and 40 mg/L may be associated with a viral infection or the onset of a bacterial infection. An inflammation caused (for example) by physical trauma can also lead to this kind of result.

POSITIVE, C-REACTIVE PROTEIN BETWEEN 40 AND 80 MG/L

In the reading window, three coloured bands appear, one by the C (Control) sign and two others by the T1 and T2 signs, as shown in the figure on the right. Consult a physician, as this result may be associated with a viral or bacterial infection or physical trauma.

POSITIVE, C-REACTIVE PROTEIN OVER 80 MG/L

Four coloured bands appear in the reading window, one by the C (Control) sign and three others by the T1, T2 and T3 signs, as shown in the figure on the right. Consult a doctor, as this result may be associated with a bacterial infection or severe inflammation.

NOTE: The colour intensity in the control and test line regions may vary (as represented in the images). The intensity of the bands is not relevant; only their presence must be considered.

NON VALID RESULT

The control band does not appear (no lines appear or one or more T bands appear without the C band). Insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the control line. Review the procedure and repeat the test with a new device and a new blood sample.

F.A.Q. - QUESTIONS AND ANSWERS

HOW DOES THE C-REACTIVE PROTEIN TEST WORK? The test detects the presence of CRP in human whole blood samples using specific antibodies and colloidal gold particles embedded in the test strip.

WHEN CAN THE TEST BE USED? The C-REACTIVE PROTEIN TEST may be used in the case of infection or inflammation symptoms such as fever, headache or weakness, or to monitor the health process following surgery and minor injuries. The test can also be used to monitor the effectiveness of a therapeutic treatment and can be performed at any time of the day.

CAN THE RESULT BE INCORRECT? The result is correct as long as the instructions are carefully followed and the test is stored under the conditions listed in the "Precautions" section. However, the result may not be correct if: the device comes into contact with other liquids before being used, if the amount of blood and/or diluent is not enough, if the number of drops dispensed in the well is incorrect or if the reading times are not strictly followed. The supplied plastic pipette allows users to be sure they have collected the right amount of blood.

HOW ACCURATE IS THE C-REACTIVE PROTEIN TEST? The C-REACTIVE PROTEIN TEST showed high concordance with analytical reference methods, with 97% accuracy in distinguishing between normal (< 10 mg/L) and abnormal values.

HOW IS THE TEST TO INTERPRET THE TEST IF THE COLOUR AND INTENSITY OF THE TEST AND CONTROL BANDS ARE DIFFERENT? The colour and the intensity of the lines are not important for the interpretation of the result. Refer

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- Hände mit Seife und warmem Wasser waschen, mit frischem Wasser abspülen und trocken lassen. Falls dies nicht möglich ist, verwenden Sie als Alternative die mitgeleistete Gaze. Hinweis: Die Verwendung von warmem Wasser erleichtert die Entnahmen von Kapillarblut.
- Legt das Material bereit: öffnen den Aluminiumbeutel, entnehmen nur die Testkassette und entsorgen den Desiccatorenbeutel. Öffnen Sie die Plastikverpackung der Pipette. Entfernen Sie die weiße Kappe, während Sie den Tropfer drehen. Drehen der weißen Kappe, achten Sie darauf, ihn nicht fallen zu lassen. **ABB. A**
- Carefully rotate the protective cap of the sterile lancet 360° without pulling it, then extract and discard the cap. **FIG. B**
- Carefully massage the finger chosen for the puncture (the side of the ring finger is recommended). It is important that the massage is done from the palm of the hand to the fingertip, to improve blood flow. Press the open end of the lancet (the side the cap has been extracted from), against the fingertip, to improve blood flow. **FIG. C**. The tip of the lancet automatically retracts after use. If the lancet does not work properly, discard it and use the second one supplied. If the second one is not required, it can be disposed of without special precautions.
- Holding the hand downwards, massage the finger until a large drop of blood forms. It is important to massage from the palm of the hand to the fingertip to improve blood flow. **FIG. D**
- Take the pipette without pressing the bulb. Two sampling methods are suggested:
 - FIG. E1:** hold the pipette horizontally without pressing the bulb place it in contact with the drop of blood, it will enter the pipette by capillarity. Move the pipette away when the black line is reached. If there is not enough blood, continue to massage the finger until the black line is reached.
 - FIG. E2:** place the pipette on a clean, flat surface with the tip protruding from the surface, then place the drop of blood in contact with the pipette, it will enter by capillarity. If blood is not sufficient, continue massaging the finger until the blood has reached the black line. Bring the pipette in contact with the finger, it will enter by capillarity. Avoid, as far as possible, to continuously move the tip of the pipette away from the finger in order to prevent the formation of air bubbles.
- Deposit the blood collected with the pipette into the previously opened vial. Press the pipette bulb 2 or 3 times to ensure that all the blood has moved into the diluent solution. Replace the screw cap onto the vial and mix well for at least 10 seconds. **FIG. F**
- Unscrew the blue cap from the dropper vial (leave the white cap tightly screwed on). **FIG. G**
- Deposit 3 drops into the well shown on the cassette (S). **FIG. H**
- Wait 5 minutes and then read the results as indicated.

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IS THE RESULT RELIABLE IF READ AFTER 6 MINUTES? No. The test must be read 5 minutes after the procedure is completed and no later than 6 minutes from this time. Results read after 6 minutes may not be correct.

WHAT SHOULD I DO IF THE RESULT IS POSITIVE? If the result is positive and the CRP is higher than 10 mg/L, refer to the "Result interpretation" section and contact your doctor as soon as possible, because antibiotic treatment may be required.

WHAT SHOULD I DO IF THE RESULT IS NEGATIVE? A negative result means that the concentration of CRP is below 10 mg/L, which means that there is neither viral nor bacterial infection in progress. If symptoms persist, seek medical advice.

HOW ACCURATE IS THE C-REACTIVE PROTEIN TEST? The C-REACTIVE PROTEIN TEST showed high