Test your health

PRIMA LAB SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND

support@primalabs.ch
primalabs.ch

EC REP Qarad EC-REP BV
Pas 257
2440 Geel - Belgium

800040_IFU_11_5.2 112023

REF 800040

Registered trademarks.

INSTRUCTIONS FOR USE
GEBRAUCHSANWEISUNG
NOTICE D'UTILISATION
ISTRUZIONI PER L'USO

0483

PRIMA
HOME TEST

BOWEL FOB TEST

Self-test for the detection of occult blood in the faeces
Selbsttest zur Bestimmung von okkultem Blut im Stuhl
Auto-test pour la détermination du sang occulte dans les selles
Self-test per la determinazione del sangue occulto nelle feci

ENGLISH

FOB (FECAL OCCULT BLOOD)

Faecal occult blood is the presence in the faeces (stool) of blood that is not clearly visible to the naked eye. Occult blood may be the consequence of gastro-intestinal disorders such as ulcers, polyps, colitis, diverticulitis, rhagades, haemorrhoids or cancer. The signs of these lesions are often silent in the first phases, so searching for occult blood in the stool is an important screening test, which allows to identify the presence of pathologies in the gastro-intestinal system at an early stage.

PRINCIPLES OF THE TEST

BOWEL FOB TEST is an immunochromatographic test that detects the smallest amounts of blood in stool samples by using highly specific monoclonal antibodies.

These properties mean that the BOWEL FOB TEST does not require particular dietary restrictions before being carried out.

REFERENCES

1. Van Rossum Leo G. et al.: "Random Comparison of Guaiac and Immunochemical Faecal Occult Blood Tests for Colorectal Cancer in a Screening Population" Gastroenterology, Volume 135, Issue 1, 82 – 90.
2. Faivre J., Dancourt V., Denis B., Dorval E., Piette C., Perrin P., Bidan J.M., Jard C., Jung S., Levillain R., Viguier J., Bretagne J.F.: "Comparison between a guaiac and three immunochemical faecal occult blood tests in screening for colorectal cancer". European Journal of Cancer. 2012; 48(16):2969-76.
3. Benton S.C., Seaman H.E., Halloran S.P.: "Faecal Occult Blood Testing for Colorectal Cancer Screening: the past or the Future". Current Gastroenterology Reports. 2015; 17(2):428
4. European Commission. European guidelines for quality assurance in colorectal cancer screening and diagnosis. First edition 2010..

CONTENT

* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the external box.

- * hermetically sealed aluminium pouch containing:
1 BOWEL FOB TEST cassette
1 desiccant bag
- * vial with red screw cap, with collection stick and dropper, containing enough BOWEL FOB TEST diluent for 1 test
- * small sheet for the collection of the faeces sample
- * 1 instructions for use leaflet

* Material required but not supplied: cotton wool, a device to measure time (i.e. timer, watch).
* Do not open the sealed aluminium pouch until just before performing the test. Take care to open it as marked.
* The desiccant bag must not be used. Dispose of it with household waste without opening it.

PRECAUTIONS

1. Read these instructions for use carefully before performing the test. The Test is reliable only if all the instructions are followed correctly.
2. Keep the Test out of the reach of children.
3. Do not use the Test after the expiry date or if the package has been damaged.
4. Store the Test components at a temperature between +4°C and +30°C. Do not freeze.
5. Use the test once only.
6. The test is for external use only. DO NOT SWALLOW.
7. In vitro diagnostic device for individual use.
8. After using, please dispose of all components according to your local waste disposal laws. Ask your pharmacist for advice.

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


ENGLISH

TEST PROCEDURE

Before performing the test, collect a faeces sample following these instructions:

COLLECTING THE SAMPLE

- 1) Check the toilet seat is dry - contact with water or detergents may prevent the faeces-collecting sheet from adhering properly and interfere with the test result.
- 2) Use the supplied special stool-collecting sheet as described below:
 - a slowly detach the two adhesive ends of the faeces-collecting sheet in the direction pointed by the arrows printed on the sheet. Open the sheet completely;
 - b attach the two adhesive surfaces of the faeces-collecting sheet near the back of the toilet seat, making sure the sheet does not get wet; **–FIG. A**
 - c deposit the sample, carefully avoiding wetting the sheet with urine, and proceed with collection as described in the next steps in this procedure.
- 3) Unscrew the vial's red top and extract the faeces-collecting stick.
- 4) Dip the stick into the faeces in 3 different points of the sample. **–FIG. B**
- 5) Screw on the cap again, immersing the stick in the faeces extraction liquid, and shake the vial for at least 10 seconds. **–FIG. C**
- 6) Detach the adhesive parts of the faeces-collecting sheet at the same time from the toilet seat and let everything fall into the toilet. To make flushing easier, before flushing, wait a few seconds for the paper to soften. Flush

EXECUTION OF THE TEST

- 7) Open the aluminium bag, remove the strip only and place it on a horizontal surface. Discard the desiccant packet.
- 8) Break the end of the red cap, freeing the dropper portion. **–FIG. D**
- 9) Dispense 3 drops of the diluted faeces sample into the well shown on the strip. **–FIG. E**
Wait 1-2 seconds between dispensing one drop and the next. Do not move the device during the test performing to avoid invalid results.
- 10) Wait 5 minutes.

RESULTS INTERPRETATION

READ THE RESULT AFTER EXACTLY 5 MINUTES

The intensity of the line colours is not relevant for the purposes of interpretation of the Test's results.

POSITIVE RESULT

Two coloured bands appear in the reading window by the T (test) and C (control) signs. The T band may be less intense (lighter) than the C line. **This result means there is blood in the faeces and you should consult a physician.** You can obtain positive results even in the presence of haemoglobin concentration in the 36-40 ng/mL range.

NEGATIVE RESULTS

A coloured band appears only under the C (control) sign. **This result means that there is no blood in the faeces or that the concentration is below the test's detection limit (40 ng/mL).**

NON VALID RESULT

If no bands appear or if there is only a line under the T (Test) sign, it means that it was not possible to read the result of the test, which must be considered as not valid. In this case it is not possible to interpret the result of the test, which must be considered not valid. It is recommended to repeat the test with a new BOWEL FOB TEST and new samples of faeces.

F.A.Q. – QUESTIONS AND ANSWERS

HOW DOES THE BOWEL FOB TEST WORK? The BOWEL FOB TEST was made to selectively identify the presence of human haemoglobin in a stool sample using a combination of specific antibodies. The lowest concentration of human haemoglobin the test can detect is 40 ng/mL, although in some cases samples with haemoglobin concentration in the 36-40 ng/mL range may yield a positive result. **WHEN CAN THE TEST BE USED?** The BOWEL FOB TEST can be carried out at any time of the day, when you have a bowel movement, and there is no need to follow a specific diet. It should not be used while menstruating, as the result obtained may not be reliable. **CAN THE RESULT BE INCORRECT?** The result is correct if the instructions are followed carefully. However, the result may not be correct if: the device comes into contact with other liquids before being used, if part of the diluent is accidentally spilt, if the number of drops dispensed in the well is higher than 4 or less than 2. **HOW DO I INTERPRET THE TEST IF THE COLOUR AND INTENSITY OF THE TEST AND CONTROL BANDS ARE DIFFERENT?** The colour and the intensity of the bands are not important for the interpretation of the result. The bands must be

full and homogeneous. The test is positive regardless of the intensity of the colour of the test band. **IS THE RESULT RELIABLE IF READ AFTER 10 MINUTES?** No. The test must be read 5 minutes after the procedure is completed (point 10). **WHAT SHOULD I DO IF THE RESULT IS POSITIVE?** If the result is positive, blood has been detected in the faeces and you must consult a physician, who will establish what to do next. **WHAT SHOULD I DO IF THE RESULT IS NEGATIVE?** A negative result means the test has not detected any blood in the stool (some kinds of lesions bleed intermittently) or that the concentration of blood is below the test's detection limit. Therefore negative results cannot guarantee there is no abnormal condition. Only a doctor can make a diagnosis. **HOW ACCURATE IS THE BOWEL FOB TEST?** The Test is very accurate. Assessment reports show a concordance rate above 92.6% (95% CI: 88.4 - 96.8%) with reference methods. Despite the test's reliability, there may be false positive and false negative results.

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