ENGLISH

CONTENT

* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the external box.

- * hermetically sealed aluminium pouch containing: 1 PROSTATE PSA TEST cassette 1 desiccant bag
* sterile lancets for self-drawing blood
* transparent plastic bag containing a pipette for collecting blood
* alcohol cleansing gauze
* 1 instructions for use leaflet

* Material required but not supplied: cotton wool, a device to measure time (i.e. timer, watch).
* Do not open the sealed aluminium bag until just before performing the test. Take care to open it as marked.
* The desiccant packet must not be used. Dispose of it with household waste without opening it.

TEST PROCEDURE

- 1) Wash hands with soap and warm water, rinse with clean water and allow to dry. If this is not possible, use the provided gauze as an alternative.
2) Prepare the necessary material as follows: open the aluminium pouch, take out only the test cassette and throw away the desiccant bag. Open the plastic packet containing the pipette.
3) Carefully rotate the protective cap of the sterile lancet 360° without pulling it. Extract and discard the released cap.
4) Carefully massage the finger chosen for the puncture (the side of the ring finger is recommended).
5) Holding the hand down, massage the finger until a large drop of blood forms.
6) Take the pipette without pressing the bulb. Two sampling methods are suggested:
7) Place the blood collected with the pipette into the well indicated on the cassette (S) by pressing the pipette bulb.
8) Unscrew the blue cap from the dropper vial (leave the white cap tightly screwed on). Deposit 2 drops into the well indicated on the cassette (S), wait for 5 seconds between the first and second drop.
9) Wait 5 minutes and read the results referring to the next section for result interpretation.

RESULTS INTERPRETATION

READ THE RESULT AFTER EXACTLY 5 MINUTES
The intensity of the line colours is not relevant for the purposes of interpretation of the Test's results.

POSITIVE RESULT

Two coloured bands appear in the reading window by the T (test) and C (control) signs. The T band may be less intense (lighter) than the C line. This result means that the value of Prostate-Specific Antigen (PSA) in the blood is higher than normal, so you must consult a physician.

NEGATIVE RESULT

A coloured band appears only under the C (control) sign. This means that the value of Prostate-Specific Antigen (PSA) in the blood is normal. A negative result does not exclude the possibility of prostatic pathologies.

NON VALID RESULT

No bands appear or there is a line only under the T (test) sign and not under the C (control) sign. In this case it is not possible to interpret the result of the test, which must be considered not valid. Repeat the test with a new blood sample.

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