ENGLISH

## OVULATION LH TEST

*Self-test for detecting the luteinizing hormone (LH) in urine*

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During a woman's fertile age (from puberty to menopause), her reproductive system will go through a

series of important structural and functional changes every month. This occurs as a result of the continuous variations in the levels of oestrogens and progesterone. These hormones are regulated by the hypothalamic-

pituitary-gonadal axis and are responsible for the menstrual cycle, which has 4 phases, i.e. menstrual phase, follicular phase, ovulation and luteal phase. During the ovulation phase, which occurs around day 12/13 of

the menstrual cycle, the level of oestrogens in the blood is so high that the pituitary gland starts to produce the luteinizing hormone (LH). This leads to follicle rupture and the release of the egg, which is ready to be

fertilised. The egg can be fertilized for the 24 hours following release.

Measuring LH levels in the urine is therefore an effective way to optimise the chances of conception.

The OVULATION LH TEST allows you to detect a level of LH above 25 mIU/ml and therefore identify the days

you are most fertile.

##### TEST PRINCIPLE

OVULATION LH TEST is an immunochromatographic test that detects the luteinizing hormone (LH) through special gold-conjugated monoclonal antibodies included in the reactive strip.

<b>REFERENCES</b>
1. Su HW, Yi YC, Wei TY, Chang TC, Cheng CM "Detection of ovulation, a review of currently available methods" <i>Bioeng Transl Med.</i> 2017 May 16;2(3):238-246.
2. Leiva RA, Bouchard TP, Abdullah SH, Ecochard R.; "Urinary Luteinizing Hormone Tests: Which concentration threshold best predicts ovulation?" <i>Front public health</i> 2017, 5: 320.

##### CONTENT

<ul style="list-style-type: none"><li>* hermetically sealed aluminium pouches containing each:</li></ul> <ul style="list-style-type: none"><li>1 device for the OVULATION LH TEST including a protective cap</li> <li>1 desiccant bag</li></ul>	<ul style="list-style-type: none"><li>* 1 instructions for use leaflet</li></ul>
---	--

\* Do not open the sealed aluminium pouch until just before performing the test. Take care to open it as marked. refer to the "content" section on the external box.

\* The desiccant bag must not be used. Dispose of it with household waste without opening it.

\* Material required but not supplied: a device to measure time (i.e. timer, watch).

##### PRECAUTIONS

1. Read these instructions for use carefully before performing the test. The Test is reliable only if all the

instructions are followed correctly.

2. Keep the Test out of the reach of children.

3. Do not use the Test after the expiry date or if the package has been damaged.

4. Follow the indicated procedure and time strictly.

5. Store the Test components at a temperature between +4°C and +30°C. Do not freeze.

6. Use the test once only.

7. The test is for external use only. DO NOT SWALLOW.

8. In vitro diagnostic device for individual use.

9. After using, please dispose of all components according to your local waste disposal laws. Ask your pharmacist for advice.

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