PRIMA®

Your partner for rapid testing

PRIMA Lab SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
support@primalabs.ch
primalabs.ch

Qarad EC-REP BV
Pas 257, 2440 Geel - Belgium

100066IP-20P 100066IP-20P_IFU_90_5.1 03/2024


Table 1. Performance / Tabelle 1. Leistung / Tableau 1. Performances / Tabella 1. Prestazioni / Tabla 1. Rendimiento / Tabela 1. Desempenho / Tabela 1. Oztalante

METHOD	PREDICATE DEVICE				TOTAL RESULT
	Results	Deficient	Insufficient	Sufficient	
VITAMIN D RAPID TEST	Deficient	4	3	0	7
	Insufficient	0	53	2	55
	Sufficient	0	0	28	28
TOTAL RESULT		4	56	30	90
ACCURACY	>99.90 %	94.60 %	93.30 %	94.40 %	

REFERENCES / BIBLIOGRAFIE / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA

- Holick MF (March 2006). High prevalence of vitamin D inadequacy and implications for health. *Mayo Clinic Proceedings*. 81 (3): 353-73.
- Eriksen EF, Glerup H (2002). Vitamin D deficiency and aging: implications for general health and osteoporosis. *Biogerontology*, 3 (1-2): 73-7.
- Grant WB, Holick MF (June 2005). Benefits and requirements of vitamin D for optimal health: a review. *Alternative Medicine Review*. 10 (2): 94-111.
- Mayad MA. Vitamin D: a rapid review. *Dermatology News*. 2009, 21:25-30.

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mens. The specimens were correctly identified >99% of the time. **INTER-ASSAY** Between-run precision has been determined by 3 independent assays on the same 4 specimens: 10 ng/mL vitamin D, 30 ng/mL vitamin D, 45 ng/mL vitamin D, 100 ng/mL vitamin D standard samples. Three different lots of the VITAMIN D RAPID TEST have been tested using these specimens. The specimens were correctly identified >99% of the time. **SENSITIVITY AND CROSS-REACTIVITY** The VITAMIN D RAPID TEST can detect levels of Vitamin D in human whole blood as low as 30 ng/mL. The addition of Vitamin A, B, C, E, K and M showed no cross-reactivity.

PRECAUTIONS

- Read these instructions for use carefully before performing the test. The test is reliable if the instructions are followed correctly (reaction times and how the sample is collected).
- Do not open the sealed aluminum pouch until just before performing the test. The desiccant bag must not be used.
- Humidity and temperature can adversely affect the results. DO NOT FREEZE.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards during all procedures and follow standard procedures for proper disposal of specimens.
- The test is for external and single use only. DO NOT SWALLOW.
- Do not use the kit after the expiry date or if the packaging and/or internal components are damaged.
- It is suggested not to use the buffer beyond 6 months after opening the vial.

CONTENT

- 20 aluminum pouches containing 1 test device for the VITAMIN D TEST and a desiccant bag;
- 20 vials with dropper containing VITAMIN D TEST DILUENT sufficient for 20 tests;
- 20 pipettes for blood sampling;
- 1 Vitamin D Color Card;
- 1 instructions for use.

TEST PROCEDURE

Prepare the necessary material as follows: open the aluminum bag, remove the cassette and discard the desiccant bag. Place the cassette on a flat, dry and clean surface. Allow the test specimen, buffer and/or controls to reach room temperature (15-30 °C) prior to testing. **1** Prick the skin of the fingertip with a sterile lancet and wipe away the first sign of blood. Obtain a drop of blood and fill the capillary pipette to the black line. **2** Transfer 20 µl of whole blood to the sample well (S) – **fig. A1**. As an alternative to the pipette, 20 µl of venous sample can be dispensed with the laboratory micropipette – **fig. A2**.

- The VITAMIN D RAPID TEST provides only a semi-quantitative analytical result. A secondary analytical method must be used to obtain a confirmed result.
- It is possible that technical or procedural errors, as well as other interfering substances in the whole blood specimen, may cause erroneous results.
- The Cut-off for the test is 30 ng/mL with a deviation range of ± 4 ng/mL.
- As with all diagnostic tests, all results must be considered with other clinical information available to the physician.

QUALITY CONTROL

An internal procedural control is included in the test. A colored line appearing in the Control (C) line region is an internal positive procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

REFERENCES

The test contains anti-25 (OH) D antibody coated particles and 25 (OH) D antigen coated on the membrane.

PERFORMANCE CHARACTERISTICS

ACCURACY The VITAMIN D RAPID TEST has been compared with Predicate Device. The following result was tabulated (table 1): Accuracy 94.4%. **INTRA-ASSAY** Within-run precision has been determined by using 3 replicates of four specimens: 10 ng/mL, 30 ng/mL, 45 ng/mL and 100 ng/mL speci-

SUFFICIENT 30-100 (ng/mL) 75-250 (nmol/L)

Two colored lines appear, one line in the Control (C) region and faint colored line appears in the Test (T) region. The line intensity in Test (T) region is equal to or lighter than 30 ng/mL depicted on Vitamin D Color Card.

EXCESS >100 (ng/mL) >250 (nmol/L)

A coloured line appears in the Control (C) region. No line appears in the Test (T) region.

INVALID RESULT

The Control (C) line does not appear. An insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the Control line (invalid result). Review the procedure and repeat the test with a new device and a new sample.

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