
**PRIMA Lab SA**  
Via Antonio Monti 7  
CH-6828 Balerna - SWITZERLAND  
support@primalabs.ch  
primallabs.ch

Qarad EC-REP BV  
Pas 257, 2440 Geel - Belgium

800060IP-20P 800060IP-20P\_IFU\_90\_55 03/2024

**SYMBOLS / SYMBOLE / SYMBOLES / SIMBOLOGIA / SÍMBOLOS / SIMBOLOGIA / SYMBOLE**

**IVD** In vitro diagnostic device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / Dispositivo de diagnóstico in vitro / Dispositivo de diagnóstico in vitro / Wyrob przeznaczony do diagnostyki in vitro

**Read the instructions before use /** **Beachten Sie die Gebrauchsanweisung /** **Consulte le mode d'emploi /** **Leggere attentamente le istruzioni per l'uso /** **Lea las instrucciones antes de usar el producto /** **Lea as instruções antes da utilização /** **Prez uvaženie prosím si prečítajte inštrukcie.**

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**CE REP** Dansu A/S Sandkæmper 1-3, Gårdsvej, DK-3660 Stejneger, Denmark

**0123** Jiangsu Suyun Medical Materials Co. Ltd No. 18 Jin Dao Road Dava Industrial Park 222002 Lianyungang, Jiangsu Province, People's Republic of China

**EC REP** Shanghai International Trading Corp. GmbH (Europe) Ellertstrasse, 80, D-20537 - Hamburg, Germany

Your partner for rapid testing

**IVD for professional use only**  
**IVD nur für den professionellen Gebrauch**  
**Div à usage professionnel uniquement**  
**IVD solo per uso professionale**  
**IVD apenas para uso profesional**  
**IVD wyłącznie do użytku profesjonalnego**

**REFERENCES / BIBLIOGRAFIE / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA**

- Anjos LM, Marcondes MB, Lima MF, Mondelli AL, Okoshi MP. Streptococcal acute pharyngitis. Rev Soc Bras Med Trop 2016; 47(4):409-13.
- Castro SA, Dorfmueller HC. A Brief Review on Group A Streptococcus pathogenesis and vaccine development. R Soc Open Sci 2022.
- Banerjee S, Ford C. Rapid Tests for the Diagnosis of Group A Streptococcal Infection: A Review of Diagnostic Test Accuracy, Clinical Utility, Safety, and Cost-Effectiveness. Ottawa (ON), Canadian Agency for Drugs and Technologies in Health; 2018 May. CADTH Rapid Response Reports.
- Llor C, Bjerrum L, Munkáçs A, Cots JM, Hernández S, Maragás A. HAPPY A/DH Investigators. Access to point-of-care tests reduces the prescription of antibiotics among antibiotic-requesting subjects with respiratory tract infections. Respir Care 2014 Dec.
- Georgers N, Nakhouli MD and John Hickner MD, MSc. Management of Adults with Acute Streptococcal Pharyngitis: Minimal Value for Backstop Strep Testing and Overage of Antibiotics. J Gen Intern Med. 2013 Jun; 28(6):830-834.
- Llor C, Bjerrum L, Munkáçs A, Cots JM, Hernández S, Maragás A. HAPPY A/DH Investigators. Access to point-of-care tests reduces the prescription of antibiotics among antibiotic-requesting subjects with respiratory tract infections. Respir Care 2014 Dec.
- Georgers N, Nakhouli MD and John Hickner MD, MSc. Management of Adults with Acute Streptococcal Pharyngitis: A Systematic Review and Meta-Analysis. PLoS One 2014; 9(11).
- Parks T, Barnett L, Jones N. Invasive streptococcal disease: a review for clinicians. British Medical Bulletin. Volume 115, Issue 1, September 2015, Pages 77-89.

**Table 1. Performance / Tabelle 1. Leistung / Tableau 1. Performances / Tabella 1. Prestazioni / Tabla 1. Rendimiento / Tabella 1. Desempenho / Tabla 1. Oztarimio**

STREP A RAPID TEST	REFERENCE METHOD	
	POSITIVE	NEGATIVE
POSITIVE	34	36
NEGATIVE	3	299
	37	201
		238

**SPECIFICITY = (199/201)\*100 = 99.00% (CI 95%: 97.63 - 100%)**  
**SENSITIVITY = (34/37)\*100 = 91.89% (CI 95%: 78,70 - 97,20%)**  
**ACCURACY = ((233/238)\*100 = 97.90% (CI 95%: 96,08 - 99,72%)**

**Table 2. Interfering substances / Tabelle 2. Beeinträchtigende Substanzen / Tabella 2. Substances interferentes / Tabella 2. Sostanze interferenti / Tabla 2. Sustancias que interfieren / Tabla 2. Substâncias interferentes / Tabla 2. Substancje zakłócające**

BACTERIA	ATCC	LOT	CONCENTRATION
<i>Bordetella pertussis</i>	9340	BP-02	10 <sup>7</sup> CFU/mL
<i>Campylobacter jejuni</i>	29428	CJ-02	10 <sup>7</sup> CFU/mL
<i>Candida albicans</i>	MUT6298	CA-01	10 <sup>7</sup> CFU/mL
<i>Clostridium perfringens</i>	13124	CP-02	10 <sup>7</sup> CFU/mL
<i>Corynebacterium diphtheriae</i>	11913	CD-02	10 <sup>7</sup> CFU/mL
<i>Enterobacter cloacae</i>	MUT6313	EC-02	10 <sup>7</sup> CFU/mL
<i>Helicobacter pylori</i>	43504	HP-01	10 <sup>7</sup> CFU/mL
<i>Hemophilus influenzae</i>	49766	HI-01	10 <sup>7</sup> CFU/mL
<i>Klebsiella pneumoniae</i>	MUT6314	KP-02	10 <sup>7</sup> CFU/mL
<i>Proteus vulgaris</i>	MUT6312	PV-02	10 <sup>7</sup> CFU/mL
<i>Pseudomonas aeruginosa</i>	MUT6310	PA-02	10 <sup>7</sup> CFU/mL
<i>Shigella flexneri</i>	12022	SF-02	10 <sup>7</sup> CFU/mL
<i>Staphylococcus aureus</i>	MUT6311	SA-02	10 <sup>7</sup> CFU/mL
<i>Streptococcus agalactiae</i>	12386	SAG-02	10 <sup>7</sup> CFU/mL
<i>Streptococcus galloyticus</i>	12388	SD-02	10 <sup>7</sup> CFU/mL
<i>Streptococcus pyogenes</i>	9808	SG-02	10 <sup>7</sup> CFU/mL
<i>Yersinia enterocolitica</i>	51871	YE-02	10 <sup>7</sup> CFU/mL

**EN - INSTRUCTIONS FOR USE STREPTOCOCCUS A**

The throat infection with group A β-hemolytic Streptococcus (also known as *Streptococcus pyogenes*) is a Gram-positive bacterium that constitutes the most frequent bacterial cause of acute pharyngitis and it is also responsible for some skin infections such as impetigo and erysipelas. Most sore throats, on the other hand, are caused by viral infections, which are cured without the need of antibiotic treatments, unnecessary in these cases. Streptococcus A infections can resolve in a few days without treatment, although doctors generally prefer to prescribe antibiotics to prevent infection-related complications that can also be serious, such as rheumatic fever or acute glomerulonephritis. Group A Streptococcus continues to be a focus of interest not only because of their causal role in acute streptococcal pharyngitis and other pyogenic infections but also because of its association with poststreptococcal sequelae, specifically, acute rheumatic fever and acute glomerulonephritis. In order to properly treat the disease using antibiotic therapy, it is important to use an accurate diagnostic method to identify the pathogenic agent. For the screening of Group A Streptococcus infection several methods are currently used including Bacitracin sensitivity tests on a sheep blood agar plate, latex agglutination and enzyme immunoassay. STREP A RAPID TEST is a rapid immunological test for the qualitative and visual detection of group A Streptococcus bacterium in throat infections. The test helps to quickly identify whether a sore throat is caused by Strep A rather than by other pathogens (usually viruses) that may not require antibiotic treatment.

**PRINCIPLE OF THE TEST**

STREP A RAPID TEST detects specific Streptococcus A antigens (antibodies) using special gold conjugate antibodies embedded to a test strip. The specimen, which is pharyngeal swab, is taken with a disposable sterile rayon swab. The specimen is placed in a plastic collection tube where the user dispenses a specific volume of 2 extraction reagents

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(in order to extract specific Streptococcus A antigens). After an incubation time, the obtained sample-solution is dispensed onto the test cassette mainly with the provided reagents are:  
- R1 solution mainly containing Sodium Nitrite and Phenol Red.  
- R2 solution mainly containing Acetic Acid.  
Once mixed, the two solutions create Nitrous Acid, acting as Streptococcus A specific antigen extractor. In case of a Strep A concentration equal or higher than 10<sup>7</sup> CFU/mL in the patient's sample, the specific antigens present in the sample solution are firstly bound by the anti-Strep A antibodies labeled with gold particles. The immunocomplexes are selectively captured by the test line (T line, made of purified anti-Strep A antibodies), thus showing a red-purple band appearance in correspondence of the T line, pointing out a positive result. In case of absence of Strep A or when its levels in patient's sample are lower than the cut-off value (10<sup>7</sup> CFU/mL), not enough specific immunocomplex is formed. Consequently, only one band, the C line, is visible indicating a negative result. To serve as a procedural control, a colored line will always appear in the Control line region, indicating that the test procedure was being properly performed and that the test components and reagents have operated as intended.

**PRECAUTIONS**

- Read these instructions for use carefully before performing the test. The test is reliable if the instructions are followed correctly (reaction times, how the sample is collected and the extraction instructions shown).
- Humidity and temperature can adversely affect the results. DO NOT FREEZE.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards during all procedures and follow standard procedures for proper disposal of specimens.
- Do not use the kit after the expiry date or if the packaging and/or internal components are damaged.
- The test is for external use only. Do not swallow solutions 1 and 2 and the other items in the box.
- If liquids come into contact with the eyes, skin or other mucosae, rinse immediately and thoroughly with water.
- Do not switch the caps of the reagent vials or mix components from different test kits.
- Not suitable for use on children under the age of 2 years.
- All kit components are single-use, with the exception of the work-station, instructions for use and solutions 1 and 2, which must be stored correctly closed with their caps.
- Do not open the sealed elements (test, swab and tongue depressor) until just before performing the test. Take care to open as marked. Dispose of the desiccant sachet inside the aluminium pouch, without opening it.

**CONTENT**

- 20 hermetically sealed aluminium pouches each containing 1 STREP A TEST cassette and 1 desiccant bag;
- 20 sterile swabs;
- 20 sterile wooden tongue depressors;
- 1 workstation;
- 20 empty plastic extraction tubes with dropper;
- 1 dropper vial containing enough STREP A TEST-R1 solution to perform 20 tests. R1 contains 0.55% Acetic Acid, Nitrite, H3O2: Harmful in case of ingestion.
- P264: Wash carefully your hands after use.
- P270: Do not eat, drink or smoke while using this product.

- 1 dropper vial containing enough STREP A TEST-R2 solution to perform 20 tests. R2 contains 0.55% Acetic Acid;
- 1 Instructions for use leaflet.
- Required material not supplied: a timer (e.g. timer, watch).

**SPECIMEN COLLECTION AND PREPARATION**

Bring tests, specimens, reagents to room temperature (15-30°C) prior to testing. To avoid cross contamination, do not allow the tips of the reagent bottles to come into contact with specimen material.  
1) Place the workstation on a clean, dry, flat surface. Place the plastic test tube in one of the compartments of the workstation – **fig. C**.  
2) Collecting the sample: Open the package containing the tongue depressor, leaving it inside the package and next to you. Open the package containing the sterile swab. Avoid touching the cotton tip and remove the swab using the plastic handle. Tilt the patient's head backwards and let open the mouth as much as possible. Pick up the tongue depressor with one hand and flatten the tongue. Use the other hand to place the swab close to the throat. Touch the back of the throat - the area around the tonsils and any redness or painful part (soft palate, uvula) with the cotton tip. Rotating the swab is recommended because it increases the amount of collected sample – **fig. B**.  
It is recommended that swab specimens be processed as soon as possible after collection.

**TEST PROCEDURE**

- After collecting the sample, insert the swab's cotton tip into the supplied plastic test tube that you had previously placed in the workstation – **fig. C**.
- Open the vial STREP A TEST - R1 by unscrewing only the blue cap (leave the white part firmly sealed), then add 4 drops into the test tube and close it with the cap – **fig. E**.
- Next, open the vial STREP A TEST - R2 by unscrewing only the green cap (leave the white part firmly sealed), then add 4 drops into the test tube and close it with the cap – **fig. E**.  
**WARNING: the addition of R2 to R1 leads to the formation of a solution for which the indications described in point 8 and 9 of the "precautions" apply.**
- Holding the swab's plastic handle, rotate the swab against the sides of the test tube 10 times to thoroughly mix the solution. Leave the swab to incubate for 2.5 minutes. – **fig. F**.
- At the end of the incubation period, remove the test tube (with the swab and extraction fluid still in it) from the workstation. Using your thumb and index finger, press the sides of the test

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**RESULTS INTERPRETATION**

READ THE RESULTS AFTER 5 MINUTES. DO NOT READ THE RESULTS AFTER 10 MINUTES. The intensity of the line colors is not relevant for the purposes of interpretation of the test's results.

**POSITIVE RESULT**

Two coloured lines appear in the reading window by the Test (T) and Control (C) signs. The T line may be less intense (lighter) than the C one. This result means that the test has detected the presence of Streptococcus A specific antigens in the sample.

**NEGATIVE RESULT**

A coloured line appears in the Control (C) region. No line appears in the Test (T) region. This result is negative and means that no Streptococcus A antigens have been detected in the sample from the pharynx or that their amount is too low to be detected.

**INVALID RESULT**

The Control (C) line does not appear. An insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the Control line (invalid result). Review the procedure and repeat the test with a new device and a new sample.

**LIMITATIONS**

- STREP A RAPID TEST should be used for the determination of Streptococcus A antigen on pharyngeal swabs only. It cannot determine either the quantitative value or the percentage increase in Streptococcus A antigen concentration.
- This test only indicates the presence of Streptococcus A antigen in the specimen due to viable or non-viable Group A Streptococcus bacteria, it does not differentiate asymptomatic carriers of group A Streptococcus from those with a symptomatic infection.
- False negative results may occur due to low antigen concentration that is not high enough to be detected by the test, or due to improper specimen collection or storage.
- An excessive amount of blood or mucus on the swab may interfere with the test results and lead to a false positive result. Avoid contact of the swab with the tongue, cheeks and teeth and any oral areas with open wounds at the time of sampling.
- As with all diagnostic tests, results should be interpreted in the light of the patient's overall clinical history.
- If the reading time (5 minutes) is not strictly respected, wrong results may be obtained.

**H302:** Harmful in case of ingestion.  
**P264:** Wash carefully your hands after use.  
**P270:** Do not eat, drink or smoke while using this product.

- 1 dropper vial containing enough STREP A TEST-R1 solution to perform 20 tests. R1 contains 0.55% Acetic Acid, Nitrite, H3O2: Harmful in case of ingestion.
- P264: Wash carefully your hands after use.
- P270: Do not eat, drink or smoke while using this product.

**PERFORMANCE CHARACTERISTICS**

Two swabs (test and reference method) were obtained from each patient and collected by physicians participating in the study. The comparator swab was analysed according to its instruction for use. The results were compared to those of the reference method and reported as positive and negative in Table 1. The performance of the test is reported in terms of Sensitivity, Specificity and Accuracy. **CROSS-REACTIVITY** Each sample was tested as indicated in Table 2, according to STREP A RAPID TEST instruction for use in force, at the indicated concentration. None of the bacteria tested in this study, at the concentration given in Table 2, have shown cross reactions or interferences with STREP A RAPID TEST.  
**MEASURING RANGE** The limit of detection study was conducted on Certified Positive ATCC Streptococcus pyogenes®. The LOD of STREP A RAPID TEST is 10<sup>7</sup> CFU/mL. No adverse effect on T-line formation (or on the results) was observed for group A Streptococcus concentrations up to 5 x 10<sup>7</sup> CFU/mL.

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