## Test your health

 **PRIMA LAB SA**  
Via Antonio Monti 7  
CH-6828 Balerna - SWITZERLAND

support@primlabsa.ch  
primlabsa.ch

 **EC REP** QbD RepS BV  
Groenenborgerlaan 16  
2610 Wilrijk - Belgium

100013\_IFU\_11\_8.0 11/2024

 **REF** 100013-X

INSTRUCTIONS FOR USE  
GEBRAUCHSANWEISUNG  
NOTICE D'UTILISATION  
ISTRUZIONI PER L'USO

## VAGINAL PH TEST

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 **Copan Italia SpA**  
Via Perotti 10  
25125 Brescia  
0123 Italy

### BIBLIOGRAPHY / BIBLIOGRAPHIE / BIBLIOGRAFIA

1. Pavletic, A. J., Hawes, S. E., Geske, J. A., Bringe, K., & Polack, S. H. (2004). Experience with routine vaginal pH testing in a family practice setting. *Infectious Diseases in Obstetrics and Gynecology*, 12(2), 63–68.
2. Huppert, J. S., Hesse, E. A., Bernard, M. C., Bates, J. R., Gaydos, C. A., & Kahn, J. A. (2012). ACCURACY AND TRUST OF SELF-TESTING FOR BACTERIAL VAGINOSIS. *The Journal of Adolescent Health : Official Publication of the Society for Adolescent Medicine*, 51(4), 400–405.

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
## ENGLISH SPECIMEN COLLECTION AND PREPARATION

The Vaginal pH Test can be performed using female vaginal swab specimen. Allow the test and the swab specimen to reach room temperature (15-30°) prior to testing. Wash hands before starting the test. Open the package, remove the swab from pouch, do not touch the head of the swab (polyester head). Separate the labia so the vagina is exposed. Insert the head of swab gently into the vagina and gently press the swab to a point, when approximately 2 inches of swab (approximately length of an AA size battery is into the vagina for approx. 10 seconds). A time of approximately 10 seconds is important to ensure that the head of swab becomes well moistened with vaginal secretions. Sufficient moisture is very important for the test evaluation. Withdraw the swab from the vagina for testing. Apply the head of swab onto pH area of the test panel for at least 5 times, make sure enough secretion is applied on the pH area of the test panel. Read the result immediately while the measurement zone is still wet. Compare the color in the pH measurement zone with the color scale on the foil pouch and note the value on the color scale that is closest to the color on the test unit. After evaluation, dispose of the used panel and swab.

## RESULTS INTERPRETATION

Read the result immediately while the measurement zone is still wet. Results are obtained by visually comparing the PH measurement zone with the color scale on the foil pouch.

< 1 min.

## NORMAL RANGE

Vaginal pH levels in the range 3.8-4.4 are within the normal range. The color of indicator is Yellow-Green.

## ABNORMAL RANGE:

Vaginal pH levels over 4.4 or under 3.8 are outside of the normal/healthy range. When pH level is over 4.4, color of the indicator is Green-Blue, when pH levels is under 3.8, it is bright Yellow.

## LIMITATIONS

- Some menopausal women may have an elevated vaginal pH, and therefore may receive an abnormal result with Vaginal pH Rapid Test Panel even though they are not suffering from a bacterial or parasitic infection.
- Certain circumstances(menstrual period, sex intercourse, use of vaginal medical products) can temporarily alter the pH level of vagina and may lead to false results. Therefore, please note the precautions when performing the test:
- Do not wash or shower the area immediately before taking the measurement.
- Be aware that urine may cause a false test result.
- Never begin any treatment before you have discussed the result of the test with a doctor.
- DO NOT use this test until 72 hours after the application of vaginal preparations like contraceptive creams and vaginal medical products (antifungal agents, vaginal suppositories, creams, gels, douching etc.).
- DO NOT use this test until 48 hours after sexual intercourse.
- DO NOT use this test until 5 days after the period is over.
- In cases of improper use, rupture of hymen may occur as it occurs when using tampons.
- This is NOT a test for sexually transmitted diseases, such as HIV, chlamydia, herpes, gonorrhea or syphilis.

## F.A.Q. – QUESTIONS AND ANSWERS

## COULD I USE THE TEST AT ANYTIME?

No. There are some restrictions on when to use this test because at certain circumstances can temporarily alter the pH level of the vagina and may lead to false results. Therefore, please find the "Limitations" section in this package insert and use the test.

## CAN THE POSITIVE RESULTS SHOW THAT THE SUBJECT HAS BACTERIAL VAGINOSIS?

Not necessarily. Some menopausal women may have an elevated vaginal pH,hence may receive a positive result even without suffering from bacterial or parasitic infection. Initiate any treatment, only after consultation with doctor.

## WHAT IS THE REASON THAT MAY RESULT IN WRONG TEST RESULTS?

Some factors can affect the test results, including inaccurate specimen collection, insufficient moistening of the swab, insufficient application of specimen on the test area.

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