Test your health

PRIMA LAB SA


Via Antonio Monti 7

CH-6828 Balerna - SWITZERLAND

support@primalabs.ch
primalabs.chQBD RepS BV - Groenenborgerlaan 16
2610 Wilrijk - Belgium

800062-1_FU_94_8.0 11/2024

REF 800062-1FSP

EN - INSTRUCTIONS FOR USE

WHAT IS HIV?

The human immunodeficiency virus (HIV) is a pathogen which attacks and suppresses the immune system by affecting in particular white blood cells. The possible transmission pathways for HIV are:

- Sexual through sexual intercourse not protected by an efficient method of prevention (e.g. condoms).
- Blood exchange of syringes, transfusions of infected blood.
- Vertical from mother to baby during pregnancy, at birth and, more rarely, through breastfeeding.

Without taking specific drugs, the HIV virus weakens the immune system leading to the acquired immunodeficiency syndrome (AIDS).

WHAT IS THE AIDS?

AIDS identifies an advanced clinical stage of HIV infection. It is a syndrome that can occur in people with HIV even several years after the infection, when the immune system's efficiency dramatically falls, and the body loses its ability to fight even the most common infections. The AIDS can be prevented by taking antiretroviral therapy at an early stage.

TEST PRINCIPLE & INTENDED USE

HIV 1/2 SELF-TEST is a self-test that is able to specifically detect the presence in a blood sample of anti-HIV 1 and/or 2 antibodies. It is the sampled drop of blood contains anti-HIV antibodies, they are captured by the test producing a red-coloured line in the Test line region (T). Otherwise, if the sample does not contain such antibodies, no colour will appear in the Test line region (T). The absence of colour in the Test line region (T) is also possible when the levels of circulating antibodies are below the limit of detection, leading to false negative results, especially for HIV-1 Group O (prevalence of HIV 1 Group O worldwide 1%). A red-coloured line will always appear in the Control line region (C) confirming that the test has been correctly operated. **IT IS RECOMMENDED TO NOT TAKE ANY DECISION OF MEDICAL RELEVANCE AFTER PERFORMING THIS TEST WITHOUT FIRST CONSULTING YOUR MEDICAL PRACTITIONER.**

F.A.Q. - QUESTIONS AND ANSWERS

Do not take any decision of medical relevance without first consulting your medical practitioner after performing the test. The test detects the presence of antibodies human immunodeficiency virus 1 and 2 (HIV-1/HIV-2) by means of specific antigens and colloidal gold nanoparticles embedded in the test strip.

WHEN CAN THE TEST BE USED? HIV 1/2 SELF-TEST can be performed at any time of the day. The test can be carried out by anyone who has been or suspects to have been exposed to the virus, people with HIV symptoms or anyone who wants to be aware of a possible HIV infection. Antibody tests can take 23 to 90 days (window period) to detect HIV infection after an exposure. **HIV 1/2 SELF-TEST is not for use in emergency situations.**

DOES HIV 1/2 SELF-TEST DETECT OTHER DISEASES? No. The HIV 1/2 SELF-TEST only detects the presence of antibodies to human immunodeficiency virus 1 and 2 (HIV-1/HIV-2) in the blood sample, which is indicative of virus infection and subsequent immune response. HIV 1/2 SELF-TEST cannot be used to detect other sexually transmitted infections.

WHAT IS THE DIFFERENCE BETWEEN HIV AND AIDS? HIV is a virus that attacks the human immune system. Under normal circumstances, immune system cells defend the body against various pathogens - including viruses. However, HIV destroys exactly those cells that should protect the body against it. AIDS is a late consequence of untreated HIV infection. Being HIV positive and having AIDS is not the same thing. The HIV multiplies and damages the immune system until AIDS develops.

CAN THE RESULT BE FALSE? The result is correct as long as the instructions are carefully followed. However, the result may not be correct if the device comes into contact with other liquids before or after, if the amount of blood and/or diluent is insufficient, if the number of drops dispensed into the well is incorrect or if the reading time of the result is not respected. The plastic pipette supplied allows to be sure that the volume of collected blood is correct. False positive results (the test incorrectly indicates a positive result even if the anti-HIV antibodies are not present) may be due to: *Mycoplasma pneumoniae*, *Borrelia*, *Hepatitis C Virus (HCV)*, *Human T-lymphotropic virus (HTLV)*, *Epstein-Barr Virus/Nuclear Antigen (EBV/EBNA)*. The main consequences of a false positive result are the person could experience a strong negative stress and would have to undergo a second unnecessary laboratory test to have the diagnosis confirmed. Due to the lack of samples with very high concentrations of antibodies, the manufacturer cannot exclude the possibility that very high concentrations of antibodies against HIV-1/2 in the blood sample could lead to false negative results (the test incorrectly indicates a negative result even if the anti-HIV antibodies are actually present).

False negative results are also possible in the case of antibodies below the limit of detection, which may be particularly relevant for HIV-1 Group O (prevalence of HIV 1 Group O worldwide 1%). The main consequences of a false negative result are a delay in receiving the proper therapy, transmission of the virus to other persons.

HOW TO INTERPRET THE TEST IF THE COLOR AND INTENSITY OF THE TEST AND CONTROL LINES ARE DIFFERENT? The colour and intensity of the lines are not relevant for the interpretation of the results. Refer to the indications and to the images in the section "Results Interpretation" to correctly evaluate the result.

IF READ AFTER 10 MINUTES IS THE RESULT RELIABLE? No. The test should be read exactly 10 minutes after the end of the procedure. Results read after 10 minutes may not be correct (false positive results may appear, see question "Can the result be false?" where the consequences of false positive results are explained).

WHAT SHOULD BE DONE IF THE RESULT IS POSITIVE? If the result is positive, it means that antibodies to HIV 1 and/or HIV 2 are present in the blood sample. You should consult your physician and report the result. obtained with this test. If the result is not sufficient for the initial diagnosis, the development of HIV infection, the result should be confirmed by your doctor. **AVOID ANY ACTIVITY THAT COULD TRANSMIT HIV TO OTHERS** until you have received the results of further diagnostic investigations in accordance with your doctor. **PROTECT YOURSELF AND OTHERS**, ask your doctor for more information on how to prevent the transmission. False positive results may be obtained due to undesired positive reaction of the test towards other pathogens.

WHAT TO DO IF THE RESULT IS NEGATIVE? If the result is negative, the test did not identify the presence of anti-HIV 1 and/or 2 antibodies. This result may mean that HIV infection has not occurred or is at an early stage where specific antibodies have not been yet developed. False negative results may be obtained with patients receiving antiretroviral or corticosteroid therapy (e.g. conditions of immunosuppression). If your result is negative but you think you have been exposed or you have symptoms, please contact your doctor to ensure that you are not in the window period (3 months), which is the time required for the development of a detectable amount of HIV antibodies. Please avoid any activity that may cause HIV transmission to others.

ALCOHOL CLEANSING GAZE / ALKOHOL-REINIGUNGSGAZE / GAZE NETTAYANTE / FALCOLO / GARZA DETERGENTE ALCOHOLICA / GAZA NASCZYDZAJĄCA ALKOHOLEM / REINIGUNGSGASS MIT ALKOHOL / ALKOHOLHOLDIG REINIGUNGSGASSE

Phenox IntraEurope Healthcare Manufacturers P.V. Ltd.
Ei-209, Shi Mahale Road, Electronic Zone, MIDC,
TTC Industrial Area, Mahape, Navi Mumbai 400 710 MH India

Adena Ltd., Tower Business Centre,
2nd Flr, Tower Street, Svitkar, BKR 4013, Malta

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