Test your health

PRIMA LAB SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
support@primalabsa.ch
primalabsa.ch

EC REP QbD Rep5 BV
Groenenborgerlaan 16
2610 Wilrijk - Belgium

800080_IFU_11_8.0 10/2024 **REF** 800080-x

Registered trademarks.

INSTRUCTIONS FOR USE GEBRAUCHSANWEISUNG NOTICE D'UTILISATION ISTRUZIONI PER L'USO

PROSTATE PSA TEST

Self-test for the determination of the Prostate-Specific Antigen (PSA) in whole blood samples

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EC REP Hualan Tianda Medical Instruments Co., Ltd., No.106 East Songjiang Road, Huaiyin, Economic/Technological Development Zone, 223300 - Hualan City, Jiangsu - P.R. China

EC REP Riomavix Sociedad Limitada
Calle de Almansa 55, 1D
Madrid 28039, Spain

CE Alcohol cleansing gauze
Alkohol-Reinigungs-gaze
Gaze nettoyante à l'alcool
Garza detergente alcolica

CE Phoenix Innovative Healthcare Manufacturers Pvt. Ltd., E1-209, Shi Mahape Road, Electronic Zone, MIDC, TTC Industrial Area, Mahape, Navi Mumbai 400 710 MH I India

EC REP Advena Ltd., Tower Business Centre, 2nd Flr, Tower Street, Swatar, BKR 4013, Malta

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CE
Legal manufacturer
Hersteller
Fabricant
Fabbricante

PROSTATE AND PSA

The Prostate Specific Antigen (PSA) is a protein produced by the epithelial cells of the prostate gland. Its function is to fluidify the ejaculate and therefore increase sperm mobility. PSA also plays a role in dissolving cervical mucus, allowing sperm to enter the uterus. Small quantities of PSA are normally found in the blood of men with a healthy prostate. Prostate pathologies (such as prostatitis, benign prostatic hyperplasia, cancer) lead to an increase in the blood levels of PSA.

PSA levels also increase naturally with age, after prostate surgery or after specific diagnostic tests (rectal examination, prostate biopsy, prostate ultrasound scan, etc.).

WHO IS THE PROSTATE PSA TEST FOR?

- For men who do not have symptoms, are aged over 45, who wish to check their prostate;
- For men with urinary tract symptoms (difficultly urinating, pain, increased frequency).

FA.Q. – QUESTIONS AND ANSWERS

HOW DOES THE PROSTATE PSA TEST WORK? The Prostate-Specific Antigen is a protein secreted normally by the prostate and circulating in our body in small quantities. Values above normal may be a sign of the presence of pathological conditions in the prostate, which may be of infectious or inflammatory origin (prostatitis), benign (benign prostatic hyperplasia) or malignant (cancer). The PROSTATE PSA TEST is a rapid qualitative test for measuring the levels of Total PSA in whole blood. The test uses a unique combination of specific antibodies that can selectively identify PSA. A positive result means PSA values are above 4 ng/mL. In some cases, samples with PSA concentrations between 3 ng/ml and 3.9 ng/ml can yield positive results.

WHEN CAN THE TEST BE USED? The PROSTATE PSA TEST can be performed at any time of the day, but in some cases it is advisable to wait for a specific period of time in order to avoid a false positive result. The blood concentration of PSA may increase in the cases listed below:

Activity	Waiting times
prolonged periods of time spent on a bicycle/motorbike	 24 hours
sexual intercourse/ejaculation	 24-48 hours
insertion of a urinary catheter	 2-3 days
transrectal ultrasound scans	 2-3 days
infections of the urinary tract	 1 week

CAN THE RESULT BE INCORRECT? The result is correct if the instructions are followed carefully. However, the result may not be correct if: the device comes into contact with other liquids before being used, if part of the diluent is accidentally spilt, if the amount of blood is not enough or the number of drops dispensed in the well is more than 3.

HOW DO I INTERPRET THE TEST IF THE COLOUR AND INTENSITY OF THE TEST AND CONTROL BANDS ARE DIFFERENT? The colour and the intensity of the bands are not important for the interpretation of the result. The bands must be full and homogeneous. The test is positive regardless of the intensity of the colour of the test band.

IS THE RESULT RELIABLE IF READ AFTER 10 MINUTES? No. The test must be read 5 minutes after the procedure is completed.

WHAT SHOULD I DO IF THE RESULT IS POSITIVE? If the result is positive, the concentration of PSA in the blood is above normal levels and you must consult a physician, who will establish what to do next. PSA IS A SPECIFIC PROSTATE MARKER BUT IT DOES NOT NECESSARILY IMPLY THE PRESENCE OF NEOPLASIA.

WHAT SHOULD I DO IF THE RESULT IS NEGATIVE? If the result is negative, the concentration of PSA is normal. A negative result does not exclude the possibility of prostatic pathologies. If the result is negative but you have urinary tract symptoms (difficulty urinating, pain, increased frequency) you should talk to your physician.

HOW ACCURATE IS THE PROSTATE PSA TEST? The test is very accurate. Assessment reports show a concordance rate of 92.34% (95% CI: 87.93-95.23%) with reference methods. Despite the reliability of the test, false positive or negative results may occur.

MEASURING RANGE	SENSITIVITY*	SPECIFICITY*	OVERALL ACCURACY*
From 4 ng/mL to 5000 ng/mL	97.22% (CI 95%: 92.15-99.05)	87.13% (CI 95%: 79.22-92.32)	92.34% (CI 95%: 87.93-95.23)

*clinical study with 209 participants enrolled

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







PRECAUTIONS

- Read these instructions for use carefully before performing the test. The Test is reliable only if all the instructions are followed correctly.
- Keep the Test out of the reach of children.
- Do not use the Test after the expiry date or if the package has been damaged.
- Follow the procedure exactly, using only the specified quantities of blood and diluent.
- Store the Test components at a temperature between +4°C and +30°C. Do not freeze.
- Use the test and lancing device once only.
- The test is for external use only. DO NOT SWALLOW.
- In vitro diagnostic device for individual use.
- Not recommended for people who take anti-coagulant medications (blood thinners) or people suffering from haemophilia.
- After using, please dispose of all components according to your local waste disposal laws. Ask your pharmacist for advice.

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
CONTENT

* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the external box.

* Material required but not supplied: cotton wool, a device to measure time (i.e. timer, watch).
* Do not open the sealed aluminium bag until just before performing the test. Take care to open it as marked.
* The desiccant packet must not be used. Dispose of it with household waste without opening it.

TEST PROCEDURE

- Wash hands with soap and warm water, rinse with clean water and allow to dry. If this is not possible, use the provided gauze as an alternative.
Note: The use of warm water facilitates capillary blood collection as it induces vasodilation.
- Prepare the necessary material as follows: open the aluminium pouch, take out only the test cassette and throw away the desiccant bag. Open the plastic packet containing the pipette. **FIG. A**
- Carefully rotate the protective cap of the sterile lancet 360° without pulling it. Extract and discard the released cap. **FIG. B**
- Carefully massage the finger chosen for the puncture (the side of the ring finger is recommended). It is important that the massage is done from the palm of the hand to the phalanx, to improve blood flow. Press the open end of the lancet (the side the cap has been extracted from), against the fingertip **FIG. C** The tip of the lancet automatically retracts after use. If the lancet does not work properly, discard it and use the second one supplied. If the second one is not required, it can be disposed of without special precautions.
- Holding the hand down, massage the finger until a large drop of blood forms. It is important to massage from the palm of the hand to the phalanx to improve blood flow. **FIG. D**
- Take the pipette **without pressing the bulb**. Two sampling methods are suggested:
FIG. E1: hold the pipette horizontally **without pressing the bulb** place it in contact with the drop of blood, it will enter the pipette by capillarity. Move the pipette away when the black line is reached. If there is not enough blood, continue to massage the finger until the black line is reached.
FIG. E2: place the pipette on a clean, flat surface with the tip protruding from the shelf, then place the drop of blood in contact with the pipette, it will enter by capillarity. If blood is not sufficient, continue massaging the finger until the blood has reached the black line.
Avoid, as far as possible, to continuously move the tip of the pipette away from the finger in order to prevent the formation of air bubbles.
- Place the blood collected with the pipette into the well indicated on the cassette (S) by pressing the pipette bulb. **FIG. F**
- Unscrew the blue cap from the dropper vial (leave the white cap tightly screwed on). Deposit **2 drops** into the well indicated on the cassette (S), wait for 5 seconds between the first and second drop. **FIG. G**
- Wait **5 minutes** and read the results referring to the next section for result interpretation.

RESULTS INTERPRETATION

READ THE RESULT AFTER EXACTLY 5 MINUTES
The intensity of the line colours is not relevant for the purposes of interpretation of the Test's results.

POSITIVE RESULT
Two coloured bands appear in the reading window by the T (test) and C (control) signs. The T band may be less intense (lighter) than the C line.
This result means that the value of Prostate-Specific Antigen (PSA) in the blood is higher than normal, so you must consult a physician.

NEGATIVE RESULT
A coloured band appears only under the C (control) sign.
This means that the value of Prostate-Specific Antigen (PSA) in the blood is normal.
A negative result does not exclude the possibility of prostatic pathologies.

NON VALID RESULT
No bands appear or there is a line only under the T (test) sign and not under the C (control) sign.
In this case it is not possible to interpret the result of the test, which must be considered not valid. Repeat the test with a new blood sample.

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