PRIMA Lab SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
support@primlabsa.ch
primlabsa.ch

EC REF QbD RepS BV,
Groenenborgerlaan 16
2610 Wilrijk - Belgium

REF 1000661P-20P 1000661P-20P_IFU_90_8.0 11/2024

Your partner for rapid testing


Table 1. Performance / Tabelle 1. Leistung / Tableau 1. Performances / Tabela 1. Prestazioni / Tabela 1. Rendimiento / Tabela 1. Desempenho / Tabela 1. Oztatante.

VITAMIN D RAPID TEST	METHOD			PREDICATE DEVICE			TOTAL RESULT
	Results	Deficient	Insufficient	Sufficient			
	Deficient	4	3	0	7		
	Insufficient	0	53	2	55		
	Sufficient	0	0	28	28		
TOTAL RESULT		4	56	30	90		
ACCURACY		>99.90 %	94.60 %	93.30 %	94.40 %		

REFERENCES / BIBLIOGRAFIE / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA

- Holick MF (March 2006). High prevalence of vitamin D inadequacy and implications for health. *Mayo Clinic Proceedings*. 81 (3): 353-73.
- Eriksen EF, Gluerp H (2002). Vitamin D deficiency and aging: implications for general health and osteoporosis. *Biogerontology*. 3 (1-2): 73-7.
- Grant WB, Holick MF (June 2005). Benefits and requirements of vitamin D for optimal health: a review. *Alternative Medicine Review*.10 (2): 94-111.
- Moyad MA. Vitamin D: a rapid review. *DermatolNurs*. 2009; 21:25-30.

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mens. The specimens were correctly identified >99% of the time.
INTER-ASSAY Between-run precision has been determined by 3 independent assays on the same 4 specimens: 10 ng/mL vitamin D, 30 ng/mL vitamin D, 45 ng/mL vitamin D, 100 ng/mL vitamin D standard samples. Three different lots of the VITAMIN D RAPID TEST have been tested using these specimens. The specimens were correctly identified >99% of the time.
SENSITIVITY AND CROSS-REACTIVITY The VITAMIN D RAPID TEST can detect levels of Vitamin D in human whole blood as low as 30 ng/mL. The addition of Vitamin A, B, C, E, K and M showed no cross-reactivity.

PRECAUTIONS

- Read these instructions for use carefully before performing the test. The test is reliable if the instructions are followed correctly (reaction times and how the sample is collected).
- Do not open the sealed aluminum pouch until just before performing the test. The desiccant bag must not be used.
- Humidity and temperature can adversely affect the results. DO NOT FREEZE.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards during all procedures and follow standard procedures for proper disposal of specimens.
- The test is for external and single use only. DO NOT SWALLOW.
- Do not use the kit after the expiry date or if the packaging and/or internal components are damaged.
- It is suggested not to use the buffer beyond 6 months after opening the vial.

CONTENT

- 20 aluminum pouches containing 1 test device for the VITAMIN D TEST and a desiccant bag;
- 20 vials with dropper containing VITAMIN D TEST DILUENT sufficient for 20 tests;
- 20 pipettes for blood sampling;
- 1 Vitamin D Color Card;
- 1 instructions for use.

TEST PROCEDURE

Prepare the necessary material as follows: open the aluminum bag, remove the cassette and discard the desiccant bag. Place the cassette on a flat, dry and clean surface.
Allow the test, specimen, buffer and/or controls to reach room temperature (15–30 °C) prior to testing.
1) Prick the skin of the fingertip with a sterile lancet and wipe away the first sign of blood. Obtain a drop of blood and fill the capillary pipette to the black line.

- Transfer 20 µl of whole blood to the sample well (S) – **fig. A1**. As an alternative to the pipette, 20 µl of venous sample can be dispensed with the laboratory micropipette – **fig. A2**.

- Add 2 drops of diluent into the buffer well (B) waiting about 10 seconds between dispensing one drop and the next. – **fig. B1, B2**.
- Start the timer. Read the results at 10 minutes.

Note: testing should be performed immediately after the whole blood have been collected.

RESULTS INTERPRETATION

READ RESULTS AT 10 MINUTES. DO NOT INTERPRET THE RESULTS AFTER 20 MINUTES.
Please refer to the illustration and compare the T line intensity with Vitamin D Color Card provided with the kit.

DEFICIENT 0-10 (ng/mL) 0-25 (nmol/L)

Two distinct colored lines appear. One is in the Control (C) region and another in the Test (T) region. The line intensity in the test region (T) is equal to or darker than 10 ng/mL line depicted on Vitamin D Color Card provided with the kit.

INSUFFICIENT 10-30 (ng/mL) 25-75 (nmol/L)

Two colored lines appear. One is in the Control (C) region and another in the Test (T) region. The line intensity in the Test (T) region is darker than the 30 ng/mL line and lighter than 10 ng/mL line, both depicted on Vitamin D Color Card provided with the kit.

SUFFICIENT 30-100 (ng/mL) 75-250 (nmol/L)

Two colored lines appear, one line in the Control (C) region and faint colored line appears in the Test (T) region. The line intensity in Test (T) region is equal to or lighter than 30 ng/mL depicted on Vitamin D Color Card.

EXCESS >100 (ng/mL) >250 (nmol/L)

A coloured line appears in the Control (C) region. No line appears in the Test (T) region.

INVALID RESULT

The Control (C) line does not appear. An insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the Control line (invalid result). Review the procedure and repeat the test with a new device and a new sample.

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