Test your health

PRIMA LAB SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
support@primalabs.ch
primalabs.ch

CE **REF** QbD RepS BV
Groenenborgerlaan 16
2610 Wilrijk - Belgium

100066_IFU_11_8.0 11/2024 **REF** 100066-x

Registered trademarks.

INSTRUCTIONS FOR USE
GEBRAUCHSANWEISUNG
NOTICE D'UTILISATION
ISTRUZIONI PER L'USO

VITAMIN D TEST

Rapid self-test for the semi-quantitative detection of vitamin D in human whole blood
Schneller Selbsttest zur semiquantitativen Bestimmung von Vitamin-D in humanem Vollblut

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ENGLISH

VITAMIN D

Vitamin D refers to a group of fat-soluble secosteroids responsible for increasing intestinal absorption of calcium, iron, magnesium, phosphate and zinc. The blood concentration of 25-hydroxy Vitamin D (including D2 and D3) is considered the best indicator of Vitamin D status. Vitamin D deficiency is now recognized as a global epidemic. The health risks associated with Vitamin D deficiency are far more severe than previously thought. Vitamin deficiency has been linked to various serious diseases: osteoporosis, osteomalacia, multiple sclerosis, cardiovascular diseases, pregnancy complications, diabetes, depression, strokes, autoimmune diseases, flu, different cancers, infectious diseases, alzheimer, obesity and higher mortality etc. Therefore, now detecting Vitamin D level is considered as "Medically Necessary Screening Test", and maintaining sufficiently high levels of Vitamin D is necessary to improve bone health, health and well-being of the body.

PRINCIPLE OF THE TEST

The VITAMIN D TEST is an immunoassay based on the principle of competitive binding. During testing, the mixture migrates upward on the membrane chromatographically by capillary action. The membrane is pre-coated with 25 (OH) D antigens on the test line region of the strip.

LIMITATIONS


- 1) The VITAMIN D TEST provides only a semi-quantitative analytical result. A secondary analytical method must be used to confirm the result.
- 2) It is possible that technical or procedural errors, as well as other interfering substances in the whole blood specimen may cause erroneous results.
- 3) As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
- 4) Other clinically available tests are required if questionable results are obtained.

BIBLIOGRAPHY

1. Holick MF (March 2006). „High prevalence of vitamin D inadequacy and implications for health.“ Mayo Clinic Proceedings. 81 (3): 353–73.
2. Eriksen EF, Glerup H (2002). Vitamin D deficiency and aging: implications for general health and osteoporosis. Biogerontology. 3 (1-2): 73-7.
3. Grant WB, Holick MF (June 2005). Benefits and requirements of vitamin D for optimal health: a review. Alternative Medicine Review.10 (2): 94-111.
4. Moyad MA. Vitamin D: a rapid review. DermatolNurs. 2009, 21:25-30.

CONTENT

* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the “content” section on the external box.

PRECAUTIONS

1. For self-testing in vitro diagnostic use only.
2. Do not eat, drink or smoke in the area where the specimens or kits are handled.
3. Store in a dry place at 2-30°C (36-86°F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. DO NOT FREEZE. Do not use after the expiration date.
4. This test kit is intended to be used as a preliminary test only and repeatedly abnormal results should be discussed with doctor or medical professional.
5. Follow the indicated time strictly, using only the specified quantities of blood and diluent.
6. Use the test only once. Do not dismantle and touch the test window of the test cassette.
7. The kit must not be frozen or used after the expiration date printed on the package.
8. Keep out of the reach of children.
9. The test is for external use only. DO NOT SWALLOW.
10. Not recommended for people who take anti-coagulant medications (blood thinners) or people suffering from haemophilia.
11. After using, please dispose of all components according to your local waste disposal laws. Ask your pharmacist for advice.
12. If you test another person, handle the blood samples as if containing a blood-borne infectious agent. Wear disposable gloves during use.

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
Alcohol cleansing gauze
Gaze nettoyante à l'alcool
Garza detergente alcolica

Huaian Tianda Medical Instruments Co., Ltd., No.106 East Songjiang Road, Huaiyin, Economic&Technological Development Zone, 223300 - Huaian City, Jiangsu - P.R. China

Phoenix Innovative Healthcare Manufacturers Pvt. Ltd., EL-209, Shil Mahape Road, Electronic Zone, MIDC, TTC Industrial Area, Mahape, Navi Mumbai 400 710 MH 1 India

CE REF

Riomavix Sociedad Limitada
Calle de Almansa 55, 10
Madrid 28039, Spain

CE REF

Advena Ltd., Tower Business Centre, 2nd Flr, Tower Street, Swatar, BKR 4013, Malta

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

DIRECTIONS FOR USE

- Wash hands with soap and warm water, rinse with clean water and allow to dry. If this is not possible, use the provided gauze as an alternative.
Note: The use of warm water facilitates capillary blood collection as it induces vasodilation.
- Prepare the necessary material as follows: open the aluminium pouch, take out only the test cassette and throw away the desiccant bag. Open the plastic packet containing the pipette. **–FIG. A**
- Carefully rotate the protective cap of the sterile lancet 360° without pulling it. Extract and discard the released cap. **–FIG. B**
- Carefully massage the finger chosen for the puncture (the side of the ring finger is recommended). It is important that the massage is done from the palm of the hand to the phalanx, to improve blood flow. Press the open end of the lancet (the side the cap has been extracted from), against the fingertip **–FIG. C** The tip of the lancet automatically retracts after use. If the lancet does not work properly, discard it and use the second one supplied. If the second one is not required, it can be disposed of without special precautions.
- Holding the hand down, massage the finger until a large drop of blood forms. It is important to massage from the palm of the hand to the phalanx to improve blood flow. **–FIG. D**
- Take the pipette **without pressing the bulb**. Two sampling methods are suggested:
–FIG. E1: hold the pipette horizontally **without pressing the bulb** place it in contact with the drop of blood, it will enter the pipette by capillarity. Move the pipette away when the black line is reached. If there is not enough blood, continue to massage the finger until the black line is reached.
–FIG. E2: place the pipette on a clean, flat surface with the tip protruding from the shelf, then place the drop of blood in contact with the pipette, it will enter by capillarity. If blood is not sufficient, continue massaging the finger until the blood has reached the black line.
Avoid, as far as possible, to continuously move the tip of the pipette away from the finger in order to prevent the formation of air bubbles.
- Deposit the blood collected into the Specimen well (S) of the cassette, by squeezing the dropper bulb. **–FIG. F**
- Wait for the blood to be totally dispensed in the well. Unscrew the cap of the buffer bottle and add 2 drops of buffer into the Buffer well (B) of the cassette and start a timer. **–FIG. G**
- Wait for the colored line(s) to appear. **Read results at 10 minutes.** Compare the T line intensity with “Vitamin D Color card” provided with the kit to get the vitamin D level in your blood. Do not interpret the result after 20 minutes.

RESULTS INTERPRETATION

Please refer to the illustration and compare the T line intensity with “Vitamin D Color card” provided with the kit. As control system, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

DEFICIENT 0-10 (ng/mL) 0-25 (nmol/L)

Two distinct colored lines appear. One is in the control region (C) and another should be in the test region (T). The line intensity in the test region (T) is equal to or darker than 10 ng/mL line depicted on Vitamin D color card provided with the kit.

INSUFFICIENT 10-30 (ng/mL) 25-75 (nmol/L)

Two colored lines appear. One is in the control region (C) and another should be in the test region (T). The line intensity in the test region (T) is darker than the 30 ng/mL line depicted on the Vitamin D color card provided with the kit and lighter than 10 ng/mL line depicted on Vitamin D color card provided with the kit.

SUFFICIENT 30-100 (ng/mL) 75-250 (nmol/L)

Two colored lines appear. One line should be always in the control region (C) and faint colored line appears in the test region (T). The line intensity in region (T) is Equal or lighter than 30 ng/mL depicted on Vitamin D Color Card.

EXCESS >100 (ng/mL) >250 (nmol/L)

One colored line appears in the control line region (C). No apparent colored line appears in the test line region (T). If the result is excess, it is recommended to consult a physician.


INVALID

Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

FA.Q. – QUESTIONS AND ANSWERS

HOW DOES THE VITAMIN D TEST WORK? In medicine, a 25-hydroxy vitamin D is the main storage form of vitamin D in the body. Therefore, the overall status of vitamin D can be determined by detecting the content of 25-hydroxy vitamin D. 25-hydroxy vitamin D level less than 30 ng/mL in case of a positive result, indicates vitamin D Deficiency or Insufficiency. vitamin D supplements can be recommended in these cases.
WHEN SHOULD THE TEST BE USED? The clinical application of vitamin D is mainly for diagnosis, treatment and monitoring of rickets (children), osteomalacia, postmenopausal osteoporosis and renal osteopathy. Vitamin D deficiency is also associated with many other diseases, including cancer, cardiovascular disease, autoimmune diseases, diabetes and depression. Monitor your vitamin D levels to determine whether to take vitamin D supplements. The Vitamin D Rapid Test can be used any time of the day.
CAN THE RESULT BE INCORRECT? The results are accurate as far as the instructions are carefully respected. Nevertheless, the result can be incorrect if the VITAMIN D TEST cassette gets wet before test performing or if the quantity of blood dispensed in the sample well is not sufficient, or if the number of buffer drops are less than 2 or more than 3. The capillary dropper provided in the box allows making sure the collected blood volume is correct. Besides, due to immunological principles involved, there exist the chances of false results in rare cases. A consultation with the doctor is always

recommended for such tests based on immunological principles.
HOW TO INTERPRET THE TEST IF THE COLOR AND THE INTENSITY OF THE LINES ARE DIFFERENT? Please refer to the illustration and compare the T line intensity with “Vitamin D Color card” provided with the kit.
IF I READ THE RESULT AFTER 20 MINUTES, WILL THE RESULT BE RELIABLE? No. The result should be read at 10 minutes after adding the buffer. The result is unreliable after 20 minutes.
WHAT DO I HAVE TO DO IF THE RESULT IS DEFICIENT OR INSUFFICIENT? If the result is deficient or insufficient, it means that the vitamin D level in blood is less than 30 ng/mL and that you should consult a physician to show the test result. Then, the physician will decide whether additional analysis should be performed.
WHAT DO I HAVE TO DO IF THE RESULT IS SUFFICIENT? If the result is sufficient, it means that the vitamin D level is higher than or equal to 30 ng/mL and is within the normal range. A case of vitamin D toxicity (hypercalcemia), though rare, cannot be excluded based on such test results. However, if the symptoms persist, it is recommended to consult a physician.
HOW ACCURATE IS VITAMIN D TEST? The test is very accurate. Accuracy was calculated through a clinical performance study on a total of 90 samples. Evaluation reports show a total concordance rate of 94.4% with reference methods (95% CI: 87.65 -97.60).

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