Test your health

PRIMA LAB SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
support@primalabsa.ch
primalabsa.ch

EC REF QbD RepS BV
Groenenborgerlaan 16
2610 Wilrijk - Belgium

100082_IFU_11_8.0 12/2024

REF 100082-x

Registered trademarks.

INSTRUCTIONS FOR USE
GEBRAUCHSANWEISUNG
NOTICE D'UTILISATION
ISTRUZIONI PER L'USO

SPERM TEST

Rapid self-test for in vitro qualitative estimation
of sperm concentration in human semen

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PRIMA
HOME TEST

PRINCIPLES OF THE TEST

SPERM TEST is biochemical assay for in vitro qualitative estimation of sperm concentration in human semen as an auxiliary aid in clinical diagnosis of the infertility and/ or pregnancy planning by self-evaluation of sperm concentration above or below the required concentration for successful pregnancy.

This product uses the inert glass fiber membrane with high water absorption and the pore size of less than 0.5µm to filter the semen. Sperm cells are trapped on the surface of the first layer of the membrane, and staining solution is used that can dye sperm cells. The darker the colour of Well A, the higher the sperm concentration. If the colour of test well A is lighter than the standard colour of reference Well B, it means that the concentration of sperm is less than 15 million/ml. If the colour of test well A is darker than the standard colour of reference Well B, it means that the sperm concentration is greater than 15 million/ml. Sperm concentration of 15 million/ml is the minimum expected sperm concentration level for pregnancy.

This kit is designed to be used for in vitro qualitative estimation of the sperm concentration of human semen. Essentially this means the test will determine if the number of sperms is at an adequate level for conception to occur with sexual intercourse, subject to female partner's ovulation in time. A low sperm concentration would indicate less likelihood of conception. It would be advisable to see your medical professional who can advise what can be done to improve the sperm concentration.

FA.Q. – QUESTIONS AND ANSWERS

WHY THE TIME IS 15 MINUTES BEFORE THE SEMEN BEING TAKEN OUT FROM THE SEMEN COLLECTING CUP FOR TEST, AND THE STORAGE TIME AFTER SAMPLING DOES NOT EXCEED 12 HOURS? The fresh semen is viscous, and normal semen needs to be incubated for 30-60 minutes at 37 °C to liquefy completely. Only in liquid state, the semen can be used for test, because the viscous semen cannot completely pass the membrane of the test well. **The sheet at bottom of the collection cup can make the semen liquefy quickly within 15 minutes.** If the storage time of sperm specimen is too long, it may cause lysis of sperms, which may affect accuracy of the results.

HOW LONG DO THE SEMEN AND STAINING SOLUTION NEED TO BE IN THE TEST WELL? In general, they may pass membrane of the test well within several seconds, if they cannot pass the membrane completely after 5 minutes, and which shows that the semen has not liquefied completely, or density of the sperm is too high, and we must repeat the test. The reasons of semen non-liquefaction may be that the enzyme in the collection cup becomes invalid or the user does not use the semen cup correctly or there is some other cup used for semen collection.

IS THE SOLUTIONS USED FOR TEST SAFE? The solution is safe, and which is a synthetic dye, intake of which in concentration less than 5mg/kg does not harm to human body. The concentration of the staining fluid is less than 10ug/ml.

CAN THE ABNORMAL RESULTS SHOW THAT THE SUBJECT HAS NO ABILITY TO HAVE CHILDREN? Sperm concentration is one of the several semen analysis tests. There are other factors that should be taken into account such as motility. Therefore, a medical expert's opinion in such cases would be highly recommended.

WHAT IS THE REASON THAT MAY RESULT IN WRONG TEST RESULTS? Any mistake at any point of time from sample collection to test timing to non-compliance to abstinence may result in erroneous test results.

HOW ACCURATE IS THE SPERM TEST? The Test is very accurate. Assessment reports show a concordance rate above 99% (95%: 97.0-99.9%) with reference methods.

PRECAUTIONS

Please read all the information in this package insert before performing the test. This kit can only be used as an in vitro diagnostic test using human semen as specimen and cannot be used with specimens of other body fluids. This test kit is intended to be used as a preliminary test only and repeatedly negative results should be discussed with doctor or medical professional. When adding specimens, staining solution and washing solution, try to avoid any bubbles as this could adversely affect the test results. Make sure you correctly follow the "time" instructions when carrying out the test and observing the results. Do NOT remove the liquifying powder from the bottom of the collection cup since it will help the semen to liquefy. The kit should be stored at room temperature, avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use.

Once the test card's package is opened, it should be used as soon as possible, to avoid being exposed to the air for long periods which could result in the test not working correctly. The kit must not be frozen or used after the expiry date printed on the outer foil. The test kit should be stored at room temperature or refrigerated (2°-30°) in the sealed pouch to the date of expiration. DO NOT FREEZE. The test kits should be kept away from direct sunlight, moisture and heat. After use, all components can be disposed according to local regulations; all specimen should be considered potentially hazardous and handled in the same manner as an infectious agent. In case of difficulties in color identification (such as Daltonism), ask for help in test reading. External use only. Do not swallow.

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BIBLIOGRAPHY / BIBLIOGRAPHIE / BIBLIOGRAPHIE / BIBLIOGRAPHIE

- Leaver R.B.:Male infertility: an overview of causes and treatment options. Br J Nurs. 2016 Oct 13;25(18):S35-S40.
- Sikka S.C., Hellstrom W.J.: Current updates on laboratory techniques for the diagnosis of male reproductive failure. Asian J Androl. 2016 May-Jun;18(3):392-401
- Sharma R., Harlev A., Agarwal A., Esteves S.C.: Cigarette Smoking and Semen Quality: A New Meta-analysis Examining the Effect of the 2010 World Health Organization Laboratory Methods for the Examination of Human Semen. Eur Urol. 2016 Oct;70(4):635-645.
- Centola G.M.: Semen assessment. Urol Clin North Am. 2014 Feb;41(1):163-7. doi: 10.1016/j.ucl.2013.08.007. Epub 2013 Sep 13.
- Cooper T.G., Noonan E., von Eckardstein S., Auger J., Baker H.W., Behre H.M., Haugen T.B., Kruger T., Wang C., Mbizvo M.T., Vogelsson K.M.: World Health Organization reference values for human semen characteristics. Hum. Reprod Update. 2010 May-Jun;16(3):231-45.

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ENGLISH KIT CONTENT

Materials Provided for 1 test

- 1 Aluminium pouch containing:
 - 1 pipette
 - 1desiccant bag
 - 1 SPERM TEST cassette
- 1 vial with dropper tip containing 3 mL of SPERM TEST staining solution
- 1 collection cup
- 1 vial with dropper tip containing 3 mL of SPERM TEST washing solution
- 1 Instructions for use

Materials Required But Not Provided: Timer

SPECIMEN COLLECTION AND PREPARATION

Before testing, read the instructions carefully and completely

- 1) Before testing, it is important that you refrain from any sexual activity for 3-7 days. This ensures that the volume and quality of sperm is at its peak and the test will then be an accurate determination of sperm concentration.
- 2) Using masturbation, the semen should be collected directly into the collection cup. **ATTENTION: Do not remove the liquifying powder from the cup since it will help the semen to liquefy.**
- 3) Care should be taken that collected semen is not contaminated by touch of hands or tissues or any other materials.
- 4) Shake the semen evenly in the collection cup and leave it to stand for 15 minutes at room temperature until the semen liquefy. Do not use semen stored for more than 12 hours after liquefaction.

TEST PROCEDURE:

- 1) Remove the test cassette from the foil pouch and lay it horizontally on a flat surface. Using the pipette provided in the foil pouch, dispense one drop of semen into sample well A. **(FIG. A-B)**
- 2) Once the semen is soaked into Well A, add three drops of the blue staining solution to Test Well A. Let it soak for 1-2 minutes. **(FIG. C)**
- 3) Now add two drops of the transparent washing solution to Test Well A **(FIG. D)**, and let it soak for 1-2 minutes, and then read the results immediately.
- 4) Read the colour of Test Well A, comparing the colour of Test Well A to Reference Well B. The darker the colour of Well A, the higher the sperm concentration.


RESULTS INTERPRETATION

Compare the colour of test well A to reference Well B

NORMAL


The colour of test well A is darker than or the same as the standard colour of reference Well B. It means that the sperm concentration is greater than or equivalent to 15 million/ml. The likelihood of conception is high with this sperm concentration, subject to other conditions, such as ovulation time of the partner being favorable.

ABNORMAL

The colour of test well A is lighter than the standard colour of reference Well B. It means that the sperm concentration is less than 15million/ml. This is known as oligospermia (a range that is normally between 5million/ml and 15million/ml). The likelihood of conception is less with this sperm concentration and further medical consultation is recommended.

NOTE: If well A is colourless, it means the sperm concentration is less than 5 million/ml or zero. This condition is known as azoospermia. If you are unsure of the result or you feel the result is abnormal you should repeat the test using the second test that is included in the pack but make sure you do not ejaculate through any sexual activity for 6 days before carrying out the second test. If the second test is also abnormal, you should discuss the results with your doctor or medical professional.

LIMITATIONS

- 1) For in vitro qualitative estimation of sperm concentration in human semen.
- 2) Sperm concentration is just one of the important tests for fertility. But other tests of semen like motility, volume and morphology as well as ovulation in females are also important. For the cases of infertility, it is recommended that other tests are also taken in consideration.

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