PRIMA Lab SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
support@primalabsa.ch
primalabsa.ch

EC REP QbD RepS BV - Groenenborgerlaan 16
2610 Wilrijk - Belgium

REF 800016IP-20P 800016IP-20P_IFU_90_80_11/2024

Your partner for rapid testing

IVD for professional use only
IVD nur für den professionellen Gebrauch
DIV à usage professionnel uniquement
IVD solo per uso professionale
IVD sólo para uso profesional
IVD apenas para uso profesional
IVD wyłącznie do użytku profesjonalnego


Table 1. Performance / Tabelle 1. Leistung / Tableau 1. Performances / Tabella 1. Prestazioni / Tabla 1. Rendimiento / Tabela 1. Desempenho / Tabela 1. Oziatanie.

	LABORATORY REFERENCE METHOD	
HELICOBACTER PYLORI ANTIGEN RAPID TEST	Positive	Negative
Positive	40	2
Negative	2	56

Specificity = 96.6 % (CI 95%: 91.8 - 100%)
Sensitivity = 95.2% (CI 95%: 88.8 - 100%)
Accuracy = 96% (CI 95%: 92.1 - 99.8%)

REFERENCES / BIBLIOGRAFIE / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA

- Shimajima T. Stool antigen tests for the management of Helicobacter pylori infection. World Journal of Gastroenterology. 2013; 19(45):8188-91.
- Ana Isabel Lopes, Filipa F Vale, Mónica Oleastro. Helicobacter pylori infection - recent developments in diagnosis. World Journal of Gastroenterology. 2014 July 28; 20(28): 9299-9313.
- Andreas Mentis, Philippe Lehours, and Francis Mégraud. Epidemiology and Diagnosis of Helicobacter pylori infection. Helicobacter. 2015 Sep; 20 Suppl 1:1-7.
- Zagari RM, Rabitti S, Eusebi LH, Bazzoli F. Treatment of Helicobacter pylori infection: A clinical practice update. Eur J Clin Invest. 2018 Jan;48(1). doi: 10.1111/eci.12857.

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CONTENT

- 20 hermetically sealed aluminum bag containing 1 HELICOBACTER PYLORI ANTIGEN TEST cassette and 1 desiccant bag;
 - 20 vials with red screw cap, with collection stick and dropper, containing enough HELICOBACTER PYLORI ANTIGEN TEST diluent for 20 tests;
 - 1 instructions for use leaflet.
- Material required but not supplied: a device for measuring time (e.g. stopwatch, clock).

PERFORMANCE CHARACTERISTICS

To assess the performance of HELICOBACTER PYLORI ANTIGEN RAPID TEST, a clinical evaluation was performed by comparing the results obtained against Helicobacter pylori Test ELISA. 100 samples were tested with both methods (Table 1).
Specificity = 96.6 % (CI 95%: 91.8 - 100%)
Sensitivity = 95.2% (CI 95%: 88.8 - 100%)
Accuracy = 96% (CI 95%: 92.1 - 99.8%)
MEASURING RANGE The study concerning the Limit Of Detection (LOD) was conducted with the Helicobacter Antigen Biospecific Inc. The LOD of HELICOBACTER PYLORI ANTIGEN RAPID TEST is 10 ng/mL. No adverse effect on T-line formation (prozone effect) was observed for antigen value up to 340.000 ng/mL.
CROSS-REACTIVITY The following compounds were tested with HELICOBACTER PYLORI ANTIGEN RAPID TEST and no interference was observed:

Table 2. Cross-reactants tested.

SUBSTANCE	CONCENTRATION ASSAYED
<i>Streptococcus Pyogenes</i>	1.32 x 10 ⁹ CFU/mL
<i>Streptococcus Agalactiae</i>	1.48 x 10 ⁹ CFU/mL
<i>Streptococcus Dysgalactiae</i>	2 x 10 ⁹ CFU/mL
<i>Pseudomonas Aeruginosa</i>	5.77 x 10 ⁹ CFU/mL
<i>Staphylococcus Aureus</i>	1.89 x 10 ⁹ CFU/mL
<i>Proteus Vulgaris</i>	4.88 x 10 ⁹ CFU/mL
<i>Corynebacterium Diphtheriae</i>	1.49 x 10 ⁹ CFU/mL
<i>Yersinia Enterocolitica</i>	1.35 x 10 ⁹ CFU/mL
<i>Escherichia coli</i>	1.43 x 10 ⁹ CFU/mL
<i>Campylobacter Jejuni</i>	6.8 x 10 ⁹ CFU/mL
<i>Enterobacter Cloacae</i>	8.57 x 10 ⁹ CFU/mL
<i>Klebsiella Pneumoniae</i>	1.81 x 10 ⁹ CFU/mL
<i>Shigella Flexneri</i>	1.3 x 10 ⁹ CFU/mL
<i>Clostridium Perfringens</i>	7.7 x 10 ⁹ CFU/mL
<i>Candida Albicans</i>	1.53 x 10 ⁹ CFU/mL

SAMPLE COLLECTION AND PREPARATION

- Allow the test, sample and buffer to reach room temperature (15-30°C) before testing.
- Unscrew the top of the sample collection device containing the sample collection probe. Note: Avoid dilution of stool specimens with urine or water from the toilet bowl. Perform the test as soon as possible after sample collection.
 - Collect the stool sample with the top of the collection device by plunging it into 3 different places of the same stool sample –**fig. A** and put it in the collection vial.
Do NOT put the collection stick back into the tube in between collecting the 3 samples.
 - Screw on the cap again, immersing the stick in the stool extraction liquid, and shake the vial for at least 10 seconds. –**fig. B**
 - Proceed with test procedure.

TEST PROCEDURE

- Open the aluminium bag, remove the test cassette and place it on a clean and level surface. Throw away the desiccant bag.
- Break the end of the red cap, exposing the dropper portion. –**fig. C**
- Dispense 3 drops of the diluted stool sample in the sample well (S) shown on the test cassette. Wait 1-2 seconds between dispensing one drop of diluted stool sample. –**fig. D**
- Wait 10 minutes.

RESULTS INTERPRETATION

READ THE RESULTS AT EXACTLY 10 MINUTES.
The intensity of the line colors is not relevant for the purposes of interpretation of the Test's results.

POSITIVE RESULT

Two coloured lines appear in the reading window by the Test (T) and Control (C) signs. The T line may be less intense (lighter) than the C one.
This result means there is *Helicobacter pylori* Antigen in the faeces.

NEGATIVE RESULT

A coloured line appears in the Control (C) region. No line appears in the Test (T) region. This result means that there is no *Helicobacter pylori* Antigen in the faeces or that its concentration is below the test's detection limit (10 ng/mL).

INVALID RESULT

The Control (C) line does not appear. An insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the Control line (invalid result). Review the procedure and repeat the test with a new device and a new sample.

QUALITY CONTROL

An internal procedural control is included in the test cassette. A coloured line appearing in the Control line region (C) is considered an internal positive procedural control, confirming sufficient specimen volume and correct procedural technique.

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