
PRIMA Lab SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
support@primlabsa.ch
primlabsa.ch

EC REF QbD RepS BV - Groenenborgerlaan 16
2610 Wilrijk - Belgium

REF 800025IP-20P 800025IP-20P_IFU_90_80 11/2024

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- MacKenzie E.L., Iwasaki K., Tsuji Y. (2008) Intracellular Iron Transport and Storage: From Molecular Mechanisms to Health Implications. Antioxid Redox Signal. 2008 Jun; 30(6): 997-1030.
- Chiancone E., Stefanini S. and Antonini E. (1980), "Ferritin: structural and functional aspects in Radioimmunoassay of Hormones, Proteins and Enzymes". Proc. Int. Symposium. Excerpta Medica Amsterdam : 197-203.
- Cullis J.O, Fitzsimons E.J, Griffiths WJ, Tsochatzis E, Thomas DW. British Society for Haematology (2018). "Investigation and management of a raised serum ferritin. Br J Haematol. 181(3):331-340.
- Simon T.L., Garry P.J., Holper E.M. (1981) "Iron stores in blood donors". JAMA 245 n°20 : 2038-2043.
- Morse E.E., Cable R., Pisciotto P., Kakaya R., Kiraly T. (1987). "Evaluation of iron status in women identified by copper sulfate screening as ineligible to donate blood". Transfusion 27 n°3 : 238-241.
- Jacob RA, Sandstead HH, Klevay LM, Johnson LK (1980). "Utility of serum ferritin as a measure of iron deficiency in normal males undergoing repetitive phlebotomy." Blood.56(5):786-91.
- Coenen J.L.L.M., Van Dieijen-Visser M.P., Van Pelt J., et al. (1991) "Measurement of serum Ferritin used to predict concentrations of iron in bone marrow in anemia of chronic disease". Clin. Chem. 37 : 560-563.

Your partner for rapid testing

IVD for professional use only
IVD nur für den professionellen Gebrauch
DIV à usage professionnel uniquement
IVD solo per uso professionale
IVD sólo para uso profesional
IVD apenas para uso profissional
IVD wyłącznie do użytku profesjonalnego

- Witte D.L., Dick F.R., Goeken J. et al. (1985) "C-reactive protein (CRP) aids interpretation of serum Ferritin (FRTN) (Abstract). Clin. Chem. 31 : 1011.
- Witte D.L., Angstadt D.S., Davis S.H., Schrantz R.D. (1988) "Predicting bone marrow iron stores in anemic patients in a community hospital using Ferritin and erythrocyte sedimentation rate. Am. J. Clin. Pathol. 90 : 85- 87.
- Alan E. Mast, Morey A. Blinder, Ann M. Granowski, Cara Chumley, Mitchell G. Scott. (1998) "Clinical utility of the soluble transferrin receptor and comparison with serum ferritin in several populations". The American Association for Clinical Chemistry.
- Weiss G., Goodnough L.T. (2005). "Anemia of chronic disease. N Engl J Med. 352(10):1011-23.
- Munoz M. et al. 2017. Current misconceptions in diagnosis and management of iron deficiency. Blood Transfus. Sep.15(5):422-437.

Table 1. Performance / Tabelle 1. Leistung / Tableau 1. Performances / Tabella 1. Prestazioni / Tabla 1. Rendimiento / Tabla 1. Desempenho / Tabla 1. Zdatanie

IRON FER RAPID TEST	Laboratory Reference Method	
	Positive	Negative
Positive	20	0
Negative	0	60

Specificity = 100 % (CI 95%: 93.98 - 100%)
Sensitivity = 100 % (CI 95%: 83.89 - 100%)
Accuracy = 100 % (CI 95%: 95.42 - 100%)

Table 2. Interfering substances / Tabelle 2. Beeinträchtigende Substanzen / Tableau 2. Substances interférentes / Tabella 2. Sostanze interferenti / Tabla 2. Sustancias que interfieren / Tabla 2. Substâncias interferentes / Tabla 2. Substancje zakłócające

SUBSTANCE	CONCENTRATION ASSAYED
Bilirubin	4.74 mg/dL (81.0 mmol/L)
Cholesterol	1248 mg/dL(32.25 mmol/L) 480 mg/dL (12.4 mmol/L)
Triglycerides	443 mg/dL (5 mmol/L) 868 mg/dL (9.8 mmol/L)
Testosterone	1200 ng/dL
EDTA	up to 10 mg/mL
Vacutainer K2EDTA (ready to use tube for blood collection)	Original tube
*HAMA Positive (Plasma)	up to 181 ng/mL
*HAMA Positive (Serum)	up to 504 ng/mL
Isopropanol 70%	up to / hasta
Hand Creams	Test in which cream is put onto the finger skin (as per claimed recommended use)

EN - INSTRUCTIONS FOR USE

IRON AND ANEMIA

Ferritin is an ubiquitous protein whose role is to store Iron, the essential element for oxygen transportation and delivery in the blood. Its structure consists of 24 subunits that form a spherical cavity in which Fe(II) is sequestered and accumulated in the cells. Although the highest concentrations of ferritin are found in liver, spleen and bone marrow, a small fraction of body ferritin is also found circulating in serum. Clinical studies have shown that serum concentration of ferritin correlates with the body's iron storage and its dosage can help to diagnose anemia: in particular, 1 ng/mL reduction in serum ferritin corresponds to an approximated loss of 5 mg of iron stored in the tissues. In case of increased body requirement, iron stores are slowly mobilized even before anemia becomes overt and ferritin levels decrease accordingly; therefore, detection of low serum ferritin level is useful to indirectly diagnose the progressive depletion of iron stores in the body. In particular, serum ferritin levels below 30 ng/mL, which is the cut-off of IRON FER RAPID

IRON FER RAPID TEST

Rapid immuno-chromatographic test for the detection of ferritin in whole blood or serum samples
Immunochromatographischer Schnelltest zum Nachweis von Ferritin in Vollblut- oder Serumproben
Test immunochromatografico rapido per la rilevazione della ferritina in campioni di sangue intero o di siero
Prueba immunocromatográfica rápida para la detección de ferritina en muestras de sangre total o de suero
Teste imunocromatográfico rápido para detecção de ferritina em amostras de sangue total ou soro
Szybki test immunochromatograficzny do wykrywania ferrytyny w próbkach krwi pełnej lub surowicy

TEST, correlate with a significant deficit in body iron stores.

PRINCIPLE OF THE TEST

IRON FER RAPID TEST is a lateral flow test device based upon a principle of solid-phase immunochromatography (also called lateral flow technology) employing a competitive technique. The detection of ferritin takes place on a test strip embedded in a plastic cassette. The methodology implied is a direct sandwich EIA test also named G.A.I.A. (Gold adsorbed Immuno Assay), in which the target analyte (ferritin) is recognized by a specific antibody (anti-ferritin), conjugated with gold particles, from one side and by a specific monoclonal antibody anti-ferritin present in the region of the Test line (T line). In case of the levels of ferritin in the blood are above 30 ng/mL, all the ferritin present in the sample is firstly bound by the monoclonal antibody anti-ferritin on the gold nanoparticles creating the immunocomplexes, which are selectively captured by the Test line (where another monoclonal anti ferritin is striped), thus showing a red-purple band in the region of the T line, indicating a negative result. Consequently, two bands are visible (T line – C line), showing a negative result. In case of ferritin level lower than the cut-off value (30 ng/mL) only one band is visible in the C line region, showing a positive result.

PRECAUTIONS

- Read these instructions for use carefully before performing the test. The test is reliable if the instructions are followed correctly (reaction times and how the sample is collected).
- Do not open the sealed aluminum pouch until just before performing the test. The desiccant bag must not be used.
- Humidity and temperature can adversely affect the results. DO NOT FREEZE.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards during all procedures and follow standard procedures for proper disposal of specimens.
- The test is for external and single use only. DO NOT SWALLOW.
- Do not use the kit after the expiry date or if the packaging and/or internal components are damaged.
- Store the buffer correctly closed with its cap.

LIMITATIONS

- As with any diagnostic method, the physician must evaluate the test result in light of other clinical evidence.
- People suffering from haemophilia or taking anti-coagulants could have an excessive bleeding after the fingertip prick.
- Despite the reliability of the test, false positive or negative results can not be excluded.
- Certain serum samples with a high concentration of Rheumatoid Factor may lead to inappropriate positive results. These cases should be identified before testing.

CONTENT

- 20 hermetically sealed aluminum bags containing 1 IRON FER TEST cassette and 1 desiccant bag;
 - 1 vial with dropper tip containing sufficient IRON FER TEST diluent for 20 tests;
 - 1 instructions for use leaflet.
- Material required but not supplied: a device for measuring time (e.g. stopwatch, clock).

TEST PROCEDURE

Allow the kit components and/or samples to reach room temperature (15-30 °C) before testing. Prepare the necessary material as follows: open the aluminum bag, remove the cassette and discard the desiccant bag. Place the cassette on a flat, dry and clean surface.

WHOLE BLOOD SAMPLING PROCEDURE

- Prick the skin of the fingertip with a sterile lancet and wipe away the first sign of blood.

Obtain a drop of blood and fill the capillary pipette to the black line.

- Transfer 30 µl of whole blood to the sample well (S) – **fig. A1**. As an alternative to the pipette, 30 µl of venous sample can be dispensed with the laboratory micropipette – **fig. A2**.

- Add 2 drops of diluent waiting about 5 seconds between dispensing one drop and the next. – **fig. B1, B2**.

- Start the timer. Read the results at 5 minutes.

SERUM COLLECTION PROCEDURE

- Using a laboratory micropipette, collect 17 µl of serum and transfer to the sample well (S) – **fig. A3**.
- Add 2 drops of diluent waiting about 5 seconds between dispensing one drop and the next – **fig. B3**.

- Start the timer. Read the results at 5 minutes.

RESULTS INTERPRETATION

READ THE RESULTS AT EXACTLY 5 MINUTES. The intensity of the line colors is not relevant for the purposes of interpretation of the Test's results.

POSITIVE RESULT

A coloured line appears in the Control (C) region. No line appears in the Test (T) region. This result means that the value of ferritin is below normal levels.

NEGATIVE RESULT

Two coloured lines appear in the reading window by the Test (T) and Control (C) signs. The T line may be less intense (lighter) than the C one. This result means the levels of ferritin in the blood are above 30 ng/mL. In some cases, samples with ferritin levels just below this value (between 27 ng/mL and 29 ng/mL) can yield a result of this type.

INVALID RESULT

The Control (C) line does not appear. An insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the Control line (invalid result). Review the procedure and repeat the test with a new device and a new sample.

PERFORMANCE CHARACTERISTICS

To assess the performance of IRON FER RAPID TEST, a clinical evaluation was performed by comparing the results obtained against the laboratory's reference method. 80 samples were analyzed. Samples with a concentration of ferritin equal or lower than 30 ng/mL were considered as positive whilst samples with ferritin concentration equal or higher than 31 were considered negative (Table 1).

Specificity = 100 % (CI 95%: 93.98 - 100%)
Sensitivity = 100 % (CI 95%: 83.89 - 100%)
Accuracy = 100 % (CI 95%: 95.42 - 100%)

LIMIT OF DETECTION

The study concerning the Limit Of Detection (LOD) was conducted with the NIBSC standard 94/572. The LOD of IRON FER RAPID TEST is 30 ng/mL of ferritin.

INTERFERING SUBSTANCES

The following compounds were tested with IRON FER RAPID TEST and no interference was observed (Table 2).

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