PRIMA Lab SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
support@primalabsa.ch
primalabsa.ch

EC REP QbD RepS BV - Groenenborgerlaan 16
2610 Wilrijk - Belgium

REF 800040IP-20P 800040IP-20P_IFU_90_80 11/2024


Table 1. Performance / Tabelle 1. Leistung / Tableau 1. Performances / Tabela 1. Desempenho / Tabela 1. Działanie

LABORATORY REFERENCE METHOD		
BOWEL FOB RAPID TEST	Positive	Negative
Positive	37	1
Negative	10	100

Specificity = 99 % (CI 95%: 97% – 100%)
Sensitivity = 78.7 % (CI 95%: 67% – 90.4%)
Accuracy = 92.6% (CI 95%: 88.4 – 96.8%)

REFERENCES / BIBLIOGRAFIE / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA

- Van Rossum Leo G. et al.: "Random Comparison of Guaiac and Immunochemical Fecal Occult Blood Tests for Colorectal Cancer in a Screening Population" *Gastroenterology*, Volume 135, Issue 1, 82 – 90
- Faivre J., Dancoart V., Denis B., Darval E., Plette C., Perrin P., Bidon J.M., Jurd C., Jung S., Levillain R., Viguer J., Bretagne J.F.: "Comparison between a guaiac and three immunochemical faecal occult blood tests in screening for colorectal cancer". *European Journal of Cancer*, 2012; 48(16):2369-76.
- Benton S.C., Seaman H.E., Halloran S.P.: "Faecal Occult Blood Testing for Colorectal Cancer Screening: the past or the Future". *Current Gastroenterology Reports*, 2015; 17(2):428.
- European Commission. European guidelines for quality assurance in colorectal cancer screening and diagnosis. First edition 2010.

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EN - INSTRUCTIONS FOR USE

FOB (FECAL OCCULT BLOOD)

Faecal occult blood is the presence in the faeces (stool) of blood that is not clearly visible to the naked eye. Occult blood may be the consequence of gastro-intestinal disorders such as ulcers, polyps, colitis, diverticulitis, rhagades, haemorrhoids or cancer. The signs of these lesions are often silent in the first phases, so searching for occult blood in the stool is an important screening test, which allows to identify the presence of pathologies in the gastro-intestinal system at an early stage.

PRINCIPLE OF THE TEST

BOWEL FOB RAPID TEST is a lateral flow test device based upon a principle of solid-phase immunochromatography (also called lateral flow technology). The detection of fecal occult blood takes place on a test strip embedded in a plastic cassette. The methodology implied is a direct sandwich EIA test, in which the target analyte (hemoglobin) is recognized by specific antibodies, one which is coated to colloidal gold particles and the other one is sprayed onto the membrane. In case of presence of hemoglobin in a patient's feces, the immunocomplexes are then captured by the Test line, thus showing a red-purple band appearance in correspondence of the T line, indicating a positive result. Consequently, two bands visible (T line – C line), after 5 minutes from the diluent drops application, showing a positive result. In case of absence of hemoglobin or with hemoglobin level lower than the cut-off value (40 ng/mL) only one band visible in the C line region, showing a negative result.

PRECAUTIONS

- Read these instructions for use carefully before performing the test. The test is reliable if the instructions are followed correctly (reaction times and how the sample is collected).
- Do not open the sealed aluminum pouch until just before performing the test. The desiccant bag must not be used.
- Humidity and temperature can adversely affect the results. DO NOT FREEZE.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards during all procedures and follow standard procedures for proper disposal of specimens.
- The test is for external and single use only. DO NOT SWALLOW.
- Do not use the kit after the expiry date or if the packaging and/or internal components are damaged.
- Store the buffer correctly closed with its cap.
- Avoid cross-contamination of samples by using a new specimen collection tube for each sample.

LIMITATIONS

- As with any diagnostic method, the physician must evaluate the test result in light of other clinical evidence.
- Despite the reliability of the test, false positive or negative results can not be excluded.
- The presence of blood in stools may be due to several causes, besides colorectal bleeding such as haemorrhoids, blood in urine or stomach irritations. Bleeding from the upper part of the digestive tract (for example in the case of stomach or duodenal ulcers) may not be detected all of the time, due to further protein digestion and the difficulty of antibodies to recognise the haemoglobin antigen after proteolysis.
- All colorectal bleeding may not be due to pre-cancerous or cancerous polyps.
- Negative results do not exclude bleeding since it can be intermittent.
- Colorectal polyps at a very early stage may not bleed. This is the reason why it is safe to periodically control (once a year) people over the age of 45.

CONTENT

- 20 hermetically sealed aluminum bags containing 1 BOWEL FOB TEST cassette and 1 desiccant bag;
- 20 vials with red screw cap, with collection stick and dropper, containing enough BOWEL FOB TEST

BOWEL FOB RAPID TEST

Rapid immunochromatographic test for the detection of fecal occult blood in stool samples
Immunochromatographischer Schnelltest zum Nachweis von fäkalem okkultem Blut in Stuhlproben
Test d'immunochromatographie rapide de détection du sang occulte dans les échantillons de selles
Test immunocromatografico rapido per la rilevazione di sangue occulto fecale in campioni di feces
Test imunocromatográfico rápido para la detección de sangre oculta en heces en muestras de heces
Teste imunocromatográfico rápido para detecção de sangue oculto nas fezes em amostras de fezes
Szybki test immunochromatograficzny do wykrywania krwi utajonej w kale

diluent for 20 tests;

- 1 instructions for use leaflet.
- Material required but not supplied: a device for measuring time (e.g. stopwatch, clock).

SAMPLE COLLECTION AND PREPARATION

Allow the test, sample and buffer to reach room temperature (15-30°C) before testing.

- Unscrew the top of the sample collection device containing the sample collection probe.
Note: Avoid dilution of stool specimens with urine or water from the toilet bowl. Patients should not collect samples during their menstrual period or 3 days before and after it, or if they have bleeding hemorrhoids, blood in urine or if they experience strain during their bowel movement. Alcohol, aspirin and other medications taken in excess may cause gastrointestinal irritation, resulting in occult bleeding. Such substances should be discontinued at least 48 hours prior to testing. No dietary restrictions are necessary before testing. Perform the test as soon as possible after sample collection.
- Collect the stool sample with the top of the collection device by plunging it into 3 different places of the same stool sample –**fig. A** and put it in the collection vial.
Do NOT put the collection stick back into the tube in between collecting the 3 samples.
- Screw on the cap again, immersing the stick in the stool extraction liquid, and shake the vial for at least 10 seconds. –**fig. B**
- Proceed with test procedure.

TEST PROCEDURE

- Open the aluminum bag, remove the test cassette and place it on a clean and level surface. Throw away the desiccant bag.
- Break the end of the red cap, exposing the dropper portion. –**fig. C**
- Dispense 3 drops of the diluted stool sample in the sample well (S) shown on the test cassette. Wait 1-2 seconds between dispensing one drop of diluted stool sample. –**fig. D**
- Wait 5 minutes.

RESULTS INTERPRETATION

READ THE RESULTS AT EXACTLY 5 MINUTES.
The intensity of the line colors is not relevant for the purposes of interpretation of the Test's results.

POSITIVE RESULT

Two coloured lines appear in the reading window by the Test (T) and Control (C) signs. The T line may be less intense (lighter) than the C one. This result means there is blood in the faeces.

NEGATIVE RESULT

A coloured line appears in the Control (C) region. No line appears in the Test (T) region. This result means that there is no blood in the faeces or that its concentration is below the test's detection limit (40 ng/mL).

INVALID RESULT

The Control (C) line does not appear. An insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the Control line (invalid result). Review the procedure and repeat the test with a new device and a new sample.

QUALITY CONTROL

An internal procedural control is included in the test cassette. A coloured line appearing in the Control line region (C) is considered an internal positive procedural Control, confirming sufficient specimen volume and correct procedural technique.

PERFORMANCE CHARACTERISTICS

To assess the performance of BOWEL FOB RAPID TEST, a clinical evaluation was performed by comparing the results obtained against Laboratory Reference Method. 148 samples were tested with both methods (Table 1).

Specificity = 99 % (CI 95%: 97 – 100%)
Sensitivity = 78.7 % (CI 95%: 67 – 90.4%)
Accuracy = 92.6% (CI 95%: 88.4 – 96.8%)

MEASURING RANGE

The study concerning the Limit Of Detection (LOD) was conducted with the Human Haemoglobin Calibrator SENTINEL CH. The LOD of BOWEL FOB RAPID TEST is 40 ng/mL. No adverse effect on T-line formation (prozone effect) was observed for FOB value up to 3 mg/mL.

CROSS-REACTIVITY

The following compounds were tested with BOWEL FOB RAPID TEST and no interference was observed. **Table 2. Cross-reactants tested.**

SUBSTANCE	CONCENTRATION ASSAYED
Bilirubin	4.74 mg/dL (81.0 mmol/L)
Albumin	15 g/dL (2.25 mmol/L)
Glucose	600 mg/dL (32 mmol/L)
Urea	198.2 mg/dL (33 mmol/L)
Uric acid	48.58 mg/dL (2.89 mmol/L)
Bovine – pig – rabbit – horse – chicken – sheep Haemoglobin	up to 2 mg/mL

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