
Test your health

PRIMA LAB SA
Via Antonio Monti 7
CH-6828 Balerna - SWITZERLAND
support@primallabs.ch
primallabs.ch

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REF CE2H88000KIT01E5

REVISIONS

REVISION	DATE	CHANGES
1.2	10/2025	Add statement where the SSP is available

BIBLIOGRAPHY / LITERATUURVERZAMELNIS / RÉFÉRENCES / BIBLIOGRAFIA / BIBLIOGRAFIE

- Lindfors, K., Ciacci, C., Kurppa, K. et al. Coeliac disease. Nat Rev Dis Primers 5, 3 (2019).
- Al-Toma et al. European Society for the Study of Coeliac Disease (EESCD) guideline for coeliac disease and other gluten-related disorders. United European Gastroenterol J. 2019 Jun;7(5):583-613.
- King, et al. 2020. Incidence of Celiac Disease is Increasing Over Time: A Systematic Review and Meta-Analysis. The American Journal of Gastroenterology 115 (4): 507-25.
- Perrinello G., Da Re M., Grizzo F., Carnelliti S., Cozzi M., Marinoni F., Villalta D. Diagnostic accuracy of a novel point-of-care test for simultaneous detection of anti-transglutaminase IgA and anti-deamidated gliadin IgG antibodies. J Clin Lab Anal. 2024

SYMBOLS / SYMBOLE / SYMBOLES / SIMBOLOGIA / SIMBOLEN

IVD	In vitro diagnostic device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / In-vitrodiaagnostisch hulpmiddel
	Read the instructions before use / Beachten Sie die Gebrauchsanweisung / Consulter le mode d'emploi / Leggere attentamente le istruzioni per l'uso / Lees de instructies voor gebruik
STERILE R	Sterilised using irradiation / Sterilisation durch Bestrahlung / Méthode de stérilisation utilisant l'irradiation / Sterilizzato mediante radiazioni / Gesteriliseerd door middel van bestraling
	Expiry date (last day of the month) / Verwendbar bis (letzter Tag des Monats) / Utiliser jusqu'à (dernier jour du mois) / Utilizzare entro (ultimo giorno del mese) / Vervaldatum (laatste dag van de maand)
	Legal manufacturer / Hersteller / Fabricant / Fabricante / Juridische fabrikant
	CE marking / CE-Kennzeichnung / Marquage CE / Marchio CE / CE-markering
	Temperature limits / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Temperatuurgrenzen
	Do not reuse / Nicht wiederverwenden / Ne pas réutiliser / Non riutilizzare / Niet opnieuw
	Sufficient for <-> tests / Ausreichend für <-> Tests / Suffisant pour <-> tests / Sufficiente per <-> test / voldoende voor <-> tests
REF	List number / Katalognummer / Code produit / Codice prodotto / Lijstnummer
LOT	Lot number / Chargenbezeichnung / Code du lot / Lotto prodotto / Lotnummer
	Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandatario nella Comunità Europea / Gennettighe vertegenwoordiger in de Europese Gemeenschap
UDI	Unique Device Identification / Eindeutige Geräteidentifikation / Identification unique du dispositif / Identificativo unico del dispositivo / Unieke apparaatidentificatie
	Importer / Importateur / Importateur / Importatore / Importeur

MEDICAL DEVICES INSIDE THE KIT / MEDIZINISCHE GERÄTE INNERHALB DES KITS / DISPOSITIFS MÉDICAUX À L'INTÉRIEUR DU KIT / DISPOSITIVI MEDICI INTERNI AL KIT / MEDISCHE HULPMIDDELEN IN DE KIT

STERILE R	Lancet / Lanzetta / Bloc autopiqueur / Lancetta / Lancet
	Huian Tianda Medical Instruments Co., Ltd. No.198 East Songjiang Road, Huayin, Economic/Technology Development Zone, 223000 - Huaiyin City, Jiangsu - P. R. China
	Romaxco Sociedad Limitada Calle de Almansa 55, 1D, Madrid 28039, Spain
	Phosma Innovative Healthcare Manufacturers Pvt. Ltd. E-209, SHI Mahape Road, Electronic Zone, MIDC, TTC Industrial Area, Mahape, Navi Mumbai 400 710 MH I India
	Advena Ltd., Tower Business Centre, 2nd Flr, Tower Street, Swatara, BKR 4013, Malta

Alcohol cleansing gauze / Alkohol-Reinigungsgaze / Gaze nettoyante à l'alcool / Garza detergente alcolica / Reinigingsgaze met alcohol

0483

CELIAC DISEASE SCREENING TEST

Rapid immunochromatographic manual self-test for the qualitative detection of anti-Deamidated Gliadin Peptide IgG and anti-tissue Transglutaminase IgA in human capillary blood sample for the screening of celiac disease

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T1	T2
measurement range limit minimum: 5 U/mL maximum: 302 U/mL*	measurement range limit minimum: 12 U/mL maximum: 128 U/mL*
SENSITIVITY** 84.30% (95% CI: 76.77% - 89.71%)	SENSITIVITY** 98.31% (95% CI: 94.03% - 99.53%)
SPECIFICITY** 90.10% (95% CI: 82.73% - 94.53%)	SPECIFICITY** 98.02% (95% CI: 93.07% - 99.46%)
ACCURACY** 86.94% (95% CI: 81.87% - 90.75%)	ACCURACY** 98.17% (95% CI: 95.40% - 99.29%)

* Maximum limit of the ThermoFisher Phadia ImmunoCAP EIA DGP-IgG and tTG-IgA reference system at which the PRIMA Test could be verified to be working perfectly.

** Clinical study conducted on 266 real samples

CROSS-REACTIVITY AND INTERFERING SUBSTANCES TESTED

To determine the effects of possible interferences, the substances listed in Tab. 1 were added to human serum samples, and the concentrations listed below were therefore added to physiological concentrations (if applicable).

Table 1.

SUBSTANCE	CONCENTRATION
Cholesterol	139 mg/dL
Glucose	262 mg/dL
Total bilirubin	4.00 mg/dL
Direct bilirubin	2.21 mg/dL
Triglycerides	254 mg/dL
Uric acid	10.20 mg/dL
Serum albumin	3.98 g/dL
Rheumatoid factor	250 IU/mL
Immunoglobulin A (IgG)	500 mg/dL
Immunoglobulin A (IgA)	50 mg/dL
Immunoglobulin M (IgM)	50 mg/dL
Isoopropanol	25%
Ethanol	25%
Human anti-mouse antibodies (HAMA)	40 ng/mL

No interference was found associated with CELIAC DISEASE SCREENING TEST.

Note: The result is reliable as long as the instructions are carefully followed and the test is stored under the specified conditions. However, false positive or negative results may be obtained, especially if the device comes into contact with other liquids before being used, if the amount of blood and/or diluent is insufficient, if the number of drops added into the well is incorrect or if the reading time of the result is not respected. The black line on the plastic pipette provided allows you to be sure that the volume of blood collected is correct.

CAN I TAKE THE TEST IF I AM ON A GLUTEN-FREE DIET? No, it is not recommended. On gluten-free diet, both the symptoms and the level of antibodies associated with the disease decrease until they are no longer detectable. Therefore, this may result in celiac testing negative.

WHAT TO DO IN CASE OF A POSITIVE RESULT? In the case of a positive result, you should consult your doctor, who will take the appropriate steps for a definitive diagnosis. Do not stop taking gluten until it is confirmed by a specialist.

WHAT TO DO IN CASE OF A NEGATIVE RESULT? If the result is negative, it means that celiac-specific antibodies are not present in the sample or that their level is below the limit detectable by the device. If the symptoms associated with the disease persist, it is still recommended that you see a doctor.

CAN I USE THE CELIAC DISEASE SCREENING TEST TO TEST WHETHER I AM GLUTEN INTOLERANT? No, the test results give no indication of possible gluten intolerance. Gluten intolerance is a pathological condition that differs from celiac disease, although the symptoms may be similar, and its diagnosis involves ruling out celiac disease first.

If you need additional information, assistance or would like to report a problem, please contact us at support@primallabs.ch. In case of serious incidents occurred in relation to the device please inform the local health authority and manufacturers.

A summary of safety and performance (SSP) for this device is available at <https://ec.europa.eu/tools/eudamed>. This is the SSP location after the launch of European Database on Medical Devices. Search for the device using the UDI-DI provided on the outer packaging of the device.

WARNING: Read the entire procedure carefully before performing the test!

- Wash your hands with soap and warm water, rinse with clean water and dry. Note: The use of warm water facilitates capillary blood sampling as it induces vasodilation.
- Prepare the necessary material as follows: open the aluminium sachet, only take out the cassette and throw away the desiccant sachet. Open the plastic sachet containing the pipette and use the provided gauze to clean the puncture site. Allow the area to dry without pulling it. Pull out and discard the released cap. -FIG. B
- Position the hand downward and gently massage from the palm to the chosen phalanx repeatedly: middle or ring fingers are recommended - avoid the fifth finger. The puncture should be made slightly off-center, in the fleshy area of the thickest part of the phalanx. Press the lancet perpendicular to the fingertip, on the side from which the cap was removed, while keeping it firmly between the fingers. -FIG. C
- The tip of the lancet automatically retracts safely after use. If the lancet does not work properly, use the second one provided. If the latter is not necessary, it can be disposed of without special precautions. Note: make a single puncture. Do not perform an immediate double sticking at the same site.
- Wipe away the first drop of blood with the cleaning gauze and allow another drop to form, by keeping your hand downward and continue to massage the finger until the specimen volume is adequate. -FIG. D
- Pick up the pipette without pressing the bulb. Two withdrawal methods are suggested: -FIG. E: hold the pipette at a slight incline without squeezing the bulb and place it in contact with the blood droplet, it will enter on its own by capillary action. Remove the pipette when it reaches the black line. If there is not enough blood, continue massaging the finger until the line is reached. -FIG. F: place the pipette on a flat, clean surface, allowing the tip to protrude, at which point the drop of blood should come into contact with the pipette, it will enter on its own by capillary action. Move your finger away when it reaches the black line. If there is not enough blood, continue massaging the finger until the line is reached. Avoid, as far as possible, breaking contact between the blood and the pipette so as to prevent the formation of air bubbles.
- Deposit the blood collected with the pipette into the sample well (S), by squeezing the bulb of the pipette. -FIG. F
- UnscREW the blue cap on the dropper vial (leave the white cap tightly screwed on). Deposit 1 drop of diluent into the well indicated on the cassette (S). -FIG. G
- Wait 10 minutes and then interpret the results as follows. Do not move the device while waiting.

LIMITATIONS

- Do not use under the age of 16 years.
- Not intended for use in urgency/emergency situations.
- People suffering from haemophilia or taking anticoagulants could have an excessive bleeding after the finger-prick.
- Situations that can lead to a "false positive" result: Crohn's disease, ulcerative colitis, systemic lupus erythematosus, chronic liver disease, in subjects with selective IgA deficiency with total IgA titres below 40 mg/dL, the test could give a false positive result due to non specific interactions among device components.
- Situations that can lead to a "false negative" result in people with celiac disease: starting a gluten-free diet before testing can lead to negative results. Furthermore, it is reported that up to 2% of celiac patients have no specific circulating antibodies, defining a seronegative celiac condition.

PRECAUTIONS

- In vitro diagnostic device for individual use.
- Read these instructions carefully before testing. The test is

INTERPRETATION OF THE RESULTS

READ THE RESULT AT 10 MINUTES.
DO NOT READ THE RESULTS AFTER 15 MINUTES.

POSITIVE*

In addition to the coloured line appearing at the C (Control) mark, the result is considered positive if:

- one coloured line appears near the T1 mark (positive for anti-DGP IgG), or

- one coloured line appears near the T2 mark (positive for IgA anti-tTG), or

- two coloured lines appear near the T1 and T2 marks (positive for both anti-DGP IgG and IgA anti-tTG)

Consult a doctor, as this result indicates probable celiac disease. *NOTE: The colour intensity in the control and test line regions may vary (as can be seen in the pictures). Therefore, any colour intensity in the region of the T1 and T2 lines must be considered positive. The user is advised to pay close attention in detecting even weak intensity test bands, which must in any case be interpreted as a positive result.

NEGATIVE

Only one coloured line appears at the C (Control) mark. This result indicates that no anti-DGP IgG or IgA anti-tTG antibodies are present in the blood or that their level is below the limit detectable by the test.

INVALID

No control line appears at the C (Control) mark. The most likely reasons for an invalid result are: insufficient sample volume or the procedure has not been performed correctly.

Review the procedure and repeat the test with a new device and a new sample.

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