## Test your health

**PRIMA LAB SA**  
Viale Serfontana, 10  
CH-6834 Morbio Inferiore - SWITZERLAND  
support@primalabs.ch  
primalabs.ch

**EC REP** QbD RepS BV  
Groenenborgerlaan 16  
2610 Wilrijk - Belgium

100013\_FU\_11\_9.0 03/2026

**REF** 100013-x

Registered trademarks.

**PRIMA**  
HOME TEST

INSTRUCTIONS FOR USE  
GEBRAUCHSANWEISUNG  
NOTICE D'UTILISATION  
ISTRUZIONI PER L'USO

## VAGINAL PH TEST

Rapid self-test for the semiquantitative determination of pH in female vaginal swab specimens.

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**CE**  
0123  
Copen Italia SpA  
Via Perotti 10  
25125 Brescia  
Italy

### BIBLIOGRAPHY / BIBLIOGRAPHIE / BIBLIOGRAFIE / BIBLIOGRAFIA

- Pavletic, A. J., Hawes, S. E., Geske, J. A., Bringe, K., & Palack, S. H. (2004). Experience with routine vaginal pH testing in a family practice setting. *Infectious Diseases in Obstetrics and Gynecology*, 12(2), 63–68.
- Huppert, J. S., Hesse, E. A., Bernard, M. C., Bates, J. R., Gaydos, C. A., & Kahn, J. A. (2012). ACCURACY AND TRUST OF SELF-TESTING FOR BACTERIAL VAGINOSIS. *The Journal of Adolescent Health : Official Publication of the Society for Adolescent Medicine*, 51(4), 400–405.

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
### INTENDED USE

The VAGINAL Ph TEST is a rapid dry chemical method for the semiquantitative detection of pH in female vaginal swab specimens to aid in the diagnosis of bacterial vaginosis.

### SUMMARY

Vaginal infections are quite common and often a recurring problem among women of all age groups. An indication of abnormal acidity of the vaginal discharge can help evaluate whether the vaginal symptoms are likely caused by an infection that may require follow-up with a healthcare professional. An acidic vaginal pH value of 3.8 to 4.5 is a basic requirement for the optimal functioning of the body's own system of protecting the vagina. This system can effectively avoid colonization by pathogenic germs and the occurrence of vaginal infections. Under the influence of lacto bacteria, the glycogen released from the upper vaginal epithelium degrades to form the sugar glucose and then lactic acid. This results in an acid environment forming in the vagina, with a pH value that is usually 3.8 to 4.5. This acidic vaginal environment forms effective protection against infections, because the vast majority of bacteria that causes disease cannot thrive under such conditions. Some vaginal infections known to be associated with abnormal vaginal acidity have no non-prescription treatments, and must be treated by a healthcare professional. Some cases of abnormal discharge with abnormal acidity may be caused by a combination of infections. Vaginal infections may be very risky for pregnant women and dangerous for the mother and her fetus, thus symptoms and results obtained with this test should always be discussed with a healthcare professional.

### PRINCIPLE OF THE TEST


A pH paper is used to determine pH levels in specimen obtained through a vaginal swab. As soon as the swab containing the specimen from vagina comes in contact with the pH paper encased in the test panel, a color change occurs. Different pH levels can result in different grades and shades of color with each color grade/ shade specifying a particular pH level on the color scale. The color obtained on the pH paper after application of the vaginal swab is then compared with the color scale onto the test aluminum pouch. The result is considered to be equivalent to the numerical value of the color on the color scale closest in appearance to the color obtained with the test.

NOTE: The initial colour of the pH paper is pale green.

### KIT CONTENT

\* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the external box.

- \* Aluminium pouches with printed color scale containing 1 VAGINAL PH TEST

- \* sterile vaginal swabs

- 1 Instructions for use

Materials Required But Not Provided: toilet paper

### PRECAUTIONS

- For self-test in vitro diagnostic use only.
- Do not eat, drink or smoke in the area where the specimens and kits are handled.
- After use, all components can be disposed according to local regulations; all specimen should be considered potentially hazardous and handled in the same manner as an infectious agent.
- Humidity and temperature can adversely affect results.
- Do not use test if pouch is damaged.
- Use each test panel only once.
- Only use the swab provided in the kit.
- Only open the package of swab before specimen collection.
- The test is stable through the expiration date printed on the sealed pouch.
- The test must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.
- Remove the test panel from the seal pouch and use it within one hour. Best result will be obtained if the test is performed immediately after opening the foil pouch.
- Place the test panel on a clean and balance desk or take it by your free hand.
- Store as packaged in the sealed pouch at room temperature or refrigerated (2-30°C).
- External use only. DO NOT SWALLOW.
- In case of difficulties in color identification (such as Daltonism), ask for help in test reading.
- DO NOT use this test until 72 hours after the application of vaginal preparations like contraceptive creams and vaginal medical products (antifungal agents, vaginal suppositories, creams, gels, douching etc.).
- DO NOT use this test until 48 hours after sexual intercourse.
- DO NOT use this test until 5 days after the period is over.

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
**ENGLISH SPECIMEN COLLECTION AND PREPARATION**

The Vaginal pH Test can be performed using female vaginal swab specimen. Allow the test and the swab specimen to reach room temperature (15-30°) prior to testing. Wash hands before starting the test. Open the package, remove the swab from pouch, do not touch the head of the swab (polyester head). Separate the labia so the vagina is exposed. Insert the head of swab gently into the vagina and gently press the swab to a point, when approximately 2 inches of swab (approximately length of an AA size battery is into the vagina for approx. 10 seconds). A time of approximately 10 seconds is important to ensure that the head of swab becomes well moistened with vaginal secretions. Sufficient moisture is very important for the test evaluation. Withdraw the swab from the vagina for testing. Apply the head of swab onto pH area of the test panel for at least 5 times, make sure enough secretion is applied on the pH area of the test panel. Read the result immediately while the measurement zone is still wet. Compare the color in the pH measurement zone with the color scale on the foil pouch and note the value on the color scale that is closest to the color on the test unit. After evaluation, dispose of the used panel and swab.

**RESULTS INTERPRETATION**

Read the result immediately while the measurement zone is still wet. Results are obtained by visually comparing the PH measurement zone with the color scale on the foil pouch.

**< 1 min.**

**NORMAL RANGE**

Vaginal pH levels in the range 3.8-4.4 are within the normal range. The color of indicator is Yellow-Green.

**ABNORMAL RANGE:**

Vaginal pH levels over 4.4 or under 3.8 are outside of the normal/healthy range. When pH level is over 4.4, color of the indicator is Green-Blue, when pH levels is under 3.8, it is bright Yellow.

**LIMITATIONS**

- Some menopausal women may have an elevated vaginal pH, and therefore may receive an abnormal result with Vaginal pH Rapid Test Panel even though they are not suffering from a bacterial or parasitic infection.
- Certain circumstances (menstruation period, sex intercourse, use of vaginal medical products) can temporarily alter the pH level of vagina and may lead to false results. Therefore, please note the precautions when performing the test:
- Do not wash or shower the area immediately before taking the measurement.
- Be aware that urine may cause a false test result.
- Never begin any treatment before you have discussed the result of the test with a doctor.
- DO NOT use this test until 72 hours after the application of vaginal preparations like contraceptive creams and vaginal medical products (antifungal agents, vaginal suppositories, creams, gels, douching etc.).
- DO NOT use this test until 48 hours after sexual intercourse.
- DO NOT use this test until 5 days after the period is over.
- In cases of improper use, rupture of hymen may occur as it occurs when using tampons.
- This is NOT a test for sexually transmitted diseases, such as HIV, chlamydia, herpes, gonorrhea or syphilis.

**F.A.Q. – QUESTIONS AND ANSWERS****COULD I USE THE TEST AT ANYTIME?**

No. There are some restrictions on when to use this test because at certain circumstances can temporarily alter the pH level of the vagina and may lead to false results. Therefore, please find the "Limitations" section in this package insert and use the test.

**CAN THE POSITIVE RESULTS SHOW THAT THE SUBJECT HAS BACTERIAL VAGINOSIS?**

Not necessarily. Some menopausal women may have an elevated vaginal pH, hence may receive a positive result even without suffering from bacterial or parasitic infection. Initiate any treatment, only after consultation with doctor.

**WHAT IS THE REASON THAT MAY RESULT IN WRONG TEST RESULTS?**

Some factors can affect the test results, including inaccurate specimen collection, insufficient moistening of the swab, insufficient application of specimen on the test area.

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