Test your health

PRIMA LAB SA
Viale Serfontana, 10
CH-6834 Morbio Inferiore - SWITZERLAND

support@primlabsa.ch
primlabsa.ch

EC REP QbD Rep5 BV
Groenenborgerlaan 16
2610 Wilrijk - Belgium

100058_IFU_11_9.0 03/2026

REF 100058-x

Registered trademarks.

INSTRUCTIONS FOR USE
GEBRAUCHSANWEISUNG
NOTICE D'UTILISATION
ISTRUZIONI PER L'USO

0483

URINARY TRACT INFECTIONS TEST

Self-Test for the detection of leukocytes, blood, nitrites and proteins in urine

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PRIMA
HOME TEST

ENGLISH THE URINARY TRACT

A urinary infection represents the most common disease of the urinary tract which includes the urethra, the bladder, the ureter and the kidneys. Men, women and children are likely to experience a urinary infection. It's mostly women who suffer from urinary infections, since the short urethra favours the penetration of germs. However, elderly males are also affected if they have an enlarged prostate which obstructs the urine flow. In healthy people, urine is sterile (i.e. it doesn't contain any micro-organisms). One of the best ways to keep your urinary tract sterile is to empty your bladder completely at regular intervals. Generally, an infection starts in the urethra and may then spread into the upper urinary tract as far as the kidneys. The symptoms vary considerably: burning when emptying the bladder, or a strong urge to urinate. The urine may also be cloudy or have a strong odour.

The URINARY TRACT INFECTIONS TEST contains firm plastic strips onto which several separate reagent areas are affixed. The test is for the qualitative and semi-quantitative detection of the following analytes in urine: Leukocytes, Blood, Nitrite and Protein by comparing the color blocks on the color chart for results.

TEST PRINCIPLE

Leucocytes: This test reveals the presence of granulocyte esterases. The esterases cleave a derivatized pyrazole amino acid ester to liberate derivatized hydroxy pyrazole. This pyrazole then reacts with a diazonium salt to produce a beige-pink to purple color.

Blood: This test is based on the peroxidase-like activity of hemoglobin which catalyzes the reaction of diisopropylbenzene dihydroperoxide and 3,3',5,5'-tetramethylbenzidine. The resulting color ranges from orange to green to dark blue.

Nitrite: This test depends upon the conversion of nitrate to nitrite by the action of Gram negative bacteria in the urine. In an acidic medium, nitrite in the urine reacts with p-arsanilic acid to form a diazonium compound. The diazonium compound in turn couples with 1 N-(1-naphthyl)- ethylenediamine to produce a pink color.

Protein: This reaction is based on the phenomenon known as the "protein error" of pH indicators where an indicator that is highly buffered will change color in the presence of proteins (anions) as the indicator releases hydrogen ions to the protein. Colors range from yellow to yellow-green for negative results and green to green-blue for positive results.

F.A.Q. - QUESTIONS AND ANSWERS

WHAT SHOULD I DO IF MY TEST RESULT IS POSITIVE? Remember that a positive result doesn't mean that all four substances have been detected in your urine. Even if your result is positive for just one of them, it is most likely that something is wrong in your urine, even if the reason may not be a urinary infection. Get in touch promptly with your own doctor, who will be able to give a more accurate diagnosis. When you visit your doctor, please take these instructions with you so that he/she will be better informed as to the type of test you have performed.

WHAT SHOULD I DO IF MY TEST RESULT IS NEGATIVE? Remember that your test result is only negative if the result on the test field for all four substances is negative. But if you still feel the signs of a UTI or have any other symptom, then contact your own doctor to arrange a more thorough examination.

WHAT IS THE ACCURACY OF URINARY TRACT INFECTIONS TEST? One study carried out in parallel with another test for the urine showed an accuracy of 97.6% for leukocytes, 96% for the blood, 100% for nitrites and 88% for proteins.

CONTENT

THE BOX CONTAINS ALL NECESSARY TO PERFORM 3 TESTS:

- 3 sealed aluminium pouches containing:
 - 1 URINARY TRACT INFECTIONS TEST strip
 - 1 desiccant bag

- 1 color chart

- 1 instructions for use leaflet

* Additionally required: a device to measure time (i.e. timer, watch), dry, residue-free container.

* The desiccant should be disposed of with your normal household waste.

PRECAUTIONS

1. Read the instructions before performing the test. In order to obtain a reliable result, the instructions should be followed carefully.
2. Keep out of the reach of children.
3. Do not use after the expiry date or if the pouch is damaged or if the strip is discolored.
4. Follow the indicated time strictly.
5. Store between +2 °C and +30 °C, DO NOT FREEZE.
6. Use the test only once. Do not dismantle and touch the reagent areas of the test strip.
7. For external use only.
8. In-vitro diagnostic device for self-testing.
9. After use, all components can be disposed according to local regulations; all specimen should be considered potentially hazardous and handled in the same manner as an infectious agent.
10. In case of difficulties in color identification (such as Daltonism), ask for help in test reading
11. The test strip should remain in the sealed pouch until use.

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



ENGLISH

PROCEDURE

ATTENTION: It's recommended to take a sample of urine for the test in the early morning since it's the most concentrated. The urine used for the test should not come into contact with water from the toilet or any disinfectant or cleaning substances.

For women only: the test should not be performed during or for three days after your menstrual period. The urine sample should not be contaminated with vaginal fluids since this may produce a misleading result.

Do not make any important medical decision without first referring to your doctor.

COLLECT URINE: Collect part of the urine using a clean cup without any residual detergents. Make sure to fill up the cup with urine. - **FIG. A**

Attention: The use of urine preservatives is not recommended. If testing cannot be done within an hour after voiding, refrigerate the specimen immediately and let it return to room temperature before testing. Prolonged storage of unpreserved urine at room temperature may result in microbial proliferation with resultant changes in pH. A shift to alkaline pH may cause false positive results with the protein test area. Urine containing glucose may decrease in pH as organisms metabolize the glucose. Contamination of the urine specimen with skin cleansers containing chlorhexidine may affect protein test results.

PERFORMING THE TEST:

1) Open the aluminium pouch where indicated and take out the test strip. Do not touch the test fields and dispose of the attached desiccant with your normal household waste. Once opened the pouch, it is recommended to perform the test immediately - **FIG. B**

2) Dip the test strip in the urine sample. **ATTENTION:** Make sure that all four test fields are immersed for about **1-2 seconds**. - **FIG. C**

3) Then remove the test strip and wipe off any surplus urine against the rim of the container or with an absorbent material (e.g. a paper towel) to avoid mixing chemicals from adjacent reagent areas.

Place the test strip horizontally on the opened aluminium pouch so as to avoid spilling any urine on any other surfaces. - **FIG. D Wait for 2 min** (do not read results after 3 minutes).

RESULT INTERPRETATION

READ THE RESULT AFTER 2 MINUTES.

READ THE RESULT SEPARATELY FOR EACH PARAMETER.

There are more possible tones and color intensities. Color changes on the edges of the test pads or color changes after more than 3 minutes have to be ignored.


NEGATIVE RESULT

The Test field for **LEUKOCYTES** stayed whitish.
The Test field for **BLOOD** stayed mustard yellow
The Test field for **NITRITE** stayed white.
The Test field for **PROTEIN** stayed yellowish.


POSITIVE RESULT FOR LEUKOCYTES

If the color of the test field has changed to **purple**, then leukocytes have been found in your urine. (see **FIG. E** for the negative result).


POSITIVE RESULT FOR BLOOD

If the color of the test field has changed to **green** (or some **green spots** appear on the background), then blood has been found in your urine. (see **FIG. E** for the negative result).

POSITIVE RESULT FOR NITRITE

If the color of the test field has changed to **pink**, then nitrites have been found in your urine. (see **FIG. E** for the negative result).

POSITIVE RESULT FOR PROTEIN

If the color of the test field has changed to **green**, then proteins have been found in your urine. (see **FIG. E** for the negative result).

PERFORMANCE

Analytical sensitivity: **Leucocytes:** Detects leukocytes as low as 9 white blood cells Leu/ μ L in clinical urine. **Blood:** Detects free hemoglobin as low as 0.018 mg/dL or 5 Ery/ μ L in urine specimens with ascorbic acid content of < 50 mg/dL. **Nitrite:** Detects sodium nitrite as low as 0.05 mg/dL in urine with a low specific gravity and less than 30 mg/dL ascorbic acid. **Protein:** Detects albumin as low as 7.5 mg/dL (0.075 g/L).

TECHNICAL NOTES ON PARAMETERS

The Test detects **LEUKOCYTES, BLOOD, NITRITE** and/or **PROTEIN** in urine.

LEUKOCYTES: When taking cephalixin and cephalothin, or high concentration of oxalic acid may also cause test results to be artificially low. Tetracycline may cause decreased reactivity, and high levels of the drug may cause a false negative reaction. High urinary protein may diminish the intensity of the reaction color. The presence of leukocytes in urine is an important symptom of an inflammation of the kidneys and the urinary tract. **BLOOD:** A uniform green color conversion indicates the presence of hemoglobin or hemolyzite erythrocytes; scattered or compacted green spots indicate intact erythrocytes. **NITRITE:** Gram-negative bacteria in urine convert nitrate from food into nitrite. Nitrite reacts with a chemical in the test field and leaves a pink shade. The test result may be distorted if urine does not stay for long in the bladder, due to hunger, a vegetable-free diet or antibiotic treatment and medicines containing phenazopyridine. Comparing the test on a white background may aid in the detection of low nitrite levels, which might otherwise be missed. **PROTEIN:** An indicator on the test field reacts with protein in the urine, changing its colour to green. They may be found where there is inflammation of the bladder or prostate or bleeding in the urinary tract. Infusions containing polyvinylpyrrolidone or medicines containing phenazopyridine may yield a false positive result

Chemical components in the test fields must be viewed as potentially dangerous substances, although they present no hazard provided that all test components are used in accordance with these instructions.

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