Test your health

PRIMA LAB SA
Viale Serfontana, 10
CH-6834 Morbio Inferiore - SWITZERLAND
support@primlabs.ch
primlabs.ch

EC REP QbD Rep5 BV
Groenenborgerlaan 16
2610 Wilrijk - Belgium

100066_IFU_11_9.0_03/2026 REF 100066-x

Registered trademarks

INSTRUCTIONS FOR USE
GEBRAUCHSANWEISUNG
NOTICE D'UTILISATION
ISTRUZIONI PER L'USO

VITAMIN D TEST

Rapid self-test for the semi-quantitative detection of vitamin D in human whole blood

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Huain Tianda Medical Instruments Co., Ltd., No.106 East Songjiang Road, Huaiyin, Economic/Technological Development Zone, 223300 - Huain City, Jiangsu - P. R. China

Riomaviv Sociedad Limitada
Calle de Almansa 55, 10
Madrid 28039, Spain

Alcohol cleansing gauze
Alkohol-Reinigungsgaze
Gaze nettoyante à l'alcool
Garza detergente alcolica

Phoenix Innovative Healthcare Manufacturers Pvt. Ltd., El-209, Shri Mahape Road, Electronic Zone, MIDC, TTC Industrial Area, Mahape, Navi Mumbai 400 710 MH 1 India

Advena Ltd., Tower Business Centre, 2nd Flr, Tower Street, Swatara, BKR 4013, Malta

SYMBOLS / SIMBOLE / SYMBOLES / SIMBOLOGIA

IVD In vitro diagnostic device In-vitro-Diagnostikum Dispositif medical de diagnostic in vitro Dispositivo medico-diagnostico in vitro	Temperature limits Temperaturbegrenzung Limites de temperature Limiti di temperatura	REF List number Katalognummer Code produit Codice prodotto
Read the instructions before use Beachten Sie die Gebrauchsanweisung Consulte le mode d'emploi Leggere attentamente le istruzioni per l'uso	Do not reuse Nicht wiederverwenden Ne pas réutiliser Non riutilizzare	LOT Lot number Chargenbezeichnung Code du lot Lotto Prodotto
STERILE R Sterilized using irradiation Sterilisation durch Bestrahlung Méthode de stérilisation utilisant l'irradiation Sterilizzato mediante radiazione	Sufficient for <v> tests Ausreichend für <v> Tests Suffisant pour <v> tests Sufficiente per <v> test	CE marking CE-Kennzeichnung Marquage CE Marchio CE
Expiry date (last day of the month) Verwendbar bis (letzter Tag des Monats) Utiliser jusqu'à (dernier jour du mois) Utilizzare entro (ultimo giorno del mese)	Importer Importeur Importatore Importatore	Legal manufacturer Hersteller Fabricant Fabricante
Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandatario nella Comunità Europea	Content / Inhalt Contenu / Contenuto	

ENGLISH VITAMIN D

Vitamin D refers to a group of fat-soluble secosteroids responsible for increasing intestinal absorption of calcium, iron, magnesium, phosphate and zinc.

The blood concentration of 25-hydroxy Vitamin D (including D2 and D3) is considered the best indicator of Vitamin D status. Vitamin D deficiency is now recognized as a global epidemic. The health risks associated with Vitamin D deficiency are far more severe than previously thought. Vitamin deficiency has been linked to various serious diseases: osteoporosis, osteomalacia, multiple sclerosis, cardiovascular diseases, pregnancy complications, diabetes, depression, strokes, autoimmune diseases, flu, different cancers, infectious diseases, alzheimer, obesity and higher mortality etc. Therefore, now detecting Vitamin D level is considered as "Medically Necessary Screening Test", and maintaining sufficiently high levels of Vitamin D is necessary to improve bone health, health and well-being of the body.

PRINCIPLE OF THE TEST

The VITAMIN D TEST is an immunoassay based on the principle of competitive binding. During testing, the mixture migrates upward on the membrane chromatographically by capillary action. The membrane is pre-coated with 25 (OH) D antigens on the test line region of the strip.

LIMITATIONS


- 1) The VITAMIN D TEST provides only a semi-quantitative analytical result. A secondary analytical method must be used to confirm the result.
- 2) It is possible that technical or procedural errors, as well as other interfering substances in the whole blood specimen may cause erroneous results.
- 3) As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
- 4) Other clinically available tests are required if questionable results are obtained.

BIBLIOGRAPHY

1. Holick MF (March 2006). „High prevalence of vitamin D inadequacy and implications for health.“ *Mayo Clinic Proceedings*. 81 (3): 353–73.
2. Eriksen EF, Glerup H (2002). Vitamin D deficiency and aging: implications for general health and osteoporosis. *Biogerontology*. 3 (1-2): 73-7.
3. Grant WB, Holick MF (June 2005). Benefits and requirements of vitamin D for optimal health: a review. *Alternative Medicine Review*. 10 (2): 94-111.
4. Moyad MA. Vitamin D: a rapid review. *DermatolNurs*. 2009, 21:25-30.

CONTENT

* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the external box.

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
PRECAUTIONS

1. For self-testing in vitro diagnostic use only.
2. Do not eat, drink or smoke in the area where the specimens or kits are handled.
3. Store in a dry place at 2-30°C (36-86°F), avoiding areas of excess moisture. If the foil packaging is damaged or has been opened, please do not use. The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. DO NOT FREEZE. Do not use after the expiration date.
4. This test kit is intended to be used as a preliminary test only and repeatedly abnormal results should be discussed with doctor or medical professional.
5. Follow the indicated time strictly, using only the specified quantities of blood and diluent.
6. Use the test only once. Do not dismantle and touch the test window of the test cassette.
7. The kit must not be frozen or used after the expiration date printed on the package.
8. Keep out of the reach of children.
9. The test is for external use only. DO NOT SWALLOW.
10. Not recommended for people who take anti-coagulant medications (blood thinners) or people suffering from haemophilia.
11. After using, please dispose of all components according to your local waste disposal laws. Ask your pharmacist for advice.
12. If you test another person, handle the blood samples as if containing a blood-borne infectious agent. Wear disposable gloves during use.

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







DIRECTIONS FOR USE

- 1) Wash hands with soap and warm water, rinse with clean water and allow to dry. If this is not possible, use the provided gauze as an alternative.
2) Prepare the necessary material as follows: open the aluminium pouch, take out only the test cassette and throw away the desiccant bag.
3) Carefully rotate the protective cap of the sterile lancet 360° without pulling it.
4) Carefully massage the finger chosen for the puncture...
5) Holding the hand down, massage the finger until a large drop of blood forms.
6) Take the pipette without pressing the bulb.
7) Deposit the blood collected into the Specimen well (S) of the cassette, by squeezing the dropper bulb.
8) Wait for the blood to be totally dispensed in the well.
9) Wait for the colored line(s) to appear. Read results at 10 minutes.

RESULTS INTERPRETATION

Please refer to the illustration and compare the T line intensity with "Vitamin D Color card" provided with the kit. As control system, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

DEFICIENT 0-10 (ng/mL) 0-25 (nmol/L) Two distinct colored lines appear. One is in the control region (C) and another should be in the test region (T). The line intensity in the test region (T) is equal to or darker than 10 ng/mL line depicted on Vitamin D color card provided with the kit.

INSUFFICIENT 10-30 (ng/mL) 25-75 (nmol/L) Two colored lines appear. One is in the control region (C) and another should be in the test region (T). The line intensity in the test region (T) is darker than the 30 ng/mL line depicted on the Vitamin D color card provided with the kit and lighter than 10 ng/mL line depicted on Vitamin D color card provided with the kit.

SUFFICIENT 30-100 (ng/mL) 75-250 (nmol/L) Two colored lines appear. One line should be always in the control region (C) and faint colored line appears in the test region (T). The line intensity in region (T) is equal or lighter than 30 ng/mL depicted on Vitamin D Color Card.

EXCESS >100 (ng/mL) >250 (nmol/L) One colored line appears in the control line region (C). No apparent colored line appears in the test line region (T). If the result is excess, it is recommended to consult a physician.

INVALID

Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

F.A.Q. - QUESTIONS AND ANSWERS

HOW DOES THE VITAMIN D TEST WORK? In medicine, a 25-hydroxy vitamin D is the main storage form of vitamin D in the body. Therefore, the overall status of vitamin D can be determined by detecting the content of 25-hydroxy vitamin D. 25-hydroxy vitamin D level less than 30 ng/mL in case of a positive result, indicates vitamin D deficiency or insufficiency.
WHEN SHOULD THE TEST BE USED? The clinical application of vitamin D is mainly for diagnosis, treatment and monitoring of rickets (children), osteomalacia, postmenopausal osteoporosis and renal osteopathy.
CAN THE RESULT BE INCORRECT? The results are accurate as far as the instructions are carefully respected. Nevertheless, the result can be incorrect if the VITAMIN D TEST cassette gets wet before test performing or if the quantity of blood dispensed in the sample well is not sufficient, or if the number of buffer drops are less than 2 or more than 3.

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