PRIMA®

PRIMA Lab SA
Viale Serfontana, 10
CH-6834 Morbio Inferiore - SWITZERLAND
support@primalabsa.ch
primalabsa.ch

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
POSITIVE RESULT

Two coloured lines appear in the reading window by the Test (T) and Control (C) signs. The T line may be less intense (lighter) than the C one. This result means there is Helicobacter pylori Antigen in the faeces.

NEGATIVE RESULT

A coloured line appears in the Control (C) region. No line appears in the Test (T) region. This result means that there is no Helicobacter pylori Antigen in the faeces or that its concentration is below the test's detection limit (10 ng/mL).

INVALID RESULT

The Control (C) line does not appear. An insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the Control line (invalid result). Review the procedure and repeat the test with a new device and a new sample.

QUALITY CONTROL

An internal procedural control is included in the test cassette. A coloured line appearing in the Control line region (C) is considered an internal positive procedural control, confirming sufficient specimen volume and correct procedural technique.

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TEST PROCEDURE

- 5) Open the aluminium bag, remove the test cassette and place it on a clean and level surface. Throw away the desiccant bag.
6) Break the end of the red cap, exposing the dropper portion. -fig. C
7) Dispense 3 drops of the diluted stool sample in the sample well (S) shown on the test cassette. Wait 1-2 seconds between dispensing one drop of diluted stool sample. -fig. D
8) Wait 10 minutes.

RESULTS INTERPRETATION

READ THE RESULTS AT EXACTLY 10 MINUTES. The intensity of the line colors is not relevant for the purposes of interpretation of the Test's results.

CONTENT

- 20 hermetically sealed aluminium bag containing 1 HELICOBACTER PYLORI ANTIGEN TEST cassette and 1 desiccant bag;
• 20 vials with red screw cap, with collection stick and dropper, containing enough HELICOBACTER PYLORI ANTIGEN TEST diluent for 20 tests;
• 1 instructions for use leaflet.
Material required but not supplied: a device for measuring time (e.g. stopwatch, clock).

PERFORMANCE CHARACTERISTICS

To assess the performance of HELICOBACTER PYLORI ANTIGEN RAPID TEST, a clinical evaluation was performed by comparing the results obtained against Helicobacter pylori Test ELISA. 100 samples were tested with both methods (Table 1). Specificity = 96.6 % (CI 95%: 91.8 - 100%) Sensitivity = 95.2% (CI 95%: 88.8 - 100%) Accuracy = 96% (CI 95%: 92.1 - 99.8%) MEASURING RANGE The study concerning the Limit Of Detection (LOD) was conducted with the Helicobacter Antigen Biospecific Inc. The LOD of HELICOBACTER PYLORI ANTIGEN RAPID TEST is 10 ng/mL. No adverse effect on T-line formation (prozone effect) was observed for antigen value up to 340.000 ng/mL. CROSS-REACTIVITY The following compounds were tested with HELICOBACTER PYLORI ANTIGEN RAPID TEST and no interference was observed:

Table 2. Cross-reactants tested.

Table with 2 columns: SUBSTANCE and CONCENTRATION ASSAYED. Lists various bacteria and their concentrations tested against the kit.

SAMPLE COLLECTION AND PREPARATION

Allow the test, sample and buffer to reach room temperature (15-30°C) before testing.

- 1) Unscrew the top of the sample collection device containing the sample collection probe. Note: Avoid dilution of stool specimens with urine or water from the toilet bowl. Perform the test as soon as possible after sample collection.
2) Collect the stool sample with the top of the collection device by plunging it into 3 different places of the same stool sample -fig. A and put it in the collection vial. Do NOT put the collection stick back into the tube in between collecting the 3 samples.
3) Screw on the cap again, immersing the stick in the stool extraction liquid, and shake the vial for at least 10 seconds. -fig. B
4) Proceed with test procedure.

TEST PROCEDURE

- 5) Open the aluminium bag, remove the test cassette and place it on a clean and level surface. Throw away the desiccant bag.
6) Break the end of the red cap, exposing the dropper portion. -fig. C
7) Dispense 3 drops of the diluted stool sample in the sample well (S) shown on the test cassette. Wait 1-2 seconds between dispensing one drop of diluted stool sample. -fig. D
8) Wait 10 minutes.

RESULTS INTERPRETATION

READ THE RESULTS AT EXACTLY 10 MINUTES. The intensity of the line colors is not relevant for the purposes of interpretation of the Test's results.

EN - INSTRUCTIONS FOR USE

HELICOBACTER PYLORI

Helicobacter pylori is a bacterium whose ideal habitat is the human gastric mucosa. The infection is often asymptomatic, but it can sometimes cause gastritis and ulcers in the stomach or in the duodenum segment of the small intestine. The most common symptoms of these conditions are heartburn or pain in the upper part of the abdomen, especially in fasting or in the early morning, when the stomach is empty. Other associated symptoms include nausea, satiety, reflux, loss of appetite, burping, weight loss and diarrhoea. Ulcers can sometimes bleed, causing anaemia over long periods of time. In the long term, Helicobacter pylori infection is associated with an increased risk of gastric cancer. The search for Helicobacter pylori-specific antigens in the feces is a valid and non-invasive method to rapidly identify an infection.

PRINCIPLE OF THE TEST

HELICOBACTER PYLORI ANTIGEN RAPID TEST is a lateral flow test device based upon a principle of solid-phase immunochemotrophie (also called lateral flow technology). The detection of Helicobacter pylori antigens takes place on a test strip embedded in a plastic cassette. The methodology implied is a direct sandwich EIA test, in which the target analyte (Helicobacter pylori Antigen) is recognized by a couple of specific antibodies, one which is coated to colloidal gold particles and the other one is fixed onto the membrane, in correspondence of the Test line. In case of presence of Helicobacter pylori antigens in a patient's feces, the immunocomplexes are then captured by the Test line, thus showing a red-purple band appearance, indicating a positive result. Consequently, two bands are visible (T line - C line), after 10 minutes from the diluent drops application, showing a positive result. In case of absence of Helicobacter pylori antigens or with level lower than the cut-off value only one band is visible in the C line region, showing a negative result.

PRECAUTIONS

- 1. Read these instructions for use carefully before performing the test. The test is reliable if the instructions are followed correctly (reaction times and how the sample is collected).
2. Do not open the sealed aluminium pouch until just before performing the test. The desiccant bag must not be used.
3. Humidity and temperature can adversely affect the results. DO NOT FREEZE.
4. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards during all procedures and follow standard procedures for proper disposal of specimens.
5. The test is for external and single use only. DO NOT SWALLOW.
6. Do not use the kit after the expiry date or if the packaging and/or internal components are damaged.
7. Store the buffer correctly closed with its cap.
8. Avoid cross-contamination of samples by using a new specimen collection tube for each sample.

LIMITATIONS

- 1. As with any diagnostic method, the physician must evaluate the test result in light of other clinical evidence.
2. Despite the reliability of the test, false positive or negative results can not be excluded.
3. During certain antibiotic treatments, the concentration of H. pylori antigens may decrease below the detection limit of the test.
4. If the test result is negative but clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of an H. pylori infection.


Table 1. Performance / Tabelle 1. Leistung / Tableau 1. Performances / Tabella 1. Prestazioni / Tabla 1. Rendimiento / Tabela 1. Desempenho / Tabela 1. Oziatolania.

Table with 3 columns: HELICOBACTER PYLORI ANTIGEN RAPID TEST, Positive, Negative. Shows results for 100 samples.

Specificity = 96.6 % (CI 95%: 91.8 - 100%) Sensitivity = 95.2% (CI 95%: 88.8 - 100%) Accuracy = 96% (CI 95%: 92.1 - 99.8%)

REFERENCES / BIBLIOGRAFIA / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA

- 1. Shimogama T. Stool antigen tests for the management of Helicobacter pylori infection. World Journal of Gastroenterology. 2013; 19(45):8188-91.
2. Ana Isabel Lopes, Filipa F Vale, Mónica Oleastro. Helicobacter pylori infection - recent developments in diagnosis. World Journal of Gastroenterology. 2014 July 28; 20(28): 9299-9313.
3. Andreas Mentis, Philippe Lehours, and Francis Mégraud. Epidemiology and Diagnosis of Helicobacter pylori infection. Helicobacter. 2015 Sep;20 Suppl 1:1-7.
4. Zagari RM, Rabitti S, Eusebi LH, Bazzoli F. Treatment of Helicobacter pylori infection: A clinical practice update. Eur J Clin Invest. 2018 Jan;48(1). doi: 10.1111/eci.12857.

SYMBOLS / SYMBOLE / SYMBOLES / SIMBOLOGIA / SIMBOLOGIA / SIMBOLOGIA / SIMBOLE

Legend for symbols: IVD (In vitro diagnostic device), Read the instructions before use, Expiry date, Legal manufacturer, CE marking, Temperature limits, Do not reuse, Sufficient for tests, REF (List number), LOT (Lot number), EC REP (Authorized Representative), Importer.

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