Your partner for rapid testing

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Table with 4 columns: RESULTS, Negative, Positive, TOTAL. Rows for Negative and Positive results with sub-rows for different concentration ranges.

Relative Specificity = 96.0% (CI 95%:86.5%–98.9%)

Relative Sensitivity = 98.7% (CI 95%: 93.0%–99.8%)

Relative Accuracy = 97.6% (CI\*: 93.3%–99.2%)

Table 2. Interfering substances / Tabelle 2. Beeinträchtigende Substanzen / Tableau 2. Substances interférentes / Tabella 2. Sostanze interferenti / Tabla 2. Sustancias que interfieren / Tabela 2. Substâncias interferentes / Tabela 2. Substancje zakłócające

Table with 2 columns: SUBSTANCE, CONCENTRATION ASSAYED. Lists various substances like Bilirubin, Cholesterol, Triglycerides, etc.

SYMBOLS / SYMBOLE / SIMBOLES / SIMBOLOGIA / SIMBOLS / SIMBOLGIA / SIMBOLE

Icons and text for various symbols: IVD, CE, REF, LOT, etc., with their respective meanings in multiple languages.

REFERENCES / BIBLIOGRAFIE / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA

- 1. C-reactive protein concentrations as a marker of inflammation or infection for interpreting biomarkers of micronutrient status. Vitamin and Mineral Nutrition Information System. Geneva: World Health Organization; 2014. (WHO/NMH/NHD/EPG/14.7).

C-REACTIVE PROTEIN RAPID TEST


Rapid test for the semi-quantitative detection of C-Reactive Protein (CRP) in whole blood or serum samples. Schnelltest für den semi-quantitativen Nachweis von C-reaktivem Protein (CRP) in Vollblut- oder Serumproben.

LIMITATIONS

- 1. As with any diagnostic method, the physician must evaluate the test result in light of other clinical evidence.
2. People suffering from haemophilia or taking anticoagulants could have an excessive bleeding after the fingerprick.

CONTENT

- 20 hermetically sealed aluminum bags containing 1 C-REACTIVE PROTEIN TEST cassette and 1 desiccant bag;
• 20 vials with dropper tip containing enough C-REACTIVE PROTEIN TEST diluent for 20 tests;

EN - INSTRUCTIONS FOR USE

C-REACTIVE PROTEIN

C-Reactive Protein (CRP) is an acute phase protein produced mainly by the liver, whose concentration increases following injury, infection and inflammation. Although it is not a specific marker for a particular disease, CRP is currently used as a generic diagnostic indicator of infection and inflammation, as well as to monitor patients' reactions to treatment and the post-operative course.

PRINCIPLE OF THE TEST

C-REACTIVE PROTEIN RAPID TEST is a lateral flow test device based upon a principle of solid-phase immunochromatography (also called lateral flow technology) employing a sandwich technique. The detection of specific CRP levels takes place on a test strip embedded in a plastic cassette. The methodology implied is a direct sandwich EIA test also named G.A.I.A. (Gold adsorbed Immuno Assay).

PRECAUTIONS

- 1. Read these instructions for use carefully before performing the test. The test is reliable if the instructions are followed correctly (reaction times and how the sample is collected).
2. Do not open the sealed aluminum pouch until just before performing the test. The desiccant bag must not be used.

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POSITIVE, C-REACTIVE PROTEIN BETWEEN 40 AND 80 mg/L

Three coloured lines appear in the reading window, one at the C sign and two more at the T1 and T2 signs, as shown in the figure opposite. This result may be associated with a viral or bacterial infection or physical trauma.

POSITIVE, C-REACTIVE PROTEIN OVER 80 mg/L

Four coloured lines appear in the reading window, one at sign C plus three more at signs T1, T2 and T3 as shown in the figure opposite. This result can be associated with a bacterial infection or severe inflammation.

NEGATIVE, C-REACTIVE PROTEIN BETWEEN BELOW 10 mg/L

A coloured line appears in the Control (C) line region. No line appears in the Test (T1, T2, T3) regions. This result indicates that the concentration of CRP is below 10 mg/L and suggests the absence of infection and inflammation.

INVALID RESULT

The Control (C) line does not appear. An insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the Control line (invalid result). Review the procedure and repeat the test with a new device and a new sample.

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1) Prick the finger with the lancet and obtain a drop of blood. Fill the capillary pipette to the black line - fig. A1.
2) Transfer 10 µl of whole blood into the vial previously opened. Press the bulb 2 or 3 times to ensure that all the blood has moved into the diluent solution. - fig. B1.

SENUM OR PLASMA COLLECTION PROCEDURE

- 1) Using a laboratory micropipette, collect 6 µl of serum and transfer to the vial previously opened. - fig. B3
2) Replace the screw cap onto the vial and mix well for at least 10 seconds. Unscrew the blue cap from the dropper vial (leave the white cap tightly screwed on). Deposit 3 drops into the well shown on the cassette (S). - fig. C1

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