PRIMA®

PRIMA Lab SA
Viale Serfontana, 10
CH-6834 Morbio Inferiore - SWITZERLAND
support@primalabsa.ch
primalabsa.ch

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Table with 3 columns: BOWEL FOB RAPID TEST, Positive, Negative. Rows show specificity and sensitivity data for positive and negative results.

Specificity = 99 % (CI 95%: 97% - 100%)
Sensitivity = 78.7 % (CI 95%: 67% - 90.4%)
Accuracy = 92.6% (CI 95%: 88.4 - 96.8%)

REFERENCES / BIBLIOGRAFIA / BIBLIOGRAPHIE / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA

- 1. Van Rossum Leo G. et al.: "Random Comparison of Guaiac and Immunochemical Fecal Occult Blood Tests for Colorectal Cancer in a Screening Population"
2. Faivre J., Dancourt V., Denis B., Darval E., Plette C., Perrin P., Bidan J.M., Jard C., Jung S., Levillain R., Viguer J., Bretagne J.F.: "Comparison between a guaiac and three immunochemical fecal occult blood tests in screening for colorectal cancer."
3. Benton S.C., Seaman H.E., Halloran S.P.: "Faecal Occult Blood Testing for Colorectal Cancer Screening: the past or the Future."
4. European Commission. European guidelines for quality assurance in colorectal cancer screening and diagnosis.

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diluent for 20 tests;

• 1 instructions for use leaflet.
Material required but not supplied: a device for measuring time (e.g. stopwatch, clock).

SAMPLE COLLECTION AND PREPARATION

Allow the test, sample and buffer to reach room temperature (15-30°C) before testing.
1) Unscrew the top of the sample collection device containing the sample collection probe.
Note: Avoid dilution of stool specimens with urine or water from the toilet bowl. Patients should not collect samples during their menstrual period or 3 days before and after it, or if they have bleeding hemorrhoids, blood in urine or if they experience strain during their bowel movement. Alcohol, aspirin and other medications taken in excess may cause gastrointestinal irritation, resulting in occult bleeding. Such substances should be discontinued at least 48 hours prior to testing. No dietary restrictions are necessary before testing. Perform the test as soon as possible after sample collection.

2) Collect the stool sample with the top of the collection device by plunging it into 3 different places of the same stool sample -fig. A and put it in the collection vial.
Do NOT put the collection stick back into the tube in between collecting the 3 samples.
3) Screw on the cap again, immersing the stick in the stool extraction liquid, and shake the vial for at least 10 seconds. -fig. B
4) Proceed with test procedure.

TEST PROCEDURE

5) Open the aluminium bag, remove the test cassette and place it on a clean and level surface. Throw away the desiccant bag.
6) Break the end of the red cap, exposing the dropper portion. -fig. C
7) Dispense 3 drops of the diluted stool sample in the sample well (S) shown on the test cassette. Wait 1-2 seconds between dispensing one drop of diluted stool sample. -fig. D
8) Wait 5 minutes.

RESULTS INTERPRETATION

READ THE RESULTS AT EXACTLY 5 MINUTES. The intensity of the line colors is not relevant for the purposes of interpretation of the Test's results.


POSITIVE RESULT

Two coloured lines appear in the reading window by the Test (T) and Control (C) signs. The T line may be less intense (lighter) than the C one. This result means there is blood in the faeces.

NEGATIVE RESULT

A coloured line appears in the Control (C) region. No line appears in the Test (T) region. This result means that there is no blood in the faeces or that its concentration is below the test's detection limit (40 ng/mL).

INVALID RESULT

The Control (C) line does not appear. An insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the Control line (invalid result). Review the procedure and repeat the test with a new device and a new sample.

QUALITY CONTROL

An internal procedural control is included in the test cassette. A coloured line appearing in the Control line region (C) is considered an internal positive procedural control, confirming sufficient specimen volume and correct procedural technique.

PERFORMANCE CHARACTERISTICS

To assess the performance of BOWEL FOB RAPID TEST, a clinical evaluation was performed by comparing the results obtained against Laboratory Reference Method. 148 samples were tested with both methods (Table 1).

Specificity = 99 % (CI 95%: 97 - 100%)
Sensitivity = 78.7 % (CI 95%: 67 - 90.4%)
Accuracy = 92.6% (CI 95%: 88.4 - 96.8%)

MEASURING RANGE

The study concerning the Limit Of Detection (LOD) was conducted with the Human Haemoglobin Calibrator SENTINEL CH. The LOD of BOWEL FOB RAPID TEST is 40 ng/mL. No adverse effect on T-line formation (prozone effect) was observed for FOB value up to 3 mg/mL.

CROSS-REACTIVITY

The following compounds were tested with BOWEL FOB RAPID TEST and no interference was observed. Table 2. Cross-reactants tested.

Table with 2 columns: SUBSTANCE, CONCENTRATION ASSAYED. Rows include Bilirubin, Albumin, Glucose, Urea, Uric acid, and Bovine - pig - rabbit - horse - chicken - sheep Haemoglobin.

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