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
Huachengyang (Shenzhen) Technology Co., Ltd
Qianhai 1st Building A, 128th Shenzhen East Rd, Huangpu Community, Xingyao St, Bao'an, Shenzhen, Guangdong, People's Republic of China
R Sicht B.V. - Roald Dahllaan 47 - 5629 MC Eindhoven - Netherlands

Sterile tongue depressor / Steriler Zungenspatel / Abaisse-langue stérile / Abassalingua stérile / Depressor de lingua estéril / Abaxador de lingua estéril / Sterylny despoter jazyka

Applimed SA
ZI, Route de Pra de Plan 1, CH-1618 Châtel-Saint-Denis - Switzerland

Dansu A/S Sandbakken 1-3, Gårlev, DK-3660 Stenløse, Denmark

Jiangsu Suyun Medical Materials Co., Ltd
No.18 Jin Qiao Road Dapu Industrial Park
222002 Lianyungang, Jiangsu Province, People's Republic of China
Shanghai International Trading Corp. GmbH (Europe)
Erlterstrasse, 80, D-20537 - Hamburg, Germany

Your partner for rapid testing

IVD for professional use only
IVD nur für den professionellen Gebrauch
DIV à usage professionnel uniquement
IVD solo per uso professionale
IVD apenas para uso profesional
IVD wyłącznie do użytku profesjonalnego

REFERENCES / BIBLIOGRAFIE / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA / BIBLIOGRAFIA

- Anjos LM, Marcondes MB, Lima MF, Mondelli AL, Okoshi MP. Streptococcal acute pharyngitis. Rev Soc Bras Med Trop 2016; 47(4):409-13.
- Castro SA, Dorfmueller HC. A brief review on Group A Streptococcus pathogenesis and vaccine development. R Soc Open Sci 2022.
- Banerjee S, Ford C. Rapid Tests for the Diagnosis of Group A Streptococcal Infection: A Review of Diagnostic Test Accuracy, Clinical Utility, Safety, and Cost-Effectiveness. Ottawa (ON), Canadian Agency for Drugs and Technologies in Health; 2018 May. CADTH Rapid Response Reports.
- Llor C, Bjerrum L, Munkáçs A, Cots JM, Hernández S, Maragás A. HAPPY AUDIT Investigators. Access to point-of-care tests reduces the prescription of antibiotics among antibiotic-requesting subjects with respiratory tract infections. Respir Care 2014; 59(12):1911-1917.
- Georgers N, Nakhouli MD and John Hickner MD, MSc. Management of Adults with Acute Streptococcal Pharyngitis: Minimal Value for Backstop Strep Testing and Overage of Antibiotics. J Gen Intern Med 2013 Jun; 28(6):830-834.
- Llor C, Bjerrum L, Munkáçs A, Cots JM, Hernández S, Maragás A. HAPPY AUDIT Investigators. Access to point-of-care tests reduces the prescription of antibiotics among antibiotic-requesting subjects with respiratory tract infections. Respir Care 2014; 59(12):1911-1917.
- Georgers N, Nakhouli MD and John Hickner MD, MSc. Management of Adults with Acute Streptococcal Pharyngitis: A Systematic Review and Meta-Analysis. PLoS One 2014; 9(11):1-11.
- Parks T, Barrett L, Jones N. Invasive streptococcal disease: a review for clinicians. British Medical Bulletin. Volume 115, Issue 1, September 2015, Pages 77-89.

Table 1. Performance / Tabelle 1. Leistung / Tableau 1. Performances / Tabella 1. Prestazioni / Tabla 1. Rendimiento / Tabella 1. Desempenho / Tabela 1. Oztarimio

STREP A RAPID TEST	REFERENCE METHOD	
	POSITIVE	NEGATIVE
POSITIVE	34	36
NEGATIVE	3	199
	37	201

SPECIFICITY = (199/201)*100 = 99.00% (CI 95%: 97.63 - 100%)
SENSITIVITY = (34/37)*100 = 91.89% (CI 95%: 78.70 - 97.20%)
ACCURACY = ((233/238)*100 = 97.90% (CI 95%: 96.08 - 99.72%))

Table 2. Interfering substances / Tabelle 2. Beeinträchtigende Substanzen / Tabella 2. Substances interferentes / Tabella 2. Sostanze interferenti / Tabla 2. Sustancias que interfieren / Tabla 2. Substâncias interferentes / Tabela 2. Substâncias zakkócające

BACTERIA	ATCC	LOT	CONCENTRATION
<i>Bordetella pertussis</i>	9340	BP-02	ASSAYED
<i>Campylobacter jejuni</i>	29428	CJ-02	10 ⁷ CFU/mL
<i>Candida albicans</i>	MUT6298	CA-01	10 ⁷ CFU/mL
<i>Clostridium perfringens</i>	13124	CP-02	10 ⁷ CFU/mL
<i>Corynebacterium diphtheriae</i>	11913	CD-02	10 ⁷ CFU/mL
<i>Enterobacter cloacae</i>	MUT6313	EC-02	10 ⁷ CFU/mL
<i>Helicobacter pylori</i>	43504	HP-01	10 ⁷ CFU/mL
<i>Hemophilus influenzae</i>	49766	HI-01	10 ⁷ CFU/mL
<i>Klebsiella pneumoniae</i>	MUT6314	KP-02	10 ⁷ CFU/mL
<i>Proteus vulgaris</i>	MUT6312	PV-02	10 ⁷ CFU/mL
<i>Pseudomonas aeruginosa</i>	MUT6310	PA-02	10 ⁷ CFU/mL
<i>Shigella flexneri</i>	12022	SF-02	10 ⁷ CFU/mL
<i>Staphylococcus aureus</i>	MUT6311	SA-02	10 ⁷ CFU/mL
<i>Streptococcus agalactiae</i>	12386	SAG-02	10 ⁷ CFU/mL
<i>Streptococcus dysgalactiae</i>	12388	SD-02	10 ⁷ CFU/mL
<i>Streptococcus gallolyticus</i>	9808	SG-02	10 ⁷ CFU/mL
<i>Yersinia enterocolitica</i>	51871	YE-02	10 ⁷ CFU/mL

EN - INSTRUCTIONS FOR USE STREPTOCOCCUS A

The throat infection with group A β-hemolytic Streptococcus (also known as *Streptococcus pyogenes*) is a Gram-positive bacterium that constitutes the most frequent bacterial cause of acute pharyngitis and it is also responsible for some skin infections such as impetigo and erysipelas. Most sore throats, on the other hand, are caused by viral infections, which are cured without the need of antibiotic treatments, unnecessary in these cases. Streptococcus A infections can resolve in a few days without treatment, although doctors generally prefer to prescribe antibiotics to prevent infection-related complications that can also be serious, such as rheumatic fever or acute glomerulonephritis. Group A Streptococcus continues to be a focus of interest not only because of their causal role in acute streptococcal pharyngitis and other pyogenic infections, but also because of its association with poststreptococcal sequelae, specifically, acute rheumatic fever and acute glomerulonephritis. In order to properly treat the disease using antibiotic therapy, it is important to use an accurate diagnostic method to identify the pathological agent. For the screening of Group A Streptococcus infection several methods are currently used including Bacitracin sensitivity tests on a sheep blood agar plate, latex agglutination and enzyme immunoassay. STREP A RAPID TEST is a rapid immunological test for the qualitative and visual detection of group A Streptococcus bacterium in throat infections. The test helps to quickly identify whether a sore throat is caused by Strep A rather than by other pathogens (usually viruses) which may not require antibiotic treatment.

PRINCIPLE OF THE TEST
STREP A RAPID TEST detects specific Streptococcus A antigens on the sides of the test. The test is performed by dipping a test strip. The specimen, which is pharyngeal swab, is taken with a disposable sterile rayon swab. The specimen is placed in a plastic collection tube where the user dispenses a specific volume of extraction reagents

STREP A RAPID TEST

Rapid Test for the detection of group A β-hemolytic Streptococcus in pharyngeal swab
Schnelltest zum Nachweis von β-Hämolysierenden Streptokokken der Gruppe A in Rachenabstrichen
Test rapide pour la détection du Streptococcus β-hémolytique sur un tampon pharyngé
Test rapido per la rilevazione dello Streptococco β-emolitico di gruppo A em hisopo faringeo
Prueba rápida para la detección de Estreptococo β-emolitico del grupo A em hisopo faringeo
Teste rápido para detecção de Estreptococos β-hemolíticos do grupo A em swab faringeo
Szybki test do wykrywania β-hemolizujących Streptococcus grupy A w wymazie z gardła

(in order to extract specific Streptococcus A antigens). After an incubation time, the obtained sample-solution is dispensed onto the test cassette mainly with the provided reagents are:
- R1 solution mainly containing Sodium Nitrite and Phenol Red.
- R2 solution mainly containing Acetic Acid.
Once mixed, the two solutions create Nitrous Acid, acting as Streptococcus A specific antigen extractor. In case of a Strep A concentration equal or higher than 10⁷ CFU/mL in the patient's sample, the specific antigens present in the sample solution are firstly bound by the anti-Strep A antibodies labeled with gold particles. The immunocomplexes are selectively captured by the test line (T line, made of purified anti-Strep A antibodies), thus showing a red-purple band appearance in correspondence of the T line, pointing out a positive result. In case of absence of Strep A or when its levels in patient's sample are lower than the cut-off value (10⁷ CFU/mL), not enough specific immunocomplex is formed. Consequently, only one band, the C line, is visible indicating a negative result. To serve as a procedural control, a colored line will always appear in the Control line region, indicating that the test procedure was being properly performed and that the test components and reagents have operated as intended.

- PRECAUTIONS**
- Read these instructions for use carefully before performing the test. The test is reliable if the instructions are followed correctly (reaction times, how the sample is collected and the extraction instructions shown).
 - Humidity and temperature can adversely affect the results. DO NOT FREEZE.
 - Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards during all procedures and follow standard procedures for proper disposal of specimens.
 - Do not use the kit after the expiry date or if the packaging and/or internal components are damaged.
 - The test is for external use only. Do not swallow solutions 1 and 2 and the other items in the box.
 - If liquids come into contact with the eyes, skin or other mucosae, rinse immediately and thoroughly with water.
 - Do not switch the caps of the reagent vials or mix components from different test kits.
 - Not suitable for use on children under the age of 2 years.
 - All kit components are single-use, with the exception of the workstation, instructions for use and solutions 1 and 2, which must be stored correctly closed with their caps.
 - Do not open the sealed elements (test, swab and tongue depressor) until just before performing the test. Take care to open as marked. Dispose of the desiccant sachet inside the aluminium pouch, without opening it.
 - Do not use the kit after the expiry date or if the packaging and/or internal components are damaged.
 - The test is for external use only. Do not swallow solutions 1 and 2 and the other items in the box.
 - If liquids come into contact with the eyes, skin or other mucosae, rinse immediately and thoroughly with water.
 - Do not switch the caps of the reagent vials or mix components from different test kits.
 - Not suitable for use on children under the age of 2 years.
 - All kit components are single-use, with the exception of the workstation, instructions for use and solutions 1 and 2, which must be stored correctly closed with their caps.
 - Do not open the sealed elements (test, swab and tongue depressor) until just before performing the test. Take care to open as marked. Dispose of the desiccant sachet inside the aluminium pouch, without opening it.

CONTENT

- 20 hermetically sealed aluminium pouches each containing 1 STREP A TEST cassette and 1 desiccant bag;
- 20 sterile swabs;
- 20 sterile wooden tongue depressors;
- 1 workstation;
- 20 empty plastic extraction tubes with dropper;
- 1 dropper vial containing enough STREP A TEST-R1 solution to perform 20 tests. R1 contains 20.7% Sodium Nitrite, H302: Harmful in case of ingestion.
- 1 dropper vial containing enough STREP A TEST-R2 solution to perform 20 tests. R2 contains 0.55% Acetic Acid;
- 1 Instructions for use leaflet.
- Required material not supplied: a timer (e.g. timer, watch).

SPECIMEN COLLECTION AND PREPARATION
Bring tests, specimens, reagents to room temperature (15-30°C) prior to testing. To avoid cross contamination, do not allow all the tips of the reagent bottles to come into contact with specimen material.

- Place the workstation on a clean, dry, flat surface. Place the empty plastic test tube in one of the compartments of the workstation - **fig. C**
- Collecting the sample: Open the package containing the tongue depressor, leaving it inside the package and next to you. Open the package containing the sterile swab. Avoid touching the cotton tip and remove the swab using the plastic handle. Tilt the patient's head backwards and let open the mouth as much as possible. Pick up the tongue depressor with one hand and flatten the tongue. Use the other hand to place the swab close to the throat. Touch the back of the throat - the area around the tonsils and any redness or painful part (soft palate, uvula) with the cotton tip. Rotating the swab is recommended because it increases the amount of collected sample - **fig. B**

It is recommended that swab specimens be processed as soon as possible after collection.

TEST PROCEDURE

- After collecting the sample, insert the swab's cotton tip into the supplied plastic test tube that you had previously placed in the workstation - **fig. C**
- Open the vial STREP A TEST - R1 by unscrewing only the blue cap (leave the white part firmly sealed), then add 4 drops into the test tube and close it with the cap - **fig. D**
- Next, open the vial STREP A TEST - R2 by unscrewing only the green cap (leave the white part firmly sealed), then add 4 drops into the test tube and close it with the cap - **fig. E**
WARNING: the addition of R2 to R1 leads to the formation of a solution for which the indications described in point 8 and 9 of the "precautions" apply.
- Holding the swab's plastic handle, rotate the swab against the sides of the test in 10 to 15 times to thoroughly mix the solution. Leave the swab to incubate for 2.5 minutes - **fig. F**
- At the end of the incubation period, remove the test tube (with the swab and extraction fluid still in it) from the workstation. Using your thumb and index finger, press the sides of the test

tube to release as much fluid as possible from the swab's cotton tip and collect it in the test tube. Remove the swab - **fig. G**. Dispose of the swab in compliance with guidelines for the handling of infectious agents and put the test tube back into one of the compartments in the workstation.

- Add a supplied dropper to the plastic test case - **fig. H**
- Dispense 1 drop into the well in the test case - **fig. I**
NOTE: If the dispensed drop contains air bubbles, add another drop to the well.
- Read the results after 5 minutes** (some positive results can be seen after just 1 minute). Do not read after 10 minutes.

RESULTS INTERPRETATION
READ THE RESULTS AFTER 5 MINUTES. DO NOT READ THE RESULTS AFTER 10 MINUTES. The intensity of the line colors is not relevant for the purposes of interpretation of the test's results.

POSITIVE RESULT
Two coloured lines appear in the reading window by the Test (T) and Control (C) signs. The T line may be less intense (lighter) than the C one. This result means that the test has detected the presence of Streptococcus A specific antigens in the sample.

NEGATIVE RESULT
A coloured line appears in the Control (C) region. No line appears in the Test (T) region. This result is negative and means that no Streptococcus A antigens have been detected in the sample from the pharynx or that their amount is too low to be detected.

INVALID RESULT
The Control (C) line does not appear. An insufficient sample volume or incorrect procedural techniques are the most likely reasons for the absence of the Control line (invalid result). Review the procedure and repeat the test with a new device and a new sample.

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