## Test your health

<span></span>	<b>PRIMA LAB SA</b> <div>Via Antonio Monti 7 CH-6828 Balerna - SWITZERLAND</div>	<span></span>	<b>Qbd RepS BV</b> <div>Groenenborgerlaan 16 2610 Wilrijk - Belgium</div>
<span></span>	<b>support@primalabsa.ch</b> primalabsa.ch	<span></span>	
<span></span>	<b>820027_IFU_11_8.0</b> 12/2024	<span></span>	<b>820027-x</b>

Registered trademarks.

### ENGLISH WHAT IS ALBUMIN?

Albumin is the most abundant protein in the plasma and it is synthesized by the liver. It has numerous functions - it carries metabolites (such as bilirubin, free fatty acids, thyroid hormones), it maintains oncotic pressure (which is essential for the distribution of liquids in blood capillaries and tissues) and it can be used as an amino acid reserve. In normal conditions, the quantity of albumin excreted via the urine is very low. Abnormally high levels of albumin in the urine mean there is pathological albuminuria, which can be caused by conditions such as hypertension, diabetes, kidney dysfunction, amyloidosis, systemic lupus erythematosus. Other factors that can affect albumin levels in the urine include smoking, dehydration, physical exercise, urinary tract infections, heart failure, hypertriglyceridemia and obesity. Kidneys Albumin Test is a device that allows to measure the level of albumin in the urine.

### REFERENCES

- Urinary biomarkers of kidney dysfunction*. *Ann Biol Clin (Paris)*. 2015 Mar-Apr;73(2):151-7. doi: 10.1684/abc.2015.1029.
- Lopez-Giacoman, Salvador and Magdalena Madero. "Biomarkers in chronic kidney disease, from kidney function to kidney damage" World journal of nephrology* vol. 4,1 (2015): 57-73.
- Sacks, David B et al. "Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus" Diabetes care* vol. 34,6 (2011): e61-99.

#### FA.Q. – QUESTIONS AND ANSWERS

**HOW DOES THE KIDNEYS ALBUMIN TEST WORK?**

Kidneys Albumin Test specifically identifies the proteins in human urine and shows if the concentration is equal to or higher than 10 µg/ml. In some cases, samples with albumin concentration above this value (between 11µg/ml and 20 µg/ml) can yield a negative result.

**WHEN CAN THE TEST BE USED?**

The Kidneys Albumin Test is recommended if there are symptoms related to the urinary system or if hyperglycemia or hypertriglyceridemia are suspected.

**CAN THE RESULT BE INCORRECT?**

The result is correct if the instructions are followed carefully. However, there may be incorrect results if other liquids end up on the device before the test is carried out or if the number of drops dispensed into the well is incorrect; or if the result reading time is not respected.

**HOW DO I INTERPRET THE TEST IF THE COLOUR AND INTENSITY OF THE TEST AND CONTROL BANDS ARE DIFFERENT?**

The colour and the intensity of the bands are not important for the interpretation of the result. The bands must be full and homogeneous. The test is negative regardless of the intensity of the colour of the test band.

**IS THE RESULT RELIABLE IF READ AFTER 5 MINUTES?**

No. The test must be read exactly 5 minutes after the procedure is completed and no later than 6 minutes from this time.

**WHAT SHOULD I DO IF THE RESULT IS POSITIVE?**

A positive result means that the level of albumin detected in the urine is higher than normal. You should consult a physician and show him/her the results of the test obtained after exactly 5 minutes. The physician will decide which other tests to perform or refer you to a consultant.

**WHAT SHOULD I DO IF THE RESULT IS NEGATIVE?**

A negative result means that the level of albumin detected in the urine is lower than normal. If symptoms persist, consult your physician.

**CAN THE TEST BE PERFORMED WHILE I AM PREGNANT?**

Yes. If positive, contact your physician.

**CAN THE RESULTS BE POSITIVE IN NORMAL CONDITIONS?**

Yes. Albumin levels in the urine may increase after intense physical exercise or when you have a high fever. If positive, contact your physician.

**HOW ACCURATE IS THE KIDNEYS ALBUMIN TEST?**

The Test is very accurate. Assessment reports show a concordance rate of 96.97% (CI 95%: 91.47 - 98.96%) with reference methods.

<span></span>	<b>PRIMA LAB SA</b> <div>Via Antonio Monti 7 CH-6828 Balerna - SWITZERLAND</div>	<span></span>	<b>Qbd RepS BV</b> <div>Groenenborgerlaan 16 2610 Wilrijk - Belgium</div>
<span></span>	<b>support@primalabsa.ch</b> primalabsa.ch	<span></span>	
<span></span>	<b>820027_IFU_11_8.0</b> 12/2024	<span></span>	<b>820027-x</b>

#### PRECAUTIONS

- Read these instructions for use carefully before performing the test. The Test is reliable only if all the instructions are followed correctly.
- Keep the Test out of the reach of children.
- Do not use the Test after the expiry date or if the package has been damaged.
- Follow the indicated procedure and time strictly.
- Store the Test components at a temperature between +4°C and +30°C. Do not freeze.
- Use the test once only.
- The test is for external use only. **DO NOT SWALLOW.**
- In vitro diagnostic device for individual use.
- After using, please dispose of all components according to your local waste disposal laws. Ask your pharmacist for advice.

INSTRUCTIONS FOR USE  
GEBRAUCHSANWEISUNG  
NOTICE D'UTILISATION  
ISTRUZIONI PER L'USO

<span></span>	<b>PRIMA HOME TEST</b>	<span></span>	<b>Qbd RepS BV</b> <div>Groenenborgerlaan 16 2610 Wilrijk - Belgium</div>
<span></span>	<b>Self-test for detecting albumin in the urine</b> <i>Selbsttest zur Bestimmung von Albumin im Urin</i> <b>Autotest pour la détermination de l'albumine dans les urines</b> <i>Self-test per la determinazione dell'albumina nelle urine</i>	<span></span>	
<span></span>	<b>820027_IFU_11_8.0</b> 12/2024	<span></span>	<b>820027-x</b>

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



**ENGLISH**


**CONTENT**

\* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the external box.

- \* hermetically sealed aluminium pouch containing: 1 KIDNEYS ALBUMIN TEST cassette, 1 desiccant bag

Do not open the sealed aluminium bag until just before performing the test. Take care to open it as marked. The desiccant packet must not be used. Dispose of it with household waste without opening it.

- \* pipette

- 1 instructions for use leaflet

**Material required but not supplied:**

- a device to measure time (i.e. timer, watch)
- clean, dry, residue-free container.

**TEST PROCEDURE**

- 1) Wash hands with soap and warm water, rinse with clean water and allow to dry.
- 2) Prepare the necessary material as follows: open the aluminium pouch, take out only the test cassette and throw away the desiccant bag. –**FIG. A**
- 3) Collect urine in a clean clean, dry, residue-free container (not supplied) –**FIG. B**
- 4) Press the pipette bulb and dip the pipette tip into the collected urine. Release the bulb to collect the urine sample. –**FIG. C**
- 5) Gently press the pipette bulb and deposit 3 drops of sample into the well (S) indicated on the cassette, wait 1-2 seconds between dispensing one drop and the next. –**FIG. D**
- 6) Wait 5 minutes and read the results as follows.

**RESULTS INTERPRETATION**

**READ THE RESULT AFTER EXACTLY 5 MINUTES**

The intensity of the line colours is not relevant for the purposes of interpretation of the Test's results.

**NEGATIVE RESULT**

Two coloured bands appear in the reading window by the T (test) and C (control) signs. The T band may be less intense (lighter) than the C line.

This result means that the levels of albumin in the urine are normal or lower than the value used to calibrate the test (10 µg/ml).

**POSITIVE RESULT**

A coloured band appears only under the C (control) sign.

This result means the levels of albumin in the urine are higher than normal.

**NON VALID RESULT**

No bands appear or there is a line only under the T (test) sign and not under the C (control) sign.

In this case it is not possible to interpret the result of the test, which must be considered not valid. Repeat the test.

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