
Test your health

PRIMA LAB SA
Viale Serfontana, 10
CH-6834 Morbio Inferiore
SWITZERLAND
support@primalabs.ch
primalabs.ch

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CE2H88000KIT01ES
CE2H88000BAG01ED

REVISIONS

Table with 3 columns: REVISION, DATE, CHANGES. Row 1: 2.0, 04/2026, New registered office of PRIMA Lab SA, and new REF added

BIBLIOGRAPHY / LITERATURVERZEICHNIS / RÉFÉRENCES / BIBLIOGRAFIA / BIBLIOGRAFIE

- 1. Lindfors, K., Ciacci, C., Kurppa, K. et al. Coeliac disease. Nat Rev Dis Primers 5, 3 (2019).
2. Al-Toma et al. European Society for the Study of Coeliac Disease (EUSC) guideline for coeliac disease and other gluten-related disorders. United European Gastroenterol J. 2019 Jun;7(5):583-613.
3. King, et al. 2020. Incidence of Coeliac Disease Is Increasing Over Time: A Systematic Review and Meta-Analysis. The American Journal of Gastroenterology 115 (4): 507-25.
4. Perrinello G., Da Re M., Grizzo F., Carnelli S., Cozzi M., Marinoni F., Villalta D. Diagnostic accuracy of a novel point-of-care test for simultaneous detection of anti-transglutaminase IgA and anti-deamidated gliadin IgG antibodies. J Clin Lab Anal. 2024

SYMBOLS / SYMBOLE / SYMBOLES / SIMBOLOGIA / SIMBOLEEN

IVD: In vitro diagnostic device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo médico-diagnóstico in vitro / In-vitrodiaagnostisch hulpmiddel
STERILE R: Read the instructions before use / Beachten Sie die Gebrauchsanweisung / Consulter le mode d'emploi / Leggere attentamente le istruzioni per l'uso / Lees de instructies voor gebruik
i: Read the instructions before use / Beachten Sie die Gebrauchsanweisung / Consulter le mode d'emploi / Leggere attentamente le istruzioni per l'uso / Lees de instructies voor gebruik
Expiry date (last day of the month) / Verwendbar bis (letzter Tag des Monats) / Utiliser jusqu'à (dernier jour du mois) / Utilizzare entro (ultimo giorno del mese) / Vervaldatum (laatste dag van de maand)
Legal manufacturer / Hersteller / Fabricant / Fabricante / Juridische fabrikant
CE marking / CE-Kennzeichnung / Marquage CE / Marchio CE / CE-markering
Temperature limits / Temperaturbegrenzung / Limites de température / Limiti di temperatura / Temperaturlimiten
Do not reuse / Nicht wiederverwenden / Ne pas réutiliser / Non riutilizzare / Niet opnieuw
Sufficient for <= tests / Ausreichend für <= Tests / Suffisant pour <= tests / Sufficiente per <= test / Adequado para <= tests
List number / Katalognummer / Code produit / Codice prodotto / Lijstnummer
Lot number / Chargenbezeichnung / Code du lot / Lotto prodotto / Lotnummer
Authorized Representative in the European Community / Bevollmächtigter in der Europäischen Gemeinschaft / Mandatario nella Comunità Europea / Gemachtigde vertegenwoordiger in de Europese Gemeenschap
Unique Device Identification / Eindeutige Geräteidentifikation / Identification unique du dispositif / Identificativo unico del dispositivo / Unieke apparatidentificatie
Importer / Importateur / Importatore / Importeur
Intended for self-testing / Für den Selbsttest bestimmt / Destiné à l'auto-test / Destinato all'auto-test / Bestemd voor zelftest
Not for near patient testing / Produkt nicht für patientennahe Tests / Dispositivo non destiné au diagnostic près du patient / Dispositivo non destinato ad analisi decentralizzate / Geen hulpmiddel voor patientnahe testen

MEDICAL DEVICES INSIDE THE KIT / MEDICINISCHES GERÄTE INNERHALB DES KITS / DISPOSITIVOS MÉDICOS A L'INTÉRIEUR DU KIT / DISPOSITIVOS MEDICI INTERNI AL KIT / MEDISCHE HULPMIDDELEN IN DE KIT
STERILE R: Lancet / Lanette / Bloc autopiqueur / Lancetta / Lancet
Hualan Trianda Medical Instruments Co., Ltd.
No.116 East Songjiang Road, Hualin, Economic&Technological Development Zone, 223300 - Hualan City, Jiangsu - P.R. China
Riomaxiv Sociedad Limitada
Calle de Almansa 55, 1D, Madrid 28039, Spain
Alcohol cleansing gauze / Alkohol-Reinigungsgaze / Gaze nettoyante à l'alcool / Garza detergente alcolica / Reinigungs-gas mit alcohol
Pharma Innovative Healthcare Manufacturers Pvt. Ltd.
E-209, SHI Mahape Road, Electronic Zone, MIDC, TTC Industrial Area, Mahape, Navi Mumbai 400 710 MH India
Advena Ltd., Tower Business Centre,
2nd Flr. Tower Street, Swatara, BKR 4013, Malta

CELIAC DISEASE SCREENING TEST

Rapid immunochromatographic manual self-test for the qualitative detection of anti-deamidated gliadin peptide IgG and anti-tissue transglutaminase IgA in human capillary blood sample for the screening of celiac disease
Schneller immunochromatographischer manueller Selbsttest zum qualitativen Nachweis von anti-deamidiertem Gliadin-Peptid IgG und Anti-Gewebe-Transglutaminase-IgA in menschlichen Kapillarblutproben zum Screening auf Zöliakie
Auto-test immunochromatographique manuel rapide pour la détection qualitative des anti-peptides désamidés de la gliadine IgG et des anti-transglutaminases tissulaires IgA dans un échantillon de sang capillaire humain pour le dépistage de la maladie cœliaque
Self Test immunocromatografico manuale rapido per la rilevazione qualitativa di IgG anti-peptide deamidato della gliadina e IgA anti-transglutaminasi tissutale (tTG) in campioni di sangue capillare umano per lo screening della celiachia
Snelle immunochromatografische handmatige zelftest voor de kwalitatieve detectie van anti-gedemideraliseerd gliadinepeptide IgG en anti-eweefsel transglutaminase IgA in humaan capillair bloedmonster voor de screening van coeliakie

Table with 2 columns: T1, T2. Rows include measurement range limit, SENSITIVITY, SPECIFICITY, ACCURACY.

* Maximum limit of the ThermoFisher Phadia ImmunoCAP EIA DGP-IgG and tTG-IgA reference system at which the PRIMA Test could be verified to be working perfectly.

** Clinical study conducted on 266 real samples
CROSS-REACTIVITY AND INTERFERING SUBSTANCES TESTED
To determine the effects of possible interferences, the substances listed in Tab. 1 were added to human serum samples, and the concentrations listed below were therefore added to physiological concentrations (if applicable).

Table with 2 columns: SUBSTANCE, CONCENTRATION. Lists cholesterol, glucose, bilirubin, etc.

No interference was found associated with CELIAC DISEASE SCREENING TEST.

Note: The result is reliable as long as the instructions are carefully followed and the test is stored under the specified conditions. However, false positive or negative results may be obtained, especially if: the device comes into contact with other liquids before being used, if the amount of blood and/or diluent is insufficient, if the number of drops added into the well is incorrect or if the reading time of the result is not respected.

WHAT TO DO IN CASE OF A POSITIVE RESULT? In the case of a positive result, you should consult your doctor, who will take the appropriate steps for a definitive diagnosis. Do not stop taking gluten until it is confirmed by a specialist.

WHEN IS IT RECOMMENDED TO USE THE TEST? CELIAC DISEASE SCREENING TEST can be performed if there are symptoms such as diarrhoea, abdominal pain and bloating, weight loss, tiredness and/or headache. Being a heritable condition, the test is also intended for those who have family members already suffering from the disease (risk of positivity of up to 20% for first-degree relatives, even if asymptomatic).

WHAT ARE THE TEST PERFORMANCE? SENSITIVITY, SPECIFICITY
The results of the CELIAC DISEASE SCREENING TEST were evaluated by comparing it with a quantitative commercial laboratory test used as a reference method (Thermo Fisher Phadia 250 - Thermo Fisher Scientific Inc.).

WHAT ARE THE TEST PERFORMANCE? SENSITIVITY, SPECIFICITY
The results of the CELIAC DISEASE SCREENING TEST were evaluated by comparing it with a quantitative commercial laboratory test used as a reference method (Thermo Fisher Phadia 250 - Thermo Fisher Scientific Inc.).

- 1. reliable if the instructions are followed properly.
2. Keep the test out of the reach of children.
3. Do not use the test after the expiry date or if the packaging and/or internal components are damaged. Do not freeze.
4. Follow the procedure exactly, using only the specified amount of capillary blood and test diluent.
5. Only use the test and its components once.
6. The test is for external use only. DO NOT SWALLOW.
7. After use, dispose of all components in accordance with local regulations.
8. If performing the test on other people, treat the blood samples as potentially infectious.
9. Only open the sealed aluminium sachet immediately prior to performing the test, pay attention to the opening mark on the sachet. The desiccant sachet must not be used. Discard it without opening it, by depositing it together with the waste.
10. No medically significant action should be taken without consulting your doctor first. This test is for screening purposes. A definitive diagnosis can only be made by a doctor after a joint evaluation of the clinical information and the results from further tests.
11. Do not confuse/mix components of this kit with components of other kits.
12. The CELIAC DISEASE SCREENING TEST DILUENT is only suitable for use with this test.

CONTENT
*The number of test devices in the package may vary. To find out the exact number of tests it contains, please refer to the "content" section on the outside of the package that you have purchased.

- 1* hermetically sealed aluminium pouch containing: 1 CELIAC RAPID TEST, 1 desiccant bag.
1* dropper vial containing the CELIAC RAPID TEST DILUENT sufficient for 10* test.
1* capillary pipette for collecting the blood sample.
2* sterile lancets.
1* alcohol cleansing gauze.
1* instruction leaflet.

WARNING: Read the entire procedure carefully before performing the test!
1) Wash your hands with soap and warm water, rinse with clean water and dry. Note: The use of warm water facilitates capillary blood sampling as it induces vasodilation.
2) Prepare the necessary material as follows: open the aluminium sachet, only take out the cassette and throw away the desiccant sachet. Open the plastic sachet containing the pipette and use the provided gauze to clean the puncture site. Allow the area to air dry.
3) Remove the protective cap from the lancet by carefully rotating it 360° without pulling it. Pull out and discard the released cap.
4) Position the hand downward and gently massage from the palm to the chosen phalanx repeatedly: middle or ring fingers are recommended - avoid the fifth finger. The puncture should be made slightly off-center, in the fleshy area of the thickest part of the phalanx. Press the lancet perpendicular to the fingertip, on the side from which the cap was removed, while keeping it firmly between the fingers.
5) The tip of the lancet automatically retracts safely after use. If the lancet does not work properly, use the second one provided. If the latter is not necessary, it can be disposed of without special precautions. Note: make a single puncture. Do not perform an immediate double sticking at the same site.
6) Wipe away the first drop of blood with the cleaning gauze and allow another drop to form, by keeping your hand downward and continue to massage the finger until the specimen volume is adequate.
7) Pick up the pipette without pressing the bulb. Two withdrawal methods are suggested:
-FIG. E1: hold the pipette at a slight incline without squeezing the bulb and place it in contact with the blood droplet, it will enter on its own by capillary action. Remove the pipette when it reaches the black line. If there is not enough blood, continue massaging the finger until the line is reached.
-FIG. E2: place the pipette on a flat, clean surface, allowing the tip to protrude, at which point the drop of blood should come into contact with the pipette, it will enter on its own by capillary action. Move your finger away when it reaches the black line. If there is not enough blood, continue massaging the finger until the line is reached.
8) Avoid, as far as possible, breaking contact between the blood and the pipette so as to prevent the formation of air bubbles.
9) Deposit the blood collected with the pipette into the sample well (S), by squeezing the bulb of the pipette.
10) Unscrew the blue cap from the dropper vial (leave the white cap tightly screwed on). Deposit 1 drop of diluent into the well indicated on the cassette (S).
11) Wait 10 minutes and then interpret the results as follows. Do not move the device while waiting.

INTERPRETATION OF THE RESULTS
READ THE RESULT AT 10 MINUTES.
DO NOT READ THE RESULTS AFTER 15 MINUTES.

POSITIVE*
In addition to the coloured line appearing at the C (Control) mark, the result is considered positive if:

- one coloured line appears near the T1 mark (positive for anti-DGP IgG), or

- one coloured line appears near the T2 mark (positive for IgA anti-tTG), or

- two coloured lines appear near the T1 and T2 marks (positive for both anti-DGP IgG and IgA anti-tTG)

Consult a doctor, as this result indicates probable celiac disease.
*NOTE: The colour intensity in the control and test line regions may vary (as can be seen in the pictures). Therefore, any colour intensity in the region of the T1 and T2 lines must be considered positive. The user is advised to pay close attention in detecting even weak intensity test bands, which must in any case be interpreted as a positive result.

NEGATIVE
Only one coloured line appears at the C (Control) mark. This result indicates that no anti-DGP IgG or IgA anti-tTG antibodies are present in the blood or that their level is below the limit detectable by the test.

INVALID
No control line appears at the C (Control) mark. The most likely reasons for an invalid result are that either there is insufficient sample volume or the procedure has not been performed correctly.

Review the procedure and repeat the test with a new device and a new sample.

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