Test your health

PRIMA LAB SA
Viale Serfontana 10
CH-6834 Morbio Inferiore
SWITZERLAND
support@primalabsa.ch
primalabsa.ch

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BIBLIOGRAPHY / LITERATURVERZEICHNIS / RÉFÉRENCES / BIBLIOGRAFIA / BIBLIOGRAFIE

- Larsen, John, Philip Buchanan, Sarah Johnson, Sonya Godbert, and Michael Zinaman. "Human Chorionic Gonadotropin as a Measure of Pregnancy Duration." *International Journal of Gynaecology and Obstetrics: The Official Organ of the International Federation of Gynaecology and Obstetrics* 123, no. 3 (December 2013): 189–95. <https://doi.org/10.1016/j.ijgo.2013.05.028>.
- Cole, Laurence A., Stephen DuToit, and Trefor N. Higgins. "Total HCG Tests." *Clinica Chimica Acta: International Journal of Clinical Chemistry* 412, no. 23–24 (November 20, 2011): 2216–22. <https://doi.org/10.1016/j.cca.2011.08.006>.
- Pike, Joanna, Sonya Godbert, and Sarah Johnson. "Comparison of Volunteers' Experience of Using, and Accuracy of Reading, Different Types of Home Pregnancy Test Formats." *Expert Opinion on Medical Diagnostics* 7, no. 5 (September 2013): 435–41. <https://doi.org/10.1517/17530059.2013.830103>.
- Gnoth, C., e S. Johnson. 2014. "Strips of Hope: Accuracy of Home Pregnancy Tests and New Developments." *Geburtshilfe und Frauenheilkunde* 74 (7): 661–69. <https://doi.org/10.1055/s-0034-1368589>.
- Moghazy, Dalia, Omar Al-Hendy, and Ayman Al-Hendy. "Krukenberg tumor presenting as back pain and a positive urine pregnancy test: a case report and literature review." *Journal of Ovarian Research* 7.1 (2014): 36.
- El Hage, Lea, and Betül Hatipoğlu. "Elevated hCG can be a benign finding in perimenopausal and postmenopausal women." *Cleveland Clinic Journal of Medicine* 88.11 (2021): 635–639.

SYMBOLS / SYMBOLE / SYMBOLES / SIMBOLOGIA / SYMBOLEN

	In vitro diagnostic device / In-Vitro-Diagnostikum / Dispositif médical de diagnostic in vitro / Dispositivo medico-diagnostico in vitro / In-vitrodiagnostisch hulpmiddel		Read the instructions before use / Beachten Sie die Gebrauchsanweisung / Consulter le mode d'emploi / Leggere attentamente le istruzioni per l'uso / Lees de instructies voor gebruik
	Sufficient for <=> tests / Ausreichend für <=> Tests / Suffisant pour <=> tests / Sufficiente per <=> test / Voldoende voor <=> tests		Temperature limits / Tempéaraturbegrenzung / Limites de température / Limiti di temperatura / Temperatuurlimieten
	Legal manufacturer / Hersteller / Fabricant / Fabricante / Juridische fabrikant		Do not reuse / Nicht wiederverwenden / Ne pas réutiliser / Non riutilizzare / Niet opnieuw
	CE marking / CE-Kennzeichnung / Marquage CE / Marchio CE / CE-markering		List number / Katalognummer / Code produit / Codice prodotto / Lijstnummer
	Importer / Importeur / Importateur / Importatore / Importeur		Lot number / Chargenbezeichnung / Code du lot / Lotto prodotto / Lotnummer
	Expiry date (last day of the month) / Verwendbar bis (letzter Tag des Monats) / Utiliser jusqu'au (dernier jour du mois) / Utilizzare entro (ultimo giorno del mese) / Vervaldatum (laatste dag van de maand)		Unique Device Identification / Eindeutige Geräteidentifikation / Identification unique du dispositif / Identificazione Unica del Dispositivo / Unieke apparatenidentificatie
	Authorized Representative in the European Community / Bevolmachtigter in der Europäischen Gemeinschaft / Mandataria dans la Communauté européenne / Mandataria nella Comunità Europea / Gemachtigd vertegenwoordiger in de Europese Gemeenschap		Intended for self-testing / Zur Selbstanwendung / Destiné à l'auto-test / Destinato all'autotest / Voor zelftest bedoeld
	Not for near patient testing / Produkt nicht für patientennahe Tests / Dispositif non destiné au diagnostic près du patient / Dispositivo non destinato ad analisi decentrali / Geen hulpmiddel voor patiëntnabije testen		

REVISIONS

REVISION	DATE	CHANGES
2.0	03/2026	New registered office of PRIMA Lab SA

EARLY DETECTION PREGNANCY TEST

Rapid immunochromatographic and manual self-test for the qualitative detection of the human Chorionic Gonadotropin (hCG) hormone in urine as an aid to the early identification of the pregnancy status in women

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EN - INSTRUCTIONS FOR USE

THE HCG HORMONE AND PREGNANCY

During the very first phases of pregnancy, the fertilised egg is carried from the Fallopian tubes to the uterus, where it implants after about seven days. In this phase, the tissue that will eventually form the placenta starts producing human chorionic gonadotropin (hCG hormone), which is essential for the development and continuation of the pregnancy. Its early appearance and continuation of the concentration in the urine make it an excellent marker to confirm the beginning of a pregnancy.

TEST PRINCIPLE

EARLY DETECTION PREGNANCY TEST is a manual self-test immunochromatographic device. The test contains monoclonal anti-hCG antibodies conjugated to colloidal gold nanoparticles in the conjugate pad and monoclonal anti-hCG antibodies immobilized in the test line region; both are considered active components of the device. During the test, if hCG is present in the urine sample, it binds to the gold-conjugated anti-hCG antibodies, forming an antigen-antibody complex. This complex migrates chromatographically along the membrane by capillary action and is captured by the immobilized anti-hCG antibodies in the test line area, producing a visible colored line. A second line in the control region serves as a procedural control and must always appear to confirm that the test has functioned properly.

F.A.Q. – QUESTIONS AND ANSWERS

WHEN THE TEST CAN BE PERFORMED?

The test can be performed starting from the 6th day before the period delay (5 days before you expect your period) and it is suggested to use the first urine of the morning. It is also recommended not to drink too much before performing the test. In case of irregular periods, it is advisable to perform the test on the day corresponding to the length of the longest menstrual cycle in recent months. Please refer to the Table 1.

WHAT ARE THE TEST PERFORMANCES?

SENSITIVITY*	SPECIFICITY*	OVERALL ACCURACY*
>99% (CI 95%: 100-95.91)	>99% (CI 95%: 100-84.54)	>99% (CI 95%: 100-96.65)

*clinical study with 111 clinical samples

CROSS-REACTIVITY AND INTERFERING SUBSTANCES TESTED

To determine the effects of possible interferences and cross-reactivity substances, the substances listed in Tab. 2 at the concentrations listed below were therefore added to negative urine sample. No interferences was found associated with device.

Table 2

SUBSTANCE	CONCENTRATION ASSAYED
Bilirubin	4.74 mg/dl
Glucose	459 mg/dL (25.5 mmol/L)
LH	500 mIU/mL
FSH	2000 mIU/mL
TSH	500 mIU/L
Uric acid	22.5 mg/dL (1340 μmol/L)
Haemoglobin	150 g/dL
Caffeine	20 mg/dL
Protein	1 mg/mL
Ascorbic acid	20 mg/dL
Gentisic acid	20 mg/dL
O-acetylsalicylic acid	20 mg/dL
Blood	2%
Semen	2%
urine pH	from 4 to 9
hCG B-core fragment	from 0.0001 nmol/mL to 0.5 nmol/mL
hCG alpha subunit	0.04 nmol/L to 40 nmol/Lh

Note: The device specifically identifies the beta subunit of the human hCG hormone and reveals if the concentration of that hormone in the urine is equal to or above 10 mIU/mL. The test is calibrated on the WHO International Standard NIBSC 18/244. The analytical sensitivity was verified according to CLSI EP17-A2 guidelines using the classical approach. The Limit of Blank (LoB) was established from 60 measurements on blank urine samples, with no false positives observed. The Limit of Detection (LoD) was verified at 10 mIU/mL using 40 replicates. A

minimum hit rate (positive results on total results) of 90% was considered acceptable, in line with EP17-A2 (Section 7.3), which allows this criterion when Type I and II error probabilities (α and β) are set at 0.05, corresponding to a 95% confidence level.

Situations that can lead to a "false positive" result: false positive results may occur in rare cases due to elevated levels of hCG in urine associated with certain pathological conditions. Elevated urine levels of hCG-βcf have been reported in approximately 49% of lung cancer cases and in 20% to 75% of patients with various gynecological cancers prior to surgery. Additionally, bladder cancer cell lines have been shown to express hCG, which may interfere with test results. False positives may also occur in cases of malignant pancreatic and biliary disease, Krukenberg tumor, and elevated hCG levels can be a benign finding in perimenopausal and postmenopausal women. In the case of unexpected or contradictory results, it is advisable to consult the gynaecologist, as some health problems may alter the test result.

Situations that can lead to a "false negative" result: miscalculation of the expected period, excessive dilution of urine due to taking too much fluid before the test.

HOW DO I INTERPRET THE TEST IF THE COLOUR AND INTENSITY OF THE TEST AND CONTROL BANDS ARE DIFFERENT?

The test should be considered positive regardless of how dark the red/pink colour is of the test and control lines. The lines should be even and full. NB: if the control line is visible, the appearance of a white line at the T (Test) mark, associated with the absence of any red/pink colouring, should always be interpreted as a negative result.

WHAT TO DO IN CASE OF A POSITIVE RESULT?

Please consult your GP or gynaecologist to confirm the pregnancy.

WHAT TO DO IN CASE OF A NEGATIVE RESULT?

A negative result means that you are not pregnant or that you performed the test too early (hormone levels at the beginning of pregnancy may vary from woman to woman). If you suspect a pregnancy, repeat the test in 2 or 3 days or when your period is due with a new device. If you still do not get your period with a negative result, consult your doctor.

If you need additional information, assistance or would like to report a problem, please contact us at support@primalabsa.ch. In case of serious incidents occur in relation to the device please inform the local health authority and manufacturer.

PRECAUTIONS





- In vitro diagnostic device for individual use.
- Read these instructions for use carefully before performing the test. The test is reliable only if all the instructions are followed correctly.
- Keep the test out of the reach of children.
- Do not use the test after the expiry date or if the packaging and/or internal components are damaged. Do not freeze.
- The test is for external and single use only. DO NOT SWALLOW.
- After use, dispose all components in accordance with local regulations.
- No medically significant action should be taken without consulting your doctor first.
- Only open the sealed aluminium sachet prior to performing the test. The desiccant sachet must not be used. Discard it, without opening it, by depositing of it together with the waste.

CONTENT

* The number of the testing devices of the kit may vary. For the exact number of tests contained, please refer to the "content" section on the outside of the bought package.

1* hermetically sealed aluminium pouch containing 1 EARLY PREGNANCY TEST, 1 desiccant bag; 1 instruction leaflet.

NOTE: Material required but not supplied: a timer.

PROCEDURE FOR USE

ATTENTION: Read the entire procedure before performing the test.

The test must be used immediately after opening the sealed aluminium pouch.

If the device has been stored at refrigerated temperature (4 °C), allow it to reach room temperature for at least 30 minutes before use.

It is suggested to use the first urine of the morning since the concentration of hCG is higher.

- Prepare the necessary material as follows: open the aluminium pouch, only take out the test device and dispose the desiccant sachet. Remove the protective cap (1) exposing the absorbing tip (2). –**FIG.A** Rotate the device so that the results window (3) faces up.
- Place the absorbing tip right under the flow of urine for at least 10 seconds, a few more doesn't interfere. –**FIG.B.1** **WARNING:** Ensure no drops of urine splash onto the reading window. If easier, urine can be collected in a clean, dry, chemically inert container (plastic or glass) with a wide mouth and stable base (not supplied) and the entire absorbing tip can be dipped into the urine for 10s as shown in –**FIG.B.2**. Avoid containers with residues of detergents, disinfectants, or other contaminants. Do not reuse food or drink containers. Use freshly collected urine samples.
- Place the protective cap back onto the absorbing tip. –**FIG.C** and put the test device onto a dry, flat surface with the results window facing up. –**FIG.D**
- Wait 5 minutes and then read the results as follows. Do not move the device while waiting.

INTERPRETATION OF THE RESULTS

READ THE RESULT AT 5 MINUTES
DO NOT READ THE RESULT AFTER 10 MINUTES

POSITIVE (PREGNANT)

Two coloured lines appear in the results window in the T (Test) and C (Control) regions.

Consult a doctor, as this result indicates a probable pregnancy.

NOTE: The colour intensity in the control and test line regions may vary (as can be seen in the pictures). Therefore, any colour intensity in the region of the T line must be considered positive. Strong positive results appear before the fifth minute.

NEGATIVE (NOT PREGNANT)

Only one coloured line appears next to the C (Control) mark.

This means that the concentration of hCG hormone in the urine is below the limit detectable.

This means that there is not an ongoing pregnancy or it's too early to detect hCG hormone in the sample.

INVALID

No control line appears at the C (Control) mark.

The most likely reasons for an invalid result are that either there is insufficient sample volume or the procedure has not been performed correctly.

Review the procedure and repeat the test with a new device and a new sample

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